making work easy
















- DE
- EN
- FR
- IT
- ES
- PT
- TR
- PL
- UK
- RU
- ZH
- JA
- KO
- AR











Vario jet

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









Made in Germany

21-9616 07062024

Vario jet

Nr. 2961-0000 / 2961-1000

DEUTSCH

1. Einleitung

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ENGLISH

1. Introduction

We are delighted that you have decided to purchase the *Vario jet* recycling sandblaster. This unit sets new standards with regard to function, efficiency, safety and ergonomics. Please note the following tips in order to ensure a long and problem-free service life.

2. Description / Applications

The *Vario jet* recyclable sandblaster is for use in dental laboratories in order to carry out the removal of residual investment material and oxides from frameworks (mainly model casting). It boasts a touch control pad that is not affected by dirt and is sealed against dust in order to ensure long and problem-free functioning.

The recyclable sandblaster can be set to manual or automatic mode dependent on choice. The construction of the unit ensures excellent blasting efficiency within a short period of time and uses comparatively less blasting material (on average 100 - 300 g per hour, dependent on extractor unit).

The funnel-shaped sand silo capacity of more than 5 kg allows long periods free from interruption. The parameters of working pressure, time, automatic / manual and extraction On / Off + run-on time (via an external unit) can be set on a sliding scale within the total potential.

2.1 Ambient Conditions

(in accordance with DIN EN 61010-1)

The *Vario jet* may only be operated:

- Indoors;
- Up to an altitude of 2.000 m [6.500 ft] above sea level;
- At an ambient temperature of 5 - 40 °C [41 - 104 °F]*);
- At a maximum relative humidity of 80 % at 31 °C [87.8 °F], decreasing linearly to 50 % relative humidity at 40 °C [104 °F]*);
- With a mains power supply where current fluctuations do not exceed 10 % of the nominal value;
- Under level 2 contamination conditions;
- Under overvoltage category II conditions.

*) Between 5 – 30 °C [41 - 86 °F], the unit can be operated at a relative humidity of up to 80 %. At temperatures between 31 - 40 °C [87.8 - 104 °F], the humidity must decrease proportionally in order to ensure operational readiness (e.g., at 35 °C [95 °F] = 65 % humidity; at 40 °C [104 °F] = 50 % humidity). The unit may not be operated at temperatures above 40 °C [104 °F].

3. Areas of Caution

Vario jet is an electrical appliance with a danger potential. This appliance may be connected / used only by authorized and qualified persons after checking its conformity with respect to the national norms.

Never operate a sandblaster without a suitable dust extraction unit or suitable personal protective clothing as this can be detrimental to your health. The type of dust extractor must relate to the type of dust produced by your sandblasting work. It is imperative to adhere to the legal regulations (in Europe this means EN 60335-2-69).

Only connect the unit to electrical sockets that are connected to the protective conductor system.

Before operating check the electrical and pneumatic leads for damage.

If defects are found do not activate or continue to use the unit.

Before opening the unit remove the plug from the socket.

Before changing the fluorescent lamp make sure it is unplugged.

Do not push the lamp too hard (danger of breakage!). Wear gloves or use a protective cloth!

Do not sandblast unprotected skin >>> danger of injury!

Before connecting the unit check the voltage against that on the unit ID.

The unit interrupts blasting when the protective glass is open. Do not readjust the switch. Escaping blasting material may cause eye injury.

It is the responsibility of the operator that national regulations during operation and regarding a repeated safety inspection of electrical equipment are complied with. For Germany these are the regulation 3 by DGUV (German Statutory Accident Insurance) in relation with VDE 0701-0702 (Association for Electrical & Electronic Technology).

Information on REACH and SVHC is available on our website www.renfert.com, in the Support area.

 Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the device, increase the risk of serious injury and lead to damage to the environment or the device itself.

3.1 Exclusion of liability

Renfert GmbH rejects any claims for compensation as well as warranty claims in case of:

- the product is used for any other purpose than those mentioned in the instruction manual.
- the product is changed or amended in any way beyond those mentioned in the instruction manual.
- the product is repaired by a non-authorized party or the product is not fitted with original Renfert replacement parts.
- the product is used continuously despite visible safety faults.

4. Installation / Operation

4.1 Set-up and Connection

1. Remove all the parts found within the blasting chamber and proceed as follows:
2. First decide whether the unit should be mounted on the wall or simply set on the work bench. For wall mounting, a drilling template and heavy load wall plug set are supplied. Anchor the hooks on the wall. Before hanging up the unit set up the connection leads according to points 3 to 5. Mount the unit onto the hooks via the openings on the back.
3. Connect the compressed air to the unit via the blue connection hose and the connection set (Fig. 1 - 7). Ensure that the compressed air hose is not kinked (Fig. 8). Then connect the foot switch to the unit (Fig. 9).
4. Set up the electrical connection to the unit using the cable supplied (Fig. 10 - 11).

 Please pay attention to the voltage specified on the identification plate!

5. Connect the extraction hose to the extractor pipe at the back of the unit (Fig. 12). When using a single-bench extraction unit connect the electrical cable of your extraction unit to the plug on the back of the unit via the supplied adapter cable (Fig. 13-A). Set the extraction unit to continuous use mode on the main switch.

 When connecting several sandblasters please note the maximum allowable connection power for extraction unit of 1200 W!
Your extraction unit will now be driven by the blasting unit.

6. Attach the supplied blasting basket to the motor shaft by pressing down (Fig. 14). Make sure that the supplied hexagonal rubber mat is placed in the basket (Fig. 16).
7. Remove the cable fastener (used to secure during transport) from the blasting nozzle unit (Fig. 15-A).
8. Fill the unit with approx. 5 kg of blasting material of your choice by simply pouring it into the blasting chamber (Fig. 17).

 **Maximum grain size should be 500 µm [32 mesh]!**
With the use of abrasives > 250 micron (for ex. 500 micron), a significantly higher level of wear and tear must be expected.

9. Close the glass lid as the unit will only work when it is closed.
10. Check the required blasting pressure on the gauge (Fig. 18-B) and change the pressure if necessary via the pressure reducer (Fig. 18-A). Pull the knob out of the raster position to reset the blasting pressure.
11. Your unit is ready for use.
12. The use of additional compressed-air based units (jet nozzles, fine sandblasting unit, deflasking chisel etc.) within the blasting chamber is possible via the feeder (Fig. 19). The unit can be connected to the free adapter (Fig. 20-B).


5. Operation

5.1 Quick overview of controls (Fig. 18)


- On / Off switch – selection button (G)
- Selection switch for manual / automatic mode (E)
- Blasting time setting (automatic mode) – selection buttons from 1 to 60 minutes (C)
- Manual On / Off switch for external extraction unit – selection button (F)
- Setting for extraction unit run-on time – dial (Fig. 21)
- Setting for blasting pressure – dial (A) from 1 to 6 bar [14.5 - 87 psi]
- Start button for automatic mode – selection button (H)
- Stop button for automatic mode – selection button (J)
- Pause button for automatic mode – selection button (I)

5.2 Automatic mode

1. Put the cast objects into the basket and set the required blasting pressure on the control (Fig. 18-A). The blasting nozzle unit (Fig. 15-A) must be fixed to the holder.
2. Press the button (Fig. 18-G) to switch on the unit. It will go into automatic mode. This is also visible via the time displayed in the LED display.

- The time can be set from 1 to 60 minutes as required via the button (Fig. 18-C). This will be stored in the memory of the unit.
- By pressing "START" the *Vario jet* will automatically begin the blasting process i.e., the basket will start turning, the extraction unit will switch on and the unit carries on for the programmed amount of time. The time display flashes.
- The blasting process can be interrupted by pressing "HOLD" (Fig. 18-I). The remaining time will flash on the display. On re-newing the start blasting will continue for the remaining time. During the pause you can work manually using the foot switch.
- By pressing "STOP" (Fig. 18-J) the blasting process will cease regardless of whether you wanted to interrupt the process or end it completely. The timer returns back to the original time setting and is ready for use.

Tips & Tricks:

- Before blasting cast objects remove as much investment material as possible >>> reduces necessary blasting time
- Remove any casting sprues >>> 50 % improvement on blasting effect
- When blasting single objects place in casting cones >>> improves object rotation

5.3 Manual Mode

- Press the button (Fig. 18-G) to switch on the unit, then press the button (Fig. 18-E) to switch the unit into manual mode. The LED display shows only the points.
- The unit can now only be used via the foot switch. The basket will not function in manual mode.
- The extraction unit will automatically start up when blasting commences.
- By re-pressing the button (Fig. 18-E) the unit will switch back into automatic mode.
- The blasting nozzle unit (Fig. 15-A) can be removed from its holder in manual mode and be guided by hand.

5.4 Extraction Unit

- The run-on time for the connected (Fig. 13-A) extraction unit can be infinitely varied between 5 and 30 seconds via the dial on the back of the unit (Fig. 21).
- By pressing the button (Fig. 18-F) you can make the extraction unit independent from the blasting process e.g., when cleaning the unit. During the blasting process the extraction unit cannot be manually switched off for safety reasons.

- On initial delivery the extraction unit is automatically configured for Renfert's extraction units. You can operate Renfert's extraction units on the lowest setting.

Tips & Tricks:

- During the extraction unit's run-on time the set time cannot be altered.
- During the extraction of the blasting dust *Vario jet* may build-up a static charge via the extraction hose. This can be avoided, if the extraction hose contacts an earthed metal object for ex. a heating radiator at the installation.

5.5 Warning Function

The red LED (Fig. 18-D) lights up when the glass screen is open >>> **no blasting is possible!** The display flashes.

5.6 Protective Ring on Basket

At the inner side the basket is fitted with a protective rubber band which prevents interlocking of the clasp pieces in the perforated sides as well as early wear and tear of the basket. It is possible to remove the band, however it must be reckoned on a higher wear and tear of the basket.

5.7 Removal of the Basket

The basket can be easily removed for manual use of the unit. Simply pull it forward off the shaft. Pull off using both hands supported by a gentle shaking movement. Now remove the basket from the blasting chamber. To reposition press the basket onto the shaft until it clicks in place. For improved central positioned the rubber mat can be removed before repositioning. The shaft is visible through an opening in the basket (Fig. 22).

5.8 Distance between Jet Nozzle Unit and Basket

The distance between the blasting nozzle unit (Fig. 15-A) and the rotative basket, and thus to the objects to be blasted, can be set in four levels as required. The standard level positions the blasting nozzle unit in a central position. To set for specific uses pull the blasting nozzle unit out of its holder and position it according to your needs (Fig. 23).

Tips & Tricks:

The blasting nozzles can be set to four levels:

- Position 1: least distance to basket >>> for regular blasting of 1 - 2 objects
- Position 2 - 3: middle distance to basket >>> for regular blasting of 3 - 5 objects
- Position 4: greatest distance to basket >>> for regular blasting of 5 or more objects

6. Cleaning / Maintenance

From time to time your *Vario Jet automatic* will need maintenance. This consists of checking the connection leads and pipes for seals and damage as well as cleaning the plate. Use only solvent-free cleaning fluid when cleaning (e.g. soapy water).

6.1 Condensed Water

Occasionally check whether water has formed in the maintenance unit. If so, press the drainage valve on the maintenance unit (Fig. 20-A).

6.2 Removal of blasting residue

A sieve which is divided into two sections is used for retaining the blasting material residues. (Fig.24).

Removing the blasting material:

1. Remove the basket (see chapter 5.7).
2. Remove the left hand sieve. Blasting material residues can be carefully emptied into the right hand sieve. (Fig. 25).
3. Remove the right hand sieve and dispose of the blasting material residues. (Fig. 26).
4. Replace the right hand sieve.
5. Replace the left hand sieve.
6. Insert the basket.

6.3 Changing the Blasting Material

To change the blasting material stand the unit on the table and pull it slightly over the edge of the table until the drainage lid on the bottom of the unit is visible (Fig. 27).

Attention: risk of tilting over!

Place a suitable vessel underneath this opening (Fig. 28) and loosen both knurled screws to let out the sandblasting material (Fig. 29).

Tighten the lid again using the knurled screws and check that the sieve at the base of the sand silo has no impurities in it (Fig. 30) >>> if it does remove them. Fill with new sandblasting material (Fig. 17).

Tips & Tricks:

Long-term tests have shown that optimum blasting efficiency is achieved with 125 µm [115 mesh] grain size.

6.4 Changing the Nozzles

- Remove the sandblasting nozzle from the holder.
- Loosen the knurled screw (Fig. 31).
- Pull the nozzle out (Fig. 32).
- Inspect the mixing chamber for wear and remove any abrasive residue (Fig. 33).
- Insert the new nozzle.
- Retighten the knurled screw hand-tight.

Note:

A setscrew was used in previous versions instead of the knurled screw. Loosen the setscrew using the Allen key included in the accessory kit and retighten it HAND-TIGHT after insertion of the new nozzle (Fig. 31a).

6.5 Changing the Viewing Pane

Loosen the crosshead screws on the top hinge (Fig. 34). Remove the pane (Fig. 35). The new pane can now be easily mounted by feeding it through the upper holder from one side and then carefully tightening the crosshead screws.

Make sure the coated side of the pane is facing into the chamber!

Note: The slight dents in the glass coating will disappear after a short time.

Tips & Tricks:

Occasionally clean the inside of the pane using a damp cloth in order to remove the cloudiness caused by fine dust.

6.6 Changing the Fluorescent Lamp

Unplug the mains cable!

Pull down the fluorescent lamp out of its clamps (Fig. 36 - 37). Now completely remove the protective pipe from its socket (move to the right) and take out the fluorescent lamp (Fig. 38). Replace the lamp by carrying out the above process in reverse. Make sure that the fluorescent pipe is fed correctly into the rubber holder in the protective pipe (Fig. 39) and that the reflector is positioned correctly. In order to guarantee perfect functioning of the unit you should also make sure that no blasting material get into the lamp. Only use original lights with inte-grated starter and 11 Watt power (see replacement parts list)!

When changing the lamp do not bend it (danger of breakage!). Wear gloves or use a cloth to avoid glass splinters.

Do not pull the power supply cable of the lamp!

6.7 Changing the Hose within the Blasting Chamber

The blasting hoses are prone to natural wear and tear and must, therefore, occasionally be replaced. Disconnect the *Vario jet* from the compressed air supply. Loose the hose clip and insert new hoses (Fig. 40 - 41).

Only use Renfert replacement hose!

7. Spare Parts

You can find components subject to wear and the spare parts on the spare part list in the internet at www.renfert.com/p918.

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

Serial number and date of manufacturing are shown on the type plate of the unit.

8. Guarantee

When used properly Renfert guarantees all parts of the *Vario jet* - with the exception of wear and tear parts - **for 3 years**. In case of any claims for warranty original dealers' invoice is required. Parts which are subject to natural wear and tear (wear parts) and consumables are excluded from the guarantee. These parts are marked in the spare part list. The guarantee is void for damage caused by misuse, disregard of the instruction manual, cleaning, maintenance and connection, for self-repaired units or repair by untrained personnel, for use of spare parts from other manufacturers and on influences outside the normal practice or unauthorised in the instruction manual. Services under guarantee do not extend the guarantee.

9. Technical data

Working pressure:	1 - 6 bar [14.5 - 87 psi]
Connecting pressure, external:	6 - 10 bar [87 - 145 psi]
Air consumption:	200 l/min at 6 bar [7.06 cfm at 87 psi]
Voltage (pre-set by Renfert):	220 - 240 V, 50 / 60 Hz 100 - 120 V, 50 / 60 Hz
Power consumption:	30 W
Lamp power:	11 W
Cable length:	2 m [78 inch]
Ø Compressed air hose:	
interior:	4 mm [157 mil]
exterior:	6 mm [236 mil]
Ø Size of pipe union for external extraction unit:	
interior:	35 mm [1.38 inch]
exterior:	40 mm [1.57 inch]
Dimensions W x H x D:	405 x 650 x 460 mm [15.9 x 25.6 x 18.1 inch]
Blasting chamber volume:	25 l [6.6 US gal]
Weight (unfilled), approx.:	26 kg [57.3 lbs]
LpA *) (at max. volume flow):	< 72 dB(A)
This product contains a light source of energy efficiency class:	G

*) Sound pressure level in accord. with EN ISO 11202

10. Parts supplied

- 1 *Vario jet* blasting unit
- 1 Foot switch, electric
- 1 Connector lead
- 1 Connector pipe for compressed air system
- 1 Connector set
- 1 Small parts set with drilling template
- 1 Instruction manual
- 1 Magnetic foil

11. Units available

No. 2961-0000	<i>Vario jet</i> , 230 V
No. 2961-1000	<i>Vario jet</i> , 120 V

12. Accessories

Renfert supplies the following abrasive materials suitable with the *Vario jet*:

Cobra aluminium oxide

No. 1594-1105	25 µm [500 mesh], white 5 kg canister
No. 1594-1205	50 µm [270 mesh], white 5 kg canister
No. 1594-2220	50 µm [270 mesh], white 20 kg bucket
No. 1584-1005	90 µm [170 mesh], white 5 kg canister
No. 1583-1005	110 µm [150 mesh], white 5 kg canister
No. 1583-1020	110 µm [150 mesh], white 20 kg bucket
No. 1587-1005	125 µm [115 mesh], pink 5 kg canister
No. 1587-1020	125 µm [115 mesh], pink 20 kg bucket
No. 1585-1005	250 µm [60 mesh], white 5 kg canister
No. 1585-1020	250 µm [60 mesh], white 20 kg bucket

Rolloblast glass beads

No. 1594-1305	50 µm [400 - 200 mesh] 5 kg canister
No. 1594-2312	50 µm [400 - 200 mesh] 12.5 kg canister
No. 1589-1005	100 µm [170 - 100 mesh] 5 kg canister

13. Disposal instructions for countries in the EU

To conserve and protect the environment, prevent environmental pollution and improve the recycling of raw materials, the European Commission adopted a directive that requires the manufacturer to accept the return of electrical and electronic units for proper disposal or recycling.

Within the European Union units with this symbol should not therefore be disposed of in unsorted domestic waste:

For more information regarding proper disposal please apply at your local authorities.

14. Troubleshooting Guide

Fault	Cause	Remedy
Unit doesn't blast, warning lamp flashing.	<ul style="list-style-type: none"> • Glass pane is not closed. • Magnet (Fig. 42-A) is missing. • Protective switch is defect. 	<ul style="list-style-type: none"> • Close the pane. • Replace the magnet (see spare parts list). • Exchange during service.
Air flow present but no sand coming through.	<ul style="list-style-type: none"> • Amount of sand too low. • Filter on suction point of suction pipe is blocked. • Suction pipe defect. 	<ul style="list-style-type: none"> • Fill with more sand. • Clean the filter. • Exchange the suction pipe.
No air, no sand.	<ul style="list-style-type: none"> • Nozzle is blocked. 	<ul style="list-style-type: none"> • Open the blasting nozzle unit (Fig. 15-A) and clean the nozzle.
On/Off switch doesn't work.	<ul style="list-style-type: none"> • No power connection. • Unit fuse (Fig. 13-B) is defect. 	<ul style="list-style-type: none"> • Connect to the power supply. • Change fuse (Fig. 13-B) - replacement fuse in set. If it happens again >>> repair.
Basket shakes, doesn't move.	<ul style="list-style-type: none"> • Basket is stuck mechanically. 	<ul style="list-style-type: none"> • Release by removing and replacing the basket.
Blasting chamber gets too dusty.	<ul style="list-style-type: none"> • Too little suction power. • Extraction unit blocked. 	<ul style="list-style-type: none"> • Raise suction power on extraction unit. • Increase suction performance by closing the air inlet boreholes under the glass holder (Fig. 42-B) (with the supplied magnetic strip). • Clean extraction unit.
Consumption of blasting material too high.	<ul style="list-style-type: none"> • Suction power is too high. 	<ul style="list-style-type: none"> • Reduce suction power. Warning: too little suction power can lead to dirty blasting material.
Objects fall out of the basket.	<ul style="list-style-type: none"> • Too many objects. 	<ul style="list-style-type: none"> • Reduce to a maximum of 6 objects.
Blasting time too long / reduced blasting effect.	<ul style="list-style-type: none"> • The blasting cones are sheltering the objects. • Too little working pressure. • Too great distance from nozzle. • Too little blasting material in unit. • Blasting material lead pipe is blocked / defect. • Blasting nozzle is split. • Air power of the compressor is too low. 	<ul style="list-style-type: none"> • Remove blasting cones before starting. • Increase working pressure >>> recommended working pressure: 4 - 5 bar [58 - 72.5 psi]. • Move nozzle further forward. • Unit should be filled with min. 4 - 5 kg blasting material. • Remove sand and check / replace lead pipe. • Replace nozzles. • Dehydrate the compressor. Use a more powerful unit.
Extraction unit still runs even when blasting has been interrupted.	<ul style="list-style-type: none"> • A certain run-on time is normal in order to remove any dust still present in the chamber. • The switch for the external extraction unit (Fig. 18-F) was reactivated. 	<ul style="list-style-type: none"> • Regulate the run-on time between 5 and 30 seconds via the control on the back of the unit (Fig. 21). • Switch off the external extraction unit.
Heavy dust build-up during blasting.	<ul style="list-style-type: none"> • No extraction unit connected. • Extraction power too weak. 	<ul style="list-style-type: none"> • Connect an extraction unit. • Close the air holes (Fig. 42-B). • Increase extraction power.
Switch for external extraction does not work.	<ul style="list-style-type: none"> • Extraction unit is not controlled by blasting unit. 	<ul style="list-style-type: none"> • Connect the extraction unit power cable to the blasting unit.
Foot switch doesn't work.	<ul style="list-style-type: none"> • Unit is in automatic mode. 	<ul style="list-style-type: none"> • Interrupt automatic mode using HOLD button (Fig. 18-I). • Switch to manual (Fig. 18-E).

Fault	Cause	Remedy
Extraction unit doesn't work.	<ul style="list-style-type: none"> • Extraction unit is not controlled by the sandblasting unit. • Extraction unit is not switched to continuous mode. • Fuse for extraction unit (Fig. 13-C) is defect. 	<ul style="list-style-type: none"> • Connect the extraction unit power cable to the blasting unit. • Switch the extraction unit to continuous mode. • Change the fuse (Fig. 13-C) - replacement fuse found in set. If it happens again >>> repair.
Pulsating sand flow.	<ul style="list-style-type: none"> • Blasting material is used or too impure. 	<ul style="list-style-type: none"> • Renew blasting material.
Unusually large amount of blasting material in blasting sand.	<ul style="list-style-type: none"> • Sieve is in an incorrect position 	<ul style="list-style-type: none"> • Replace the sieve correctly, observe the insertion sequence: <ol style="list-style-type: none"> 1. Insert right hand sieve. 2. Insert left hand sieve.
Independent switching off while working.	<ul style="list-style-type: none"> • Too low supply voltage. • Inadequate grounding because the device is located too far from the nearest ground connection. 	<ul style="list-style-type: none"> • Checking of the power capacity of the mains supply. • Temporary switching off of other more powerful units. • Decrease the distance to the nearest ground connection

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c/o Renfert GmbH
Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany