making work easy













- de
- en
- fr
- it
- es
- pt
- tr
- ru
- pl
- zh
- ja
- ko
- ar
- uk
- hr












Vario basic

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

















Made in Germany

21-6575 08102024

Vario basic

Nr. 2960-0005 / 2960-1005

DEUTSCH

1. Einleitung

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ENGLISH

1. Introduction

We are pleased with your decision to purchase the *Vario basic*.

This unit sets new standards with regard to multifunctionality, performance, and ergonomics.

Please read the following operating instructions carefully and observe all the information they contain in order to ensure long and trouble-free service.

2. Application Area

The *Vario basic* recyclable blaster is used in dental laboratories to remove investments and oxides from cast parts (primarily model castings).

The unit can be fitted with up to 4 fine blasting tanks. Thus the application area also includes the removal of investments and oxides from crowns / bridges and inlays as well as surface treatment by polishing.

The unit is supplied without the fine blasting tanks. It has the capacity to hold up to 4 fine blasting tanks. As an alternative to a sandblasting tank, it is also possible to connect a device which runs on compressed air. This is determined by the selection dial within the blasting chamber.

If micro sandblasting tanks are to be installed from the outset, this should be done prior to the unit being filled with abrasive.

2.1 Ambient Conditions

(in accordance with DIN EN 61010-1)

The *Vario basic* may only be operated:

- Indoors;
- Up to an altitude of 2,000 m [6,500 ft] above sea level;
- At an ambient temperature range of 5 - 40 °C [41 - 104 °F];
- At a maximum relative humidity of 80 % at 31 °C [87.8 °F], dropping linearly to 50 % relative humidity at 40 °C [104 °F];
- With a mains power supply whose voltage fluctuations do not exceed 10% of the rated value;
- Under contamination level 2 conditions;
- Under over-voltage category II conditions.

*) Between 5 - 30 °C [41 - 86 °F] the unit can be operated at a relative humidity of up to 80 %. At temperatures between 31 - 40 °C [87.8 - 104 °F] the humidity must decrease proportionally in order to ensure proper operability (e.g., at 35 °C [95°F] = 65 % relative humidity, at 40 °C [104 °F] = 50 % relative humidity). The unit may not be operated at temperatures above 40 °C [104 °F].

3. Hazard Information

3.1 Employed Symbols

The following symbols are used in these instructions and on the unit:

Danger

This indicates an imminent risk of injury.

Electrical hazard

There is a risk of electrical shock.

Attention

Failure to observe the associated information can result in damage to the unit.

Note

The operating instructions contain useful information to make handling easier.

Wear eye protection.

Wear protective gloves.

For indoor use only.

Note the operating instructions.

3.2 Hazard Information

- The unit should only be operated using a mains cable with a country-specific plug system. Any modification required should only be carried out by an electrician.
- Only connect the unit to electrical sockets that are connected to the protective conductor system.
- Prior to commissioning, compare the information on the nameplate with the specifications of the local mains power network and the connection values for compressed air.
- In order to avoid health hazards, never operate the sandblasting unit without adequate dust extraction and proper personal protection. In order to protect the health of persons working in a dental laboratory, the law requires that specialised dust extractors must be used. National regulations state the MAK value - maximum concentration of suspended particles in the air. Find out about the legal national threshold value and the type of dust generated in your laboratory.
- Improper operation can result in the risk of eye and skin injuries.

- Always wear appropriate eye protection when operating the unit! Failure to wear eye protection can result in eye injuries due to airborne particles.
- Never direct the nozzles towards your eyes or uncovered areas of skin!
- Never work with the viewing pane open!
- Regularly inspect all connecting lines and hoses (e.g., power cord and compressed air hoses) for damage (e.g., kinks, cracks, porosity).
- Units with damaged connecting lines or hoses, or those which exhibit other defects may not be used.
- Always disconnect the unit from the compressed air supply prior to performing any maintenance tasks.
- Always disconnect the unit from the mains power supply prior to performing any work on electrical components.
- Do not activate the foot switch when only the compressed air inlet hose is connected. The unconnected hose can whip about, resulting in a significant risk of injury.
- It is the responsibility of the operator that national regulations during operation and regarding a repeated safety inspection of electrical equipment are complied with. For Germany these are the regulation 3 by DGUV (German Statutory Accident Insurance) in relation with VDE 0701-0702 (Association for Electrical & Electronic Technology).
- Information on REACH and SVHC is available on our website www.renfert.com, in the Support area.
- Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the device, increase the risk of serious injury and lead to damage to the environment or the device itself.

When using micro sandblasting tanks:

- Check the tank cover for a secure seat prior to starting operation. Tank covers that are not securely closed can explosively open. Ejected parts and abrasive represent a significant risk of injury.
- Always wear safety glasses when working with the tank (filling, cleaning, maintenance).
- Abrasive residue on the seal will result in leaks and premature seal wear. Always clean the threads and seal after filling and handtighten the cover.
- Solvents and tensides can result in the formation of microscopic cracks in the plastic (explosion hazard!). Use only a clean, dry cloth to wipe the tank out and off. Never write on the tanks or place adhesive labels on them.
- Regularly inspect the tanks and tank covers for damage and immediately replace any suspect components.

3.3 Liability Exclusion

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- The product is employed for any purposes other than those cited in the operating instructions;
- The product is altered in any way other than those alterations described in the operating instructions;
- The product is repaired by other than an authorized facility or if any but Renfert OEM parts are employed;
- The product continues to be employed, despite obvious safety faults or damage;
- The product is subjected to mechanical impacts or is dropped.

4. Commissioning

4.1 Wall Mounting

The holder required for wall mounting can be ordered separately as an optional accessory (refer to the accessories list). The holder includes instructions for mounting the unit on a wall.

4.2 Setup and Connection

1. Remove the foot switch stored in the blasting chamber.
2. Connect the foot switch to the sandblasting unit by completely inserting the two red and yellow pneumatic hoses into their correspondingly marked connectors until they reach their stops (Figure 1a).
3. Connect the compressed air supply to the sandblasting unit:
 - Loosen the air connection's knurled sleeve nut on the unit (Figure 1b).
 - Push the knurled sleeve nut over the included blue compressed air connecting hose (Figure 1c).
 - Push the compressed air connecting hose onto the compressed air connection and tighten the knurled sleeve nut (Figure 1d, 1e).
4. Connect the compressed air hose to the compressed air supply:
 - Loosen the compressed air connector's knurled sleeve nut and push it over the included blue compressed air connecting hose (Figure 2a).
 - Push the compressed air connecting hose onto the compressed air connection and tighten the knurled sleeve nut (Figure 2b, 2c).
 - Connect the connector to the laboratory compressed air supply (Figure 3).

Do not bend or kink the hose (Figure 4)!

4.3 Installing the Nozzle

1. Remove the base grate (Figure 5).
2. Remove the transportation safety device (cable clip) from the blasting nozzle (Figure 6).
3. Check the extraction hose for proper seat (Figure 7).

4.4 Adjusting the Operating Pressure

1. Set the selector switch to the recycling sandblasting position (Figure 8, Pos. A).
2. Adjust the blasting pressure to 1 - 6 bar [14.5 - 87 psi] (Figure 9). The current pressure setting is continuously displayed on the manometer.

4.5 Adding Abrasive

1. Remove the base grate (Figure 5).
2. Add abrasive (Figure 10). Approx. 5 kg (+/- 2 kg).
3. Reinsert the base grate (Figure 5).

i Use only clean, dry abrasive with the appropriate grain size (refer to the accessories list). We recommend an abrasive with a grain size of 125 - 250 µm [115 - 60 mesh].

! Max. grain size: 500 µm [32 mesh].
With the use of abrasives > 250 micron (for ex. 500 micron), a significantly higher level of wear and tear must be expected

4.6 Starting the Unit Illumination

1. Connect one end of the power cord to the unit (Figure 11) and plug the other end into the wall outlet (Figure 12).
2. Switch the lamp on (Figure 13).

If the membrane key is damaged or defect, please disconnect the device and send away for repair! If the membrane is removed there is the risk of an electrical shock!

4.7 Connection to the Dust Extractor

Push the extraction hose onto the extractor connection (Figure 14, Pos. a). An extraction hose is not included in the standard *Vario basic* delivery package. When using extraction units with an „automatic start function“, ensure that the *Vario basic* mains cable is plugged into the relevant socket on the extraction unit. Observe the maximum connection load on the plug socket.

5. Operation

5.1 Device Selection

The *Vario basic* is supplied without the fine blasting tanks. It has the capacity to hold up to 4 fine blasting tanks. As an alternative to a sandblasting tank, it is also possible to connect a device which runs on compressed air. This is determined by the selection dial within the blasting chamber. (Figure 8):

- (A) Recycling blasting function
- (B) Selection of various consumers (optional)
 - (B1) Top right fine blasting tank
 - (B2) Top left fine blasting tank
 - (B3) Bottom right fine blasting tank
 - (B4) Bottom left fine blasting tank

Sandblasting continues as long as the foot switch is depressed (Figure 15). A small amount of abrasive continues to flow once the foot switch is released.

Please use the provided gloves when sandblasting.

The cloth hand port cuffs do not need to be removed.

Do not operate the unit if the viewing pane is open.

5.2 Positioning the Blasting Nozzle

Nozzle distance and nozzle angle can be adjusted separately (Figure 16).

5.2.1 Adjusting the Nozzle Distance:

Three different nozzle inserts are available to provide three different nozzle distances (Figure 16, Pos. a + Figure 6).

5.2.2 Adjusting the Blasting Angle:

1. Loosen the upper and lower knurled nuts on the retainer (Figure 16, Pos. b).
2. Adjust the nozzle to the desired angle.
3. Retighten the screw (Figure 16, Pos. b).

5.3 Adjusting the Extraction Power

The suction performance of a weak extraction unit can be increased by closing the boreholes on the back side of the unit (Fig. 14, Pos. b) with the supplied magnetic foil.

5.4 Cold Silanization in the Sandblasting Unit

EN

The *Basic* series of Renfert sandblasting units is exclusively recommended by 3M ESPE for use with ROCATEC™. This bonding system is ideal if an adhesive bond is required between a composite and a dental material (metal, acrylic, porcelain, zirconia porcelain).

Advantage:

The surface film is applied cold, preventing thermal loading of metal frameworks and consequently the risk of distortion. All units in the *Basic* series can be used for ROCATEC™ due to their regular, homogeneous jet flow.

6. Connecting Additional Tanks and Other Compressed Air Devices

6.1 Installing Micro Sandblasting Tanks

The installation of the additional tanks is described in a separate set of instructions included with the tanks.

6.2 Connecting Other Devices

The various devices are connected to the connectors on the back of the unit (Figure 17):

- 1 Installation of top left blasting tank (Before Installation, remove the blank plug)
- 2 Installation of top right blasting tank (Before Installation, remove the blank plug)
- 3 Installation of bottom left blasting tank (Before Installation, remove the blank plug)
- 4 Installation of bottom right blasting tank (Before Installation, remove the blank plug)

i *Alternatively, it is possible to connect a device which runs on compressed air!*

6.3 Filling the Micro Sandblasting Tanks

Never activate the foot switch while filling the tanks.

Regardless of the pressure indicated on the manometer, the micro sandblasting tanks are only under pressure if the foot switch is activated.

1. Unscrew the tank cover (Figure 18).
2. Add abrasive only to the max. fill mark (the point at which the tank jacket thickens below the cover thread) (Figure 19).
3. Reinstall the tank cover (Figure 18).

Before closing the tank, always clean any abrasive residue from the threads on the cover and on the tank. Do not misalign the cover.

Use only clean, dry abrasive with the appropriate grain size (refer to the accessories list).

6.4 Connecting Additional Compressed Air Devices

Additional compressed air devices such as the *Basic mobil* are connected to the back of the *Vario basic* (Figure 21, Section 6.2, Figure 17).

The connected device's compressed air hose can be inserted through the socket on the side (Figure 20).

Figure 21: Example of connecting a *Basic mobil*.

7. Cleaning / Maintenance

Always disconnect the unit from the mains power and compressed air supplies prior to performing any cleaning and maintenance tasks.

7.1 Cleaning the Interior

- Use only solvent-free cleansers (e.g., soapy water).
- Remove the base grate (Figure 5) and vacuum out the blasting chamber.
- Using a dry, lint-free cloth, wipe out the insides of the tanks and the tank covers!

7.2 Inspecting the Water Separator

- Drain any condensation water through the outlet valve.
- Press up on the valve from below (Figure 22).

7.3 Blasting Nozzle Replacement

- Remove the sandblasting nozzle from the holder.
- Loosen the knurled screw (Figure 23).
- Pull the nozzle out (Figure 24).
- Inspect the mixing chamber for wear and remove any abrasive residue (Figure 25).
- Insert the new nozzle (Figure 24).
- Retighten the knurled screw hand-tight.

i *A setscrew was used in previous versions instead of the knurled screw. Loosen the setscrew using the Allen key included in the accessory kit and retighten it **HAND-TIGHT** after insertion of the new nozzle (Figure 23a).*

7.4 Changing the Abrasive

- Pull the unit out over the edge of the table (Figure 26).

Tipping hazard!
A second individual should help steady the unit.

- Position the container to catch the used abrasive (Figure 27).
- Open the abrasive drain hole (Figure 28).
- Clean / and resecure the filter (Figure 29 / Figure 7).
- Close the abrasive drain hole (Figure 28).
- Add new abrasive (Figure 10).

7.5 Replacing the Viewing Pane

- Loosen the retaining screws (Figure 30).
- Remove the old pane and install the new one (Figure 31).
- Carefully retighten the retaining screws (Figure 30).

7.6 Replacing the Blasting Hoses

Regularly inspect the blasting hoses for damage and wear. Replace damaged hoses immediately.

Disconnect the sandblasting unit from the compressed air supply when replacing the blasting hoses (Figure 3).

- Drain the abrasive out of the unit (refer to Section 7.4).
- Loosen the hose clamp and pull the hose off the connector (Figure 32, 33, 34).
- Install new hoses (Figure 32, 33, 34).

Use only Renfert OEM replacement hoses!

7.7 Replacing the Fluorescent Lamp

7.7.1 Device with LED lighting

The service life of LED lighting is far higher than that of a fluorescent tube. Therefore, the LEDs can not be changed individually, the complete lighting has to be changed (see chapter 8, Spare Parts).

7.7.2 Device with fluorescent tubes

In these devices, the fluorescent tubes can be changed.

Danger of breakage!

Do not press on the fluorescent lamp or try to bend it when replacing it. We recommend using gloves or a cloth to protect your hands in case of breakage.

- Switch off the unit.
- Pull plug out.
- Lightly press light cover on the right side and remove it (Figure 35).
- Pull fluorescent lamp horizontally out of the socket and remove it.
- Vacuum socket as thoroughly as possible.
- Insert new fluorescent lamp in the retainer and slide it horizontally into the socket (Figure 36).
- Install lamp cover in the housing on the left side, then press it to the right until it clicks into position (Figure 37).

Use only original Renfert replacement fluorescent (see chapter 8, Spare Parts).

7.8 Lamp cover

In the course of time, the lamp cover can become matt. This cover can be ordered as a spare part.

Changing the lamp cover:

- Switch off the unit.
- Pull plug out.
- Lightly press light cover on the right side and remove it.
- Insert the new lamp cover on the left side of the lighting case, then press on the right side until it clicks into place.

8. Spare Parts

You can find components subject to wear and the spare parts on the spare part list in the internet at www.renfert.com/p918.

Enter the following item number: 29600005 or 29601005.

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

Serial number and date of manufacturing are shown on the type plate of the unit.

9. Warranty

EN Provided the unit is properly used, Renfert **warrants** all components of the *Vario basic* for a period of **3 years**. The original sales receipt from the dealer must be provided when making warranty claims.

Parts which are subject to natural wear and tear (wear parts) and consumables are excluded from the guarantee. These parts are marked in the spare part list.

The warranty is voided in case of improper use; failure to observe the operating, cleaning, maintenance, and connection instructions; in case of independent repairs or repairs by unauthorized personnel; if spare parts from other manufacturers are employed, or; in case of unusual influences or influences not in compliance with the utilization instructions. Warranty service shall not extend the original warranty.

10. Technical Data

Mains voltage:	220 - 240 V, 50 / 60 Hz 100 - 120 V, 50 / 60 Hz
Power consumption:	16 W
Power cord length:	2 m [78 inches]
Connecting pressure, external:	6 - 8 bar [87 - 116 psi]
Operating pressure:	1 - 6 bar [14.5 - 87 psi]
Air consumption:	200 l/min. at 6 bar [7.06 cfm at 87 psi]
Ø Compressed air hose:	
internal:	4 mm [157 mil]
external:	6 mm [236 mil]
Ø Size of pipe union for external extraction unit:	
internal:	35 mm [1.38 inch]
external:	40 mm [1.57 inch]
Dimensions (W x H x D):	380 x 540 x 450 mm [15.0 x 21.3 x 17.7 inches]
Blasting chamber volume:	23 l [6.1 US gal]
Weight (empty), approx.:	14 kg [31 lbs]
This product contains a light source of energy efficiency class	F

11. Notes Regarding Return Shipment

- Request Renfert shipping boxes (Tel. +49 (0)7731 8208-383).
- Empty all blasting tanks.

! If a Renfert shipping box is not used, always remove the additional tanks (otherwise we can assume no liability for damage during shipping).

The proper return shipment form can be downloaded from the Renfert Website, www.renfert.com, under "Consulting / Repair Service".

12. Disposing of the Unit

The unit must be disposed of by an authorized recycling operation. The selected firm must be informed of all possibly health-hazardous residues in the unit.

12.1 Disposal of the Lighting Equipment

Before disposing of the sandblasting unit / the extraction box, please remove the lighting equipment including the lamp.

- Remove the mains plug.
- Loosen the 2 fixation screws from the lighting equipment at the back of the sandblasting unit / extraction box.
- Remove the lighting equipment including the lamp.

For disposal of the lighting equipment please observe the following advice:

To conserve and protect the environment, prevent environmental pollution and improve the recycling of raw materials, the European Commission adopted a directive that requires the manufacturer to accept the return of electrical and electronic units for proper disposal or recycling.

Within the European Union units with this symbol should not therefore be disposed of in unsorted domestic waste.

For more information regarding proper disposal please apply at your local authorities.

13. Standard Delivery

- 1 *Vario basic*
- 1 Foot switch
- 1 Connector kit
- 1 Accessories kit
- 1 Pair of sandblasting gloves
- 1 Operating instructions, with attachments
- 1 Power cord
- 1 Magnetic foil

14. Delivery Versions

- No. 2960-0005 *Vario basic*, basic unit, 230 V, 50 / 60 Hz, without sandblasting tanks
- No. 2960-1005 *Vario basic*, basic unit, 100 - 120 V, 50 / 60 Hz without sandblasting tanks

15. Accessories

Micro sandblasting tanks

No. 2954-0050	25 - 70 µm, Right
No. 2954-0051	25 - 70 µm, Left
No. 2954-0250	70 - 250 µm, Right
No. 2954-0251	70 - 250 µm, Left

Wand mounting

No. 2960-0001	Wall mounting, incl. 90° extractor bracket
---------------	--

Extractor bracket

No. 2960-0002	90° extractor bracket
---------------	-----------------------

Protective grate

No. 2960-0003	Protective grate for the viewing pane
---------------	---------------------------------------

Blasting nozzle

No. 90003-5520	Borocarbide nozzle
----------------	--------------------

Abrasive

Cobra Aluminium oxide

No. 1594-1105	25 µm [500 mesh], white 5 kg canister
No. 1594-1205	50 µm [270 mesh], white 5 kg canister
No. 1594-2220	50 µm [270 mesh], white 20 kg bucket
No. 1584-1005	90 µm [170 mesh], white 5 kg canister
No. 1583-1005	110 µm [150 mesh], white 5 kg canister
No. 1583-1020	110 µm [150 mesh], white 20 kg bucket
No. 1587-1005	125 µm [115 mesh], pink 5 kg canister
No. 1587-1020	125 µm [115 mesh], pink 20 kg bucket
No. 1585-1005	250 µm [60 mesh], white 5 kg canister
No. 1585-1020	250 µm [60 mesh], white 20 kg bucket

Rolloblast glass beads

No. 1594-1305	50 µm [400 - 200 mesh], 5 kg canister
No. 1594-2312	50 µm [400 - 200 mesh], 12.5 kg canister
No. 1589-1005	100 µm [170 - 100 mesh], 5 kg canister

16. Troubleshooting

Problem	Possible cause	Corrective action
No air and no abrasive	<ul style="list-style-type: none"> Nozzle blocked. Nozzle blocked at grip and/or blasting hose blocked. Kink(s) in the foot switch connector hose. Switch set to intermediate position. Pressure set too low. 	<ul style="list-style-type: none"> Open the blasting nozzle unit and clean the nozzle (Figure 23, 24, 25). Unscrew the nozzle from the handle and blow it out with compressed air. Without reconnecting the nozzle, activate the foot switch to blow compressed air through the hose. Inspect and correct the connecting hoses. Turn the switch until you feel it click into place. Open the pressure regulator further (recommended pressure: 4 - 5 bar).
Air only and no or very little abrasive	<ul style="list-style-type: none"> Dirty or moist abrasive. Insufficient abrasive in the unit. Intake hose blocked. Problems with the additional tank. 	<ul style="list-style-type: none"> Replace the abrasive. Add abrasive. Check that the intake hose is not kinked / pinched / blocked. The intake hose must be securely attached to the unit base (Figure 7). Refer to "Additional tank - Mounting instructions", No. 21-9774.

Problem	Possible cause	Corrective action
Reduced blasting force	<ul style="list-style-type: none"> • Function selector switch not properly set. • Nozzle worn. • Abrasive used up. • Abrasive feed hose is blocked / damaged / kinked. • Insufficient abrasive in the reservoir. 	<ul style="list-style-type: none"> • Make sure the selector switch clicks into place. • Replace the nozzle. • Replace the abrasive. • Drain the abrasive out and inspect / replace the feed hose. • Add abrasive.
Excessive abrasive consumption	<ul style="list-style-type: none"> • Extraction power set too high. 	<ul style="list-style-type: none"> • Reduce the extraction power.
Blasting chamber fills with dust.	<ul style="list-style-type: none"> • Insufficient extraction power. • Extractor dust bag full. 	<ul style="list-style-type: none"> • The suction performance of a weak extraction unit can be increased by closing the boreholes on the back side of the unit (Fig. 14, Pos. b) with the supplied magnetic foil. • Clean the extractor / replace the dust bag.
Abrasive spurts in pulses.	<ul style="list-style-type: none"> • Nozzle connectors not vertical. • Abrasive used up or too dirty. • Kink in the suction hose. 	<ul style="list-style-type: none"> • Make sure the nozzle is vertical. • Replace the abrasive. • Inspect the feed hose.
Unit continuously blasting and stops when foot switch is activated.	<ul style="list-style-type: none"> • Foot switch connected incorrectly. 	<ul style="list-style-type: none"> • Inspect and correct the foot switch connection.
Foot switch fails to switch on.	<ul style="list-style-type: none"> • No pressure at the compressed air connection. • Kink(s) in the pneumatic hoses. 	<ul style="list-style-type: none"> • Inspect the compressed air lines and compressor. • Realign the hoses.
Air leaking from the foot switch.	<ul style="list-style-type: none"> • Pneumatic hoses connected incorrectly. 	<ul style="list-style-type: none"> • Replace the connectors.
Lamp fails to go on.	<ul style="list-style-type: none"> • Power cord not connected. • Illumination defect. <p>The device with fluorescent tubes</p> <ul style="list-style-type: none"> • The fluorescent tube is not situated correctly in the socket. • The fluorescent tube is defect. 	<ul style="list-style-type: none"> • Plug the power cord into the wall outlet. • Change light using the service or see spare part list. • Check that the fluorescent tube is fitted correctly (see chapter 7.7.2). • Changing the fluorescent tubes (see chapter 7.7.2).

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declare under our sole responsibility, that the product(s)

Name	Part No.
Vario basic	2960 0005
Basic quattro IS	2959 0000 / 2959 0000 03
Basic quattro	2958 0000 / 2958 0000 03
Basic master	2948 2000 / 2948 2000 03 2948 2025 / 2948 2025 03 2948 2250 / 2948 2250 03
Basic classic	2947 1050 / 2947 1050 03 2947 1250 / 2947 1250 03 2947 2000 2947 2025 / 2947 2025 03 2947 2250
Basic eco	2949 1050 / 2949 1050 03 2949 1250 / 2949 1250 03 2949 2025 / 2949 2025 03
Basic mobil	2914 3050 2914 3250
Dustex master plus	2626 0105 / 2626 0105 03

is (are) in conformity with the relevant regulatory requirements by compliance with the UK designated standards.

UK legislation

Electromagnetic Compatibility Regulations 2016

Electrical Equipment (Safety) Regulations 2016

The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

Designated standards

EN 61326-1:2013

EN 61010-1:2010/A1:2019/AC:2019-04

EN IEC 63000:2018

Tilo Burgbacher,
Engineering Director

Hilzingen, 02.11.2022

The Technical Documentation has been retained by Renfert GmbH

**Hochaktuell und ausführlich auf ...
Up to date and in detail at ...
Actualisé et détaillé sous ...
Aggiornato e dettagliato su ...
La máxima actualidad y detalle en ...
Актуально и подробно на ...**

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Renfert GmbH, Industriegebiet, 78247 Hilzingen / Germany
Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany