making work easy

- de
- en
- fr
- it
- es
- pt
- tr
- ru
- pl
- zh
- ja
- ko
- ar

SILENT flow sensor

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Made in Germany

21-6023 22082024

making work easy

SILENT flow sensor

Inhalt

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Contents

en

1	Introduction	3
1.1	Symbols used.....	3
2	Safety	3
2.1	Intended use.....	3
2.2	Improper use.....	4
2.3	Ambient conditions for safe operation.....	4
2.4	Ambient conditions for storage and transport.....	4
2.5	Hazard and warning information.....	4
2.5.1	General information.....	4
2.5.2	Specific information.....	5
2.6	Authorized persons.....	5
2.7	Disclaimer.....	5
3	Product description	6
3.1	General description.....	6
3.2	Assemblies and functional elements.....	6
3.3	Scope of delivery.....	6
4	Startup	7
4.1	Unpacking.....	7
4.2	Connection.....	7
4.3	Electrical connection.....	8
4.4	Setting the warning threshold (minimum flow rate).....	8
4.5	Zero Adjustment (Taring).....	8
4.6	Checking leak-tightness.....	9
4.7	Checking function.....	9
5	Operation	9
5.1	Operation.....	9
5.2	Checking function.....	10
6	Cleaning / maintenance	10
6.1	Cleaning.....	10
6.2	Cleaning the suction tunnel.....	11
6.3	Effectiveness of the measurement equipment.....	11
7	Troubleshooting	11
8	Technical data	12
9	Warranty	12
10	Disposal information	12
10.1	Disposing of consumables.....	12
10.2	Disposing of the equipment.....	12
10.2.1	Disposal information for countries in the EU.....	12

1 Introduction

The SILENT flow sensor is used to monitor the flow rate at a suction point for workbenches where dust is suctioned. For dust detection units that are tested and recognized in accordance with GS-IFA-M20, a minimum flow rate is defined that must be reached for correct operation.

The minimum flow rates currently used (15 l/s, 20 l/s, 25 l/s) can be defined on the SILENT flow sensor. 25 l/s must be defined in the case of a suction point that is not certified.

If the minimum flow rate that is defined is not reached, the SILENT flow sensor outputs a warning signal. The SILENT flow sensor is not suitable for use with, e.g., sandblasters!

1.1 Symbols used

Symbols with the following meanings are used in these instructions or on the equipment:

Danger
Immediate risk of injury. Consult accompanying documents!

Electrical voltage
Danger due to electrical voltage.

Caution
Risk of damage to the equipment if the instruction is not followed.

Notice
Provides the operator with useful information that eases use.

Protection class II

The equipment complies with applicable EU directives.

This product complies with the relevant UK legislation.
See the UKCA declaration of conformity available online at www.renfert.com.

The equipment is subject to the provisions of the EU directive 2002/96/EC (WEEE directive).

► List, particular attention should be paid

- List
- Subordinate list

⇒ Instructions / appropriate action / input / operational sequence:

You are asked to carry out the specified action in the order indicated.

- ◆ Result of an action / response by the equipment / response by the program:

The equipment or program responds as a result of your actions or when a specific incident occurs.

Other symbols are explained as they occur.

2 Safety

2.1 Intended use

The area of application of the SILENT flow sensor is with extraction systems that are installed at dental technician workbenches where restorations are processed using dental handpieces.

The SILENT flow sensor can be used together with one of the following extraction units in the SILENT range:

- SILENT compact
- SILENT TS / TS2
- SILENT TC / TC2 / EC2

within the scope of measures to protect against exposure to dust *) as defined by the German ordinance on hazardous substances ("GefStoffV") / respective employers' liability insurance association.

In doing so, both the instructions for use of the SILENT flow sensor as well as those of the extraction unit must be observed, particularly with regard to low-dust disposal of the extracted material and filter replacement.

*) Note for customers in Germany: A validated and recognized dust protection measure in accordance with testing principle GS-IFA-M20 (issued by the German Institute for Work and Health IFA) shall only be deemed to have been provided when all components, i.e. the flow rate monitoring unit such as the SILENT flow sensor, the extraction unit, and the detection unit, are validated and recognized. In the case of multi-workbench extraction units and the use of Y-junctions or extractor switches, every suction channel / every suction point must be fitted with a flow rate monitoring unit.

2.2 Improper use

The SILENT flow sensor may not be operated in combination with equipment that generates a significant amount of dust such as sandblasting units, trimmers, or CAM milling systems!

Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the equipment, increase the risk of serious injury and lead to damage to the environment or the product itself.

2.3 Ambient conditions for safe operation

The unit may only be operated:

- Indoors
- Up to an altitude of 2,000 m above sea level
- At an ambient temperature of between 5 - 40 °C [41 - 104 °F] *)
- At a maximum relative humidity of 80 % at 31 °C [87.8 °F], decreasing linearly to 50% relative humidity at 40 °C [104 °F] *)
- With mains power where the voltage fluctuations do not exceed 10 % of the nominal value
- Under contamination level 2 conditions
- Under over-voltage category II conditions

*) At between 5 - 30 °C [41 - 86 °F], the unit can be operated at a humidity of up to 80 %. At temperatures between 31 - 40 °C [87.8 - 104 °F], the humidity must decrease proportionally in order to ensure operational readiness (e.g., at 35 °C [95 °F] = 65 % humidity, at 40 °C [104 °F] = 50 % humidity). The unit may not be operated at temperatures above 40 °C [104 °F].

2.4 Ambient conditions for storage and transport

The following ambient conditions must be observed during storage and transport:

- Ambient temperature of -20 - +60 °C [-4 - +140 °F]
- Maximum relative humidity 80 %

2.5 Hazard and warning information

2.5.1 General information

- ▶ If the unit is not operated in accordance with the supplied instructions for use, the intended safety can no longer be guaranteed.
- ▶ The unit may only be operated using a power cord equipped with a plug appropriate to the local power supply. Any necessary modifications may only be carried out by a qualified electrician.
- ▶ Only operate the unit if the information on the nameplate corresponds to the specifications of your local mains power supply.
- ▶ The unit may only be plugged into electrical sockets that are connected to the protective conductor system.
- ▶ The power plug must be easily accessible.
- ▶ Disconnect the unit from the mains power supply before carrying out work on the electrical parts.
- ▶ It is the responsibility of the operator to ensure that national regulations are complied with during operation and regarding recurring safety inspection of electrical equipment. In Germany, these are DGUV Regulation 3 in conjunction with VDE 0701-0702.
- ▶ Check connection cables (such as power cords), tubes and housing (such as the button pad) regularly for damage (such as kinks, cracks and porosity) or signs of ageing.
Do not operate equipment with damaged connection cables, hoses or housing parts or other defects!
- ▶ Defective equipment must be taken out of service immediately. Disconnect the power plug and secure the equipment against further use. Send the equipment in for repair!
- ▶ Only operate the equipment under supervision.
- ▶ Please observe the national accident prevention regulations!
- ▶ Information on REACH and SVHC is available on our website at www.renfert.com, in the Support area.

2.5.2 Specific information

- ▶ Only the adapter provided on delivery may be used for adaption to the local power supply.
- ▶ The German ordinance on hazardous substances ("GefStoffV") or equivalent national regulations must be observed during operation, cleaning, and maintenance.

2.6 Authorized persons

Operation and maintenance of the equipment may only be performed by trained personnel.

Any repairs not described in this user information may only be performed by a qualified electrician.

2.7 Disclaimer

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- ▶ The product is employed for any purposes other than those specified in the instructions for use.
- ▶ The product is altered in any way other than those alterations described in the instructions for use.
- ▶ The product is not repaired by an authorized dealer or if original Renfert spare parts are not used.
- ▶ The product continues to be used despite obvious safety defects or damage.
- ▶ The product is subjected to mechanical impact or is dropped.

3 Product description

en

3.1 General description

The SILENT flow sensor measures the flow rate that flows through it. The measured flow rate is compared electronically with a threshold value defined in the unit. If the flow rate is less than the defined threshold value for a certain period, this is indicated by an acoustic as well as a visual signal.

The SILENT flow sensor meets the requirements of the section on flow rate monitoring units in Annex A of the testing principles GS-IFA-M20-2012-12 issued by the IFA (German Institute for Work and Health).

3.2 Assemblies and functional elements

- | | | | |
|---|---|---|---|
| 1 | Inlet nozzle (IN) | 6 | Indicator for threshold value 20 l/s (yellow) |
| 2 | Indicators and operation | 7 | Indicator for threshold value 25 l/s (yellow) |
| 3 | Outlet nozzle (OUT) | 8 | Button |
| 4 | End bushing | 9 | Operating indicator (green) |
| 5 | Indicator for threshold value 15 l/s (yellow) | | |


Fig. 1

3.3 Scope of delivery

- 1 SILENT flow sensor
- 2 End bushing
- 1 Universal power plug set (only for item number 2914 1000)
- 4 Cable ties
- 1 Quick Start Guide

4 Startup

4.1 Unpacking

- ⇒ Remove the unit and accessories from the shipping box.
- ⇒ Check the delivery for completeness (compare with the "Scope of delivery").

4.2 Connection

⇒ Determine a location where you can install the SILENT flow sensor.
 The position where the SILENT flow sensor is installed between the extraction unit and the suction point is not important.

- ! **When using extractor switches, the SILENT flow sensor must always be placed between the switch and the suction mouth.**
- ! **The SILENT flow sensor must be operated in an upright position. The suction current must flow through the SILENT flow sensor from top to bottom. If installed on its side, dust could accumulate and distort the measurement result.**
- ! **When routing the tube, avoid sharp bends and kinks.**
- ! **Before installation, verify that the inlet nozzle IN (1) and the outlet nozzle OUT (3) are screwed in fully until stop into the SILENT flow sensor.**
- i* **If the inlet nozzle and outlet nozzle are not fully screwed in, this could result in leakage that could distort the measurement result and cause unwanted noise.**


Fig. 2

- ⇒ Cut the suction hose leading to the relevant suction point at a suitable location, e.g., using wire cutters.
- ⇒ Screw the two end bushings provided (4) on to the newly created tube ends. Direction of rotation: Counter-clockwise!
- ⇒ By rotating, insert the end bushing (4) of the tube leading from the suction point into the inlet nozzle IN (1) until stop.
- ⇒ By rotating, insert the end bushing (4) of the tube leading to the extraction unit into the outlet nozzle OUT (3) until stop.
- ⇒ Install the SILENT flow sensor using appropriate attachment material, for example the cable ties provided or screws (screws are not included with the delivery).

Typical assembly


Fig. 3

4.3 Electrical connection

If required, use one of the mains power adapters provided for connection to a building electrical socket.

⇒ Unroll the power cord and insert the power plug into a building electrical socket.

- ◆ All 4 LEDs light up, the measurement system initializes. This can take approx. 10 seconds.
- ◆ A short acoustic signal is output when initialization is complete.

Once initialization is complete:

- ◆ The green operating indicator (9) lights up.
- ◆ The yellow indicator for the defined threshold value lights up (5, 6, 7).


Fig. 4

4.4 Setting the warning threshold (minimum flow rate)

The threshold value to be set depends on the dust detection unit (suction mouth) used.

In the case of dust detection units tested in accordance with GS-IFA-M20, the minimum flow rate is indicated on the dust detection unit.

For the Renfert suction mouth (REF: 2925 0000), a threshold value of "20" must be defined. For threshold values for other certified dust detection units, contact the respective manufacturer.

If a threshold value cannot be obtained, the largest threshold value of "25" must be set accordingly on the SILENT flow sensor and an appropriately adjusted, higher suction level selected on the extraction unit.

⇒ Press the button until a confirmation tone is output.

- ◆ The green operating indicator (9, Fig. 4) flashes quickly.
- ◆ The yellow indicator for the defined threshold value lights up (5, 6, 7, Fig. 4).

To modify the threshold value:

⇒ Press the button briefly.

- ◆ Every time you press the button, a confirmation tone is output, and the next threshold value is selected.

To save the selected threshold value:

⇒ Press the button until a confirmation tone is output.

- ◆ The green operating indicator (9, Fig. 4) lights up.
- ◆ The yellow indicator for the defined threshold value lights up (5, 6, 7, Fig. 4).

If a button is not pressed for 5 seconds, the threshold value defined at that time is automatically saved; a short warning tone is output and the unit exits setting mode.

4.5 Zero Adjustment (Taring)

To ensure a correct measurement, a zero adjustment (taring) must be performed on the SILENT flow sensor during commissioning and at regular intervals.

The zero adjustment is performed if no flow rate passes through the SILENT flow sensor.

⇒ The SILENT flow sensor is installed and electrically connected.

⇒ Extraction is not in operation!

⇒ Press the button for approx. 5 sec..

- ◆ After approx. 2 sec., a single acoustic signal sounds for setting the warning threshold.

⇒ Continue to hold down the button.

- ◆ After approx. 5 sec., a double acoustic signal sounds to perform the zero adjustment.

⇒ Release the button.

- ◆ The displays flash alternately.
- ◆ The zero adjustment is performed.
- ◆ After approx. 5 sec. the SILENT flow sensor returns automatically into operating mode; the zero adjustment is completed.

4.6 Checking leak-tightness

- ◆ The SILENT flow sensor is connected to a socket.
 - ◆ The green operating indicator (9, Fig. 4) lights up.
 - ◆ The yellow indicator for the defined threshold value lights up (5, 6, 7, Fig. 4).
- ⇒ Switch suction on.
- ◆ The green operating indicator begins to flash.
 - ◆ The yellow indicator for the defined threshold value lights up.
- ⇒ Close the suction mouth with your hand so that air is no longer drawn in through the suction mouth.
- ⇒ Pay attention to noises that could indicate leakage.

Green indicator (9, Fig. 4)	Yellow indicator (5, 6, 7, Fig. 4)	Result / measure
Light is on	Light is on	Leak-tight / no leakage
Flashing	Light is on or flashing	Leak between the suction mouth and the SILENT flow sensor / check the suction route between the suction mouth and the SILENT flow sensor, in particular whether the inlet nozzle IN (1) has been fully and securely screwed into the sensor.

4.7 Checking function

After assembly, check the leak-tightness of the suction route and the function of the SILENT flow sensor.

- ◆ The SILENT flow sensor is connected to a socket.
 - ◆ The green operating indicator (9, Fig. 4) lights up.
 - ◆ The yellow indicator for the defined threshold value lights up (5, 6, 7, Fig. 4).
- ⇒ Switch suction on.
- ◆ The green operating indicator begins to flash.
 - ◆ The yellow indicator for the defined threshold value lights up.
- ⇒ Close part of the suction mouth with your hand so that the flow rate drops below the defined threshold value.
- ◆ After a short time, the indicator for the corresponding threshold value starts to flash.
 - ◆ The warning tone is output after approx. 30 seconds.
- ⇒ Pay attention to noises that could indicate leakage.
- ⇒ Release the suction mouth again and continue to suction for approx. 30 seconds.
- ◆ The indicator for the corresponding threshold value lights up steadily.
 - ◆ The warning tone is no longer output.
- ⇒ Switch off the extraction unit.

5 Operation

5.1 Operation

If the SILENT flow sensor detects a flow rate at the suction point, the green operating indicator flashes slowly (9, Fig. 4).

As long as the flow rate is higher than the defined threshold value, the corresponding indicator remains lit up (5, 6, 7, Fig. 4)

If the flow rate is less than the defined threshold value for a certain period, the corresponding indicator begins to flash and a warning tone is output every 30 seconds.

If the flow rate is too low at the suction point being monitored, safe operation is no longer guaranteed and the level of exposure to dust can be harmful to health!

In this case, the following measures must be taken:

- ⇒ Stop working at all suction points of the affected extraction unit. This applies to the individual suction channels of any multi-workbench extraction unit as well as to all suction points operated via a Y-junction or via extractor switches.
- ⇒ Determine why the flow rate is too low and resolve the issue.

Cause	Measure
Suction level is too low.	<ul style="list-style-type: none"> Choose a higher suction level
In the case of multi-workbench extraction units: the suction power is too low because an additional suction channel has been connected.	<ul style="list-style-type: none"> Stop working with the additional suction channel.
Filter blockage	<ul style="list-style-type: none"> In the case of extraction units with a pre-filter (filter bag): change the filter bag. When doing so, observe the instructions and the safety information for the corresponding extraction unit. In the case of extraction units without a filter bag: clean the filter.
Changing the filter bag / cleaning the filter do not result in the expected improvement.	<ul style="list-style-type: none"> Change the fine filter. When doing so, observe the instructions and the safety information for the corresponding extraction unit in order to ensure that changing the fine filter results in as little dust as possible.
Blockage in the suction point (suction mouth).	<ul style="list-style-type: none"> Disconnect the suction hose at the suction point and check whether there is a blockage in the suction point; remove where appropriate.
Blockage in the suction line.	<ul style="list-style-type: none"> Check the suction line at all disconnection points / service openings section by section for blockages.
Leakage in the section of the suction line between the SILENT flow sensor and the extraction unit.	<ul style="list-style-type: none"> Check that the outlet nozzle of the SILENT flow sensor has been fully screwed into the SILENT flow sensor until stop and that the end bushing of the suction tube is fully inserted into the outlet nozzle.

Determining the cause and measures can be carried out in the order shown.

After every measure, you should switch on the extraction unit again and check whether the minimum flow rate is again exceeded so that correct operation is possible again.

⇒ Work should only resume after the fault has been rectified and the minimum flow rate is exceeded again. It is also the case with multi-workbench extraction units that work must be interrupted at all connected suction points until the fault is rectified.

5.2 Checking function

Every day when beginning your work, check the function of the SILENT flow sensor as described in Section 4.6 and 4.7.

6 Cleaning / maintenance

The unit does not contain any internal parts that require maintenance.

6.1 Cleaning

Disconnect the power plug!

To clean the unit, only wipe the outside with a damp cloth.

Do not use any solvent-based or abrasive cleaning agents.

6.2 Cleaning the suction tunnel

The inside of the unit must be cleaned once a year.

Disconnect the power plug!

Do not use compressed air!

Compressed air irreparably damages the sensitive measurement equipment. Using compressed air to blow off deposits of dust is generally not permitted.

- ⇒ Withdraw the end bushing from the inlet and outlet nozzle.
- ⇒ Unscrew the inlet and outlet nozzle and wipe out using a dry cloth.
- ⇒ Clean the suction tunnel using a dry cloth or a bottle cleaning brush.
- ⇒ Screw the inlet and outlet nozzle into the SILENT flow sensor again fully until stop.
- ⇒ By rotating, insert the end bushing of the ends of the suction hose back into the inlet and outlet nozzle until stop.

6.3 Effectiveness of the measurement equipment

Check the effectiveness of the measurement equipment twice a year.

To do so, perform the steps in the following segments:

- 4.5 Zero Adjustment (Taring)
- 4.6 Checking leak-tightness
- 4.7 Checking function

7 Troubleshooting

Errors	Cause	Corrective action
The warning tone is output.	<ul style="list-style-type: none"> • The threshold value defined for the flow rate so that work can be carried out safely is not reached. 	<ul style="list-style-type: none"> • Choose a higher suction level. • Stop working. • Check the suction route, remove blockages. • Check the extraction unit, clean the filter. See also Chapter 5.1
The warning tone is output although the flow rate at full suction power is greater than the defined warning threshold.	<ul style="list-style-type: none"> • The flow through the SILENT flow sensor is in the wrong direction. 	<ul style="list-style-type: none"> • Observe the labels on the connection pieces. • Connect the suction hose leading from the suction point to the inlet nozzle (IN). • Connect the suction hose leading to the extraction unit to the outlet nozzle (OUT).
The green indicator light (9, Fig. 1) shows a volume flow (flashing) even though the extraction is not switched on.	<ul style="list-style-type: none"> • The zero adjustment has changed. 	<ul style="list-style-type: none"> • Perform a zero adjustment, see Chapter 4.5.
The acoustic warning sounds even though the volume flow is evidently sufficient.	<ul style="list-style-type: none"> • The zero adjustment has changed. 	<ul style="list-style-type: none"> • Perform a zero adjustment, see Chapter 4.5.
A short acoustic signal sounds every second during zero adjustment.	<ul style="list-style-type: none"> • The determined value for the zero adjustment is too high indicating a defective sensor. 	<ul style="list-style-type: none"> • Switch off the acoustic signal by briefly pressing the button. • Send in the SILENT flow sensor.

8 Technical data

en

Item number:	2914 0000 / 2914 1000
Nominal voltage / V:	230
Permissible mains voltage / V:	100 - 240
Mains frequency / Hz:	50 / 60
Power consumption / W:	0,6
Dimensions (width x height x depth) / mm [inch]:	108 x 223 x 94 [4.25 x 8.78 x 3.70]
Weight / kg [lbs.]:	0,8 [1.76]

9 Warranty

Renfert offers a **3-year warranty** on all parts of the equipment, provided that the equipment is properly used. Warranty claims may only be made upon presentation of the original sales receipt from the authorized dealer.

Parts that are subject to natural wear and tear (wear parts) and consumables are excluded from the warranty. These parts are marked in the spare parts list.

The warranty is voided in the case of improper use; failure to observe the operating, cleaning, maintenance, and connection instructions; in the case of independent repairs or repairs that are not carried out by an authorized dealer; if spare parts from other manufacturers are used, or in the case of unusual influences or influences not in compliance with the utilization instructions.

Warranty claims shall not extend the warranty period.

10 Disposal information

10.1 Disposing of consumables

10.2 Disposing of the equipment

The equipment must be disposed of by a specialist facility. The specialist facility must be informed of any residue in the equipment that is harmful to health.

10.2.1 Disposal information for countries in the EU

To conserve and protect the environment, prevent environmental pollution, and improve the recycling of raw materials, the European Commission has adopted a directive that requires the manufacturer to accept the return of electrical and electronic equipment for proper disposal or recycling.

For this reason, equipment with this symbol may not be disposed of in unsorted domestic waste within the European Union.

Please contact your local authorities for more information on proper disposal.

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declare under our sole responsibility, that the product(s)

Name	Part No.
SILENT flow sensor 100-240 V 50/60 Hz Steckertyp 1	29140000
SILENT flow sensor 100-240 V 50/60 Hz Universaladapter	29141000

is (are) in conformity with the relevant regulatory requirements by compliance with the UK designated standards.

UK legislation

Electromagnetic Compatibility Regulations 2016

Electrical Equipment (Safety) Regulations 2016

The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

Designated standards

EN 61010-1: 2010 + A1:2019 + A1:2019/ AC:2019

EN 61326-1: 2013

EN IEC 63000:2018

Tilo Burgbacher,
Engineering Director

Hilzingen, 09.03.2023

The Technical Documentation has been retained by Renfert GmbH

**Hochaktuell und ausführlich auf ...
Up to date and in detail at ...
Actualisé et détaillé sous ...
Aggiornato e dettagliato su ...
La máxima actualidad y detalle en ...
Актуально и подробно на ...**

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Renfert GmbH, Industriegebiet, 78247 Hilzingen / Germany
Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany