making work easy

# SILENT EC2

FAQ - frequently asked questions

Made in Germany

FAQ - Deutsch

FAQ - English

FAQ - Français

FAQ - Italiano

FAQ - Español

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# FAQ - English

For information on operating elements, see [Fig. 1](#) and [Fig. 2](#) on page 2 of this document.

For chapter indications, refer to the complete instruction manual on the internet.

EN

Trouble	Reason	Solution
<b>When the device is switched on, there is a loud vibration noise for approx. 8 seconds.</b>	<ul style="list-style-type: none"> <li>The suction unit performs an automatic filter cleaning sequence.</li> </ul>	<ul style="list-style-type: none"> <li>This is functional and requires no intervention.</li> </ul>
<b>An acoustic signal is emitted, the suction unit is switched off and a filter cleaning is performed.</b>	<ul style="list-style-type: none"> <li>The internal flow velocity limit was not reached.</li> <li>After 8 hours of operation (turbine air time) without switching off, a filter cleaning is carried out.</li> </ul>	<ul style="list-style-type: none"> <li>Repeat work after cleaning is finished.</li> <li>Switch off devices at the end of each working day at the on / off switch <a href="#">(3, Fig. 1)</a> ausschalten.</li> </ul>
<b>After switching on, the display shows "Dust drawer" <a href="#">(23, Fig. 2)</a> and a 3 x acoustic signal is heard.</b>	<ul style="list-style-type: none"> <li>The time interval for emptying the dust drawer has expired and the dust drawer has not yet been emptied.</li> <li>The emptying of the dust drawer has not yet been confirmed.</li> </ul>	<ul style="list-style-type: none"> <li>Empty dust drawer and confirm by pressing the Enter key <a href="#">(30, Fig. 2)</a>.</li> <li>Confirm emptying the dust drawer by pressing the Enter key <a href="#">(30, Fig. 2)</a>.</li> </ul>
<b>"Cal" is shown in the display.</b>	<ul style="list-style-type: none"> <li>The automatic switch-on function cannot be set. At the selected rotational speed on the handpiece, the difference between stand-by current and operating current is too low.</li> <li>The automatic switch-on function cannot be set.</li> </ul>	<ul style="list-style-type: none"> <li>Adjust the switch-on threshold to a higher speed on the handpiece.</li> <li>Operate the device manually (see chapter 5.2).</li> </ul>
<b>"Air" is shown in the display.</b>	<ul style="list-style-type: none"> <li>During self-diagnosis a lack of compressed air was detected.</li> </ul>	<ul style="list-style-type: none"> <li>Check the air pressure connection. When compressed air is available again, clear error message pressing the menu key (20).</li> </ul>
<b>"HOT" is shown in the display.</b>	<ul style="list-style-type: none"> <li>The electronics have become overheated.</li> </ul>	<ul style="list-style-type: none"> <li>Switch off the device and let it cool down.</li> <li>Provide adequate cooling, e.g. <ul style="list-style-type: none"> <li>- Chapter 4.2.</li> <li>- Use external exhaust air duct (see chapter 4.7).</li> <li>- Change the fine filter (see accessories as well as chapter 6.3).</li> </ul> </li> </ul>
<b>"ELEC" is shown in the display.</b>	<ul style="list-style-type: none"> <li>Leakage: Dust drawer not airtight.</li> <li>Leakage: Second channel open although it should be closed.</li> <li>Turbine does not run even though electronics have been switched on.</li> <li>Faults in electronics.</li> </ul>	<ul style="list-style-type: none"> <li>Check the correct fit of the dust drawer (see chapter 6.2).</li> <li>Check compressed air pinch valves.</li> <li>Clean the air filter (see chapter 6.5).</li> <li>Check that the pinch valve in the unused channel is closed.</li> <li>Check whether the turbine cable is properly plugged in (especially after changing the turbine).</li> <li>Contact Renfert / Service.</li> </ul>
<b>"FILT" is shown in the display.</b>	<ul style="list-style-type: none"> <li>The fine filter is so heavily clogged that the filter cleaning does not show sufficient effect.</li> </ul>	<ul style="list-style-type: none"> <li>Change the fine filter (see accessories as well as chapter 6.3).</li> </ul>
<b>Dust drawer over full.</b>	<ul style="list-style-type: none"> <li>Time interval to the "FULL" display of the dust drawer too high.</li> </ul>	<ul style="list-style-type: none"> <li>Set a shorter time interval (see chapter 5.4.2).</li> </ul>
<b>The signal for emptying the dust drawer appears, even though it is not yet full.</b>	<ul style="list-style-type: none"> <li>The set time interval to the "FULL" dust drawer is too low.</li> </ul>	<ul style="list-style-type: none"> <li>Set a longer time interval (see chapter 5.4.2).</li> </ul>

Trouble	Reason	Solution
<b>The suction performance is insufficient.</b>	<ul style="list-style-type: none"> <li>• Suction performance is set too low.</li> <li>• Blockage in the Suction tube.</li> <li>• High loss of suction due to unfavorable tube position.</li> <li>• Leakage: Dust drawer not airtight.</li> <li>• Leakage: Second channel open although it should be closed.</li> <li>• Fine filter blocked</li> </ul>	<ul style="list-style-type: none"> <li>• Select a higher suction level.</li> <li>• Check suction tube.</li> <li>• Read instructions on connection to the extraction point (see chapter 4.5).</li> <li>• Check the correct fit of the dust drawer (see chapter 6.2).</li> <li>• Check pinch valves in compressed air.</li> <li>• Check that the pinch valve in the unused channel is closed.</li> <li>• Switch the device off and on again, so that a filter cleaning is carried out.</li> <li>• Change fine filter (see chapter 6.3) (if filter cleaning does not improve the performance).</li> </ul>
<b>Electrical equipment connected to the appliance socket cannot be operated.</b>	<ul style="list-style-type: none"> <li>• Suction not switched on.</li> <li>• The device fuse has triggered because the power of the connected electrical device is too high.</li> </ul>	<ul style="list-style-type: none"> <li>• Switch suction on.</li> <li>• Observe the maximum connection capacity of the appliance socket.</li> <li>• Connect the electrical appliance to a separate socket and operate the extraction unit manually (see chapter 5.2).</li> </ul>
<b>Suction starts even though the electrical appliance connected to the socket is not being used.</b>	<ul style="list-style-type: none"> <li>• The automatic switch-on function has been set too low.</li> </ul>	<ul style="list-style-type: none"> <li>• Setting the automatic switch-on (see chapter 5.4.3).</li> </ul>
<b>Suction does not stop when the electrical appliance is switched off.</b>	<ul style="list-style-type: none"> <li>• The automatic switch-on function has been set too low.</li> </ul>	<ul style="list-style-type: none"> <li>• Setting the automatic switch-on (see chapter 5.4.3).</li> </ul>
<b>Suction does not start even though the electrical appliance connected to the socket is in used.</b>	<ul style="list-style-type: none"> <li>• The automatic switch-on function has been set too high.</li> </ul>	<ul style="list-style-type: none"> <li>• Setting the automatic switch-on (see chapter 5.4.3).</li> </ul>
<b>Fuse triggers when a connected electrical device is switched on.</b>	<ul style="list-style-type: none"> <li>• Power consumption of the connected device too high.</li> </ul>	<ul style="list-style-type: none"> <li>• Observe the max. connection value (see chapter 8).</li> </ul>
<b>Suction stops suddenly in continuous or automatic mode.</b>	<ul style="list-style-type: none"> <li>• Suction turbine is overheated.</li> <li>• Suction turbine is defect.</li> </ul>	<ul style="list-style-type: none"> <li>• Switch off the device and leave to cool for at least 60 min.</li> <li>• Check whether suction tube is blocked, eliminate blockage.</li> <li>• Switch the device off and then on again to perform a filter cleaning.</li> <li>Change fine filter (see chapter 6.3) if filter cleaning does not improve the performance.</li> <li>• Change Suction turbine.</li> </ul>
<b>After the suction motor was changed, the suction unit stopped functioning or stops suddenly.</b>	<ul style="list-style-type: none"> <li>• The plug from the suction motor is not sufficiently clicked in place.</li> </ul>	<ul style="list-style-type: none"> <li>• Connect the suction motor plug correctly. Ensure that the plug is fitted correctly and clicked into place by pulling it.</li> </ul>

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