making work easy






- de
- en
- fr
- it
- es
- pt
- tr
- ru
- pl
- cs
- uk
- zh
- ja
- ko
- ar
- sv
- fi




# Vibrax

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Made in Germany

21-9793 15112024

making work easy

# Vibrax

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Made in Germany

# Inhalt

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# 1 Introduction

We are pleased with your decision to purchase the Renfert Vibrax vibrator. This device sets a new standard with regard to functionality, performance and ergonomics.

**Please read the following operating instructions carefully and observe the safety information they contain in order to ensure a long and trouble-free service life.**

## 1.1 Symbols used

The following symbols are employed in these instructions and on the unit itself:

### **Danger**

This indicates a direct risk of injury.

### **Electrical current**

This indicates a hazard due to electrical current

### **Attention**

Failure to observe the associated information can result in damage to the unit.

### **Note**

*This provides the operator with useful information to make working with the unit easier.*

The device complies with the requirements of the applicable EU directives.

This product complies with the relevant UK legislation.

See UKCA Declaration of Conformity in the Internet under [www.renfert.com](http://www.renfert.com).

The device is subject to the EU directive 2002/96/EG (WEEE directive).

Only intended for indoor use.

Before opening the unit, disconnect it from the mains power supply by unplugging the power cord from the wall outlet.

Tool included.

**Observe the Operating Instructions.**

Other symbols are explained as they occur.

## 2 Safety

**Using these operating instructions as a starting point, instruct all operators of the unit with regard to the area of application, the possible hazards during operation, and the proper operation of the Vibrax.**

Please have these operating instructions readily available for the operators.

### 2.1 Proper Use

This device is intended solely for the compaction of dental investments and plasters and to remove air pockets from them.

The selective introduction of vibrations improves the flow characteristics of investments and plasters and permits bubble-free pouring of castings and cast muffles.

### 2.2 Ambient Conditions (in accordance with DIN EN 61010-1)

The unit may only be operated:

- Indoors;
- Up to an altitude of 2,000 m above sea level;
- At an ambient temperature range between 5 - 40 °C [41 - 104 °F] \*);
- At a maximum relative humidity of 80 % at 31 °C [87.8 °F], dropping in a linear manner to 50 % relative humidity at 40 °C [104 °F] \*);
- With mains power where the voltage fluctuations do not exceed 10 % of the nominal value;

- Under contamination level 2 conditions;
- Under over-voltage category II conditions.

\*) Between 5 – 30 °C [41 – 86 °F], the unit can be operated at a relative humidity of up to 80 %. At temperatures between 31 – 40 °C [87.8 – 104 °F], the humidity must decrease proportionally in order to ensure operational readiness (e.g., at 35 °C [95 °F] = 65 % humidity; at 40 °C [104 °F] = 50 % humidity). The unit may not be operated at temperatures above 40 °C [104 °F].

## 2.3 Hazard and Warning Information

**Only intended for indoor use. The unit is only designed for dry applications and may not be operated or stored outdoors or under wet conditions.**

**The unit may only be operated with a power cord equipped with a plug appropriate to the local power supply.**

**The unit may only be operated if the information on the nameplate conforms with the specifications of your local mains power supply.**

**Regularly inspect connecting lines and hoses (e.g., the power cord) for damage (e.g., kinks, cracks, porosity) or signs of aging. Units exhibiting damaged connecting lines, hoses, or other defects must be taken out of service immediately.**

**Always unplug the unit from the wall outlet before beginning any work on the unit's electrical components.**

**Do not allow the device to be operated without supervision.**

**It is the responsibility of the operator that national regulations during operation and regarding a repeated safety inspection of electrical equipment are complied with. For Germany these are the regulation 3 by DGUV (German Statutory Accident Insurance) in relation with VDE 0701-0702 (Association for Electrical, Electronic and Information Technology).**

**Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the device, increase the risk of serious injury and lead to damage to the environment or the device itself.**

**Information on REACH and SVHC is available on our website [www.renfert.com](http://www.renfert.com), in the Support area.**

## 2.4 Authorised Individuals

Only properly trained individuals may operate and service the Vibrax. Any repairs not specifically described in these operating instructions may only be performed by a qualified electrician.

## 2.5 Liability Exclusion

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- **The product is employed for any purposes other than those cited in the operating instructions;**
- **The product is altered in any way other than those alterations described in the operating instructions;**
- **The product is repaired by other than an authorized facility or if any but Renfert OEM parts are employed;**
- **The product continues to be employed, despite obvious safety faults or damage;**
- **The product is subjected to mechanical impacts or is dropped.**

# 3 Product description

## 3.1 Operating Elements

Refer to Figures 1 and 2.

- (A) ON switch and oscillation mode switch Coarse/Fine
- (B) Power level switch
- (C) Rubber mat
- (D) Mains connection
- (E) Device fuse

## 3.2 Standard Delivery

- 1 Vibrax
- 1 Rubber mat
- 1 Power cord
- 1 Operating instructions
- 1 SW4 Allen key (in the control lever, Figure 3)

## 3.3 Accessories

1830-0001 Vibration ball

# 4 Setup and Commissioning

The unit may only be operated if the information on the nameplate conforms with the specifications of your local mains power supply.

## 4.1 Setup

Select a stable base suitable to bear the weight of the Vibrax. Preferably, the Vibrax should be set up on a solid, heavy slab.

Please install the unit in order to reach the mains connection easily.

## 4.2 Rubber Mat Installation

Before you can begin working with the Vibrax you must first secure the included rubber mat on the vibrator table.

**Ensure that the rubber mat locks into place with its groove below the edge of the vibration plate over its entire circumference (Figure 6).**

## 4.3 Mains Connection

The device may only be commissioned if the information on the nameplate corresponds to the specifications of your local electric power supply utility!

- Move switch (A) to its vertical (OFF) position (Figure 1).
- Plug the power cord into the mains connection (D) (Figure 2).
- Plug the other end of the power cord into the mains outlet.

The Vibrax is now ready to use.

# 5 Operation

Refer to Figure 4.

Switch (A) is used to switch the Vibrax ON and OFF.

Switching the device on simultaneously selects the vibrator activation.

Switch settings:

OFF

ON, fine oscillation mode

ON, coarse oscillation mode

Using the 4 power level switch settings in combination with the 2 selectable wave forms, you can adjust the vibration intensity to optimally match the materials being processed.

You can change the setting of the power level switch (B) at any time, without having to first switch the Vibrax off.

### Tip:

Until you become familiar with the Vibrax and its operation, we recommend turning the power level switch to its lowest setting before switching the Vibrax on.

## 6 Cleaning / Maintenance

The Vibrax is maintenance-free.

EN

### 6.1 Cleaning the Rubber Mat

If and when it becomes dirty, the rubber mat can be easily removed for cleaning:

- Remove the mat (Figure 5).
- Clean the mat with water and then wipe it dry.
- If necessary, apply a thin coat of silicone to the mat.
- Reinstall the mat. Ensure that the rubber mat locks into place with its groove below the edge of the vibration plate over its entire circumference (Figure 6).

**Tip:**

Regularly apply a thin coat of silicone to the mat. This will extend the mat's service life and will make it easier to remove hardened residue.

### 6.2 Cleaning the Housing

**Always disconnect the device from the power supply before cleaning the housing!**

**Never clean the device under running water and never submerge it.**

**NEVER use cleansers containing solvents to clean the device!**

Wipe the housing with a damp cloth.

The control levers can be removed to clean the housing:

- Remove the Allen key from its holder on the control lever (Figure 3).
- Loosen the fastening screw (Figure 7).
- Pull the control levers off (Figure 8).
- Wipe the housing with a damp cloth.
- Reinsert the control levers.
- Retighten the fastening screw (Figure 9).
- Return the Allen key to its holder.

### 6.3 Fuse Replacement

**Always disconnect the device from the power supply before replacing the fuse.**

When replacing fuses, only use fuses with the same rating as the original.

Refer to the "Technical Specifications" section.

**NEVER use fuses with a higher rating!**

The fuse holder (E) is located on the rear of the unit next to the power cord socket (Figure 2).

- Unplug the power cord.
- Release and pull the fuse holder out (Figure 10).
- Remove the faulty fuse and replace it (Figure 11).
- Push the fuse holder back in place until its locks on both sides.

### 6.4 Spare Parts

You can find components subject to wear and the spare parts on the spare part list in the internet at [www.renfert.com/p918](http://www.renfert.com/p918).

Enter the following item number: 14600000 oder 14601000.

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

Serial number and date of manufacturing are shown on the type plate of the unit.

## 7 Error List

Error	Possible cause	Corrective action
The device generates no vibrations.	<ul style="list-style-type: none"> <li>Power cord not properly plugged in at the device.</li> <li>Blown fuse.</li> <li>No power at the wall outlet.</li> </ul>	<ul style="list-style-type: none"> <li>Check that the power cord is securely seated in the socket.</li> <li>Replace the fuse (refer to Sec. 3.4).</li> <li>Test the wall outlet for power.</li> </ul>
No difference in the vibration intensity between the switch settings Fine / Coarse.	<ul style="list-style-type: none"> <li>Faulty diode (short circuit).</li> </ul>	<ul style="list-style-type: none"> <li>Return the device for service.</li> </ul>
No function in coarse oscillation mode.	<ul style="list-style-type: none"> <li>Faulty diode (interruption).</li> <li>Faulty switch.</li> <li>Faulty switch connection wiring.</li> </ul>	<ul style="list-style-type: none"> <li>Return the device for service.</li> </ul>
No function in fine oscillation mode.	<ul style="list-style-type: none"> <li>Faulty switch.</li> <li>Faulty switch connection wiring.</li> </ul>	<ul style="list-style-type: none"> <li>Return the device for service.</li> </ul>
No effect at the different power level settings.	<ul style="list-style-type: none"> <li>Faulty switch.</li> <li>Faulty switch connection wiring.</li> </ul>	<ul style="list-style-type: none"> <li>Return the device for service.</li> </ul>
Device extremely loud during operation.	<ul style="list-style-type: none"> <li>Loose internal components.</li> <li>Loose rubber mat.</li> </ul>	<ul style="list-style-type: none"> <li>Return the device for service.</li> <li>Attach the rubber mat correctly (see Chap. 1.2).</li> </ul>

## 8 Technical Data

Mains voltage:	<b>230 V, 50 Hz / 220 V, 60 Hz</b> <b>120 V, 60 Hz</b> <b>100 V, 50 Hz / 60 Hz</b>
Power consumption:	185 VA (230 V) 170 VA (120 V / 220 V) 190 VA (100 V)
Device fuse *):	2 x 1,6 A (T)
Dimensions (W x H x D):	275 x 140 x 220 mm (10.8 x 5.5 x 8.7 inch)
Weight:	6,4 kg [14.1 lbs]

\*) Fuse in accordance with IEC60127

## 9 Warranty

Provided the unit is properly used, Renfert **warrants** the all components of the *Vibrax* for a **period of 3 years**.

Warranty claims may only be made upon presentation of the original sales receipt from the authorized dealer.

Parts which are subject to natural wear and tear (wear parts) and consumables are excluded from the guarantee. These parts are marked in the spare part list.

The warranty is voided in case of improper use; failure to observe the operating, cleaning, maintenance, and connection instructions; in case of independent repairs or repairs by unauthorized personnel; if spare parts from other manufacturers are employed, or; in case of unusual influences or influences not in compliance with the utilization instructions.

Warranty service shall not extend the original warranty.

# 10 Disposal Information

EN

## 10.1 Instrument Disposal

The instrument must be disposed of by a specialist facility. This specialist facility must be informed of any hazardous residue in the instrument.

### 10.1.1 Disposal Information for EU Nations

In order to maintain and protect the environment, to prevent environmental contamination and to improve the reutilisation of raw materials (recycling), the European Commission has established guidelines according to which manufacturers must take back electrical and electronic devices in order to forward them to a regulated disposal or recycling system.

**Within the European Union, devices identified with this symbol may therefore not be disposed of with unsorted municipal solid waste:**

For more information regarding proper disposal please apply at your local authorities.

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declare under our sole responsibility, that the product(s)

Name	Part No.
Vibrax	1830 0000 / 1830 0000 03

is (are) in conformity with the relevant regulatory requirements by compliance with the UK designated standards.

#### UK legislation

*Electromagnetic Compatibility Regulations 2016*

*Supply of Machinery (Safety) Regulations 2008*

*The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012*


#### Designated standards

EN 61326-1:2013

EN ISO 12100:2010

EN 61010-1:2010/A1:2019/AC:2019-04

EN IEC 63000:2018

Tilo Burgbacher,  
Engineering Director

Hilzingen, 02.11.2022

The Technical Documentation has been retained by Renfert GmbH

**Hochaktuell und ausführlich auf ...  
Up to date and in detail at ...  
Actualisé et détaillé sous ...  
Aggiornato e dettagliato su ...  
La máxima actualidad y detalle en ...  
Актуально и подробно на ...**

**[www.renfert.com](http://www.renfert.com)**

Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany  
Tel.: +49 7731 82 08-0 • Fax: +49 7731 82 08-70  
[www.renfert.com](http://www.renfert.com) • [info@renfert.com](mailto:info@renfert.com)

Renfert USA • 3718 Illinois Avenue • St. Charles IL 60174/USA  
Tel.: +1 6307 62 18 03 • Fax: +1 6307 62 97 87  
[www.renfert.com](http://www.renfert.com) • [info@renfertusa.com](mailto:info@renfertusa.com)  
USA: Free call 800 336 7422

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**Renfert GmbH, Industriegebiet, 78247 Hilzingen / Germany**