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

High Performance Compressed Air Filter

Nr. 2940-0000

ENGLISH

1. Description / Applications

You have just purchased a high quality product by Renfert. In order to be able to enjoy your high performance compressed air filter for a long time to come, please follow the instructions below.

The high performance compressed air filter should keep the compressed air free from condensation, oil residue and other dirt. Only clean compressed air can guarantee trouble free operation and an extended durability of the connected blasting units.

2. Hazard Information

2.1 Employed Symbols

The following symbols are used in both these instructions as well as on the unit itself:

Hazard

There is an imminent risk of injury.

Attention

Failure to observe this information may result in damage to the unit.

2.2 Hazard Information

- Make sure the filter containers do not come into contact with diluting agents, trichloride, kerosene or other aromatic carbon hydrides (solvents, etc. can damage or even destroy the seals and viewing windows).

- Use only neutral cleaners for cleaning.
- Check container during every maintenance period.
- Containers showing signs of mechanical damage must be replaced by new containers.
- Suitable only for use in compressed air systems. Do not use if the pressure and temperature exceed the values indicated in the technical data.
- Disposal of damaged components as per regulations.

2.3 Liability Exclusion

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- The product is employed for any purposes other than those cited in the operating instructions;
- The product is altered in any way other than those alterations described in the operating instructions;
- The product is repaired by other than an authorized facility or if any but Renfert OEM parts are employed;
- The product continues to be employed, despite obvious safety faults;
- The product is exposed to mechanical shocks or is being dropped.

3. Installation / Operation

The high performance compressed air filter has a standard distributor (D) for the connection of three blasting units. The connection to the compressed air net should be made with the hose and nipple connection (F) on the filter. The air flows in the direction of the throughput direction arrows visible on the casing.

! Use only within the allowed operational limits (see technical data)!

Mount unit in a protected yet easily accessible location. Avoid mechanical loading of input lines and screw fittings.

Install the unit in the vicinity of the units to be provided with the compressed air.

The waste condensation will automatically escape upon start-up.

4. Handling

Drainage is done automatically.

5. Troubleshooting Guide

Fault	Cause	Remedy
Air is oily.	Filter is clogged.	Change filter.
Performance is too low.	Filter is clogged.	Change filter.

6. Cleaning / Maintenance

Maintenance of the high performance compressed air filter should be carried out at least once a year. Unusually arduous operational conditions can shorten this recommended maintenance interval. Wear and tear of the active carbon is realized by a discolouring (from black to reddish).

Before starting maintenance, disconnect the high performance filter from the pressure system and ventilate completely.

To facilitate maintenance, the filter can be removed from the installation system. To do so, unscrew the clamping ring (C) and remove the filter toward the bottom from the assembly frame (E). Unscrew the filter container (G) by turning clockwise and remove along with the O-ring (K). Now both filter elements (H + J) can be unscrewed.

Unscrew the transparent filter container (A) by turning counter-clockwise and loosen the active coal filter (B) along with the O-ring (L).

The air flow resistance of the filter element increases according to the amount of dust particles retained. The filter elements must be replaced if they are dirty or damaged. Cleaning of the elements is not possible. Check parts for damage, foreign particles and impurities and reassemble in reversed order. Make sure that the clamping ring (C) is properly placed onto the assembly frame and that the throughput direction arrows on the unit and mounting frame point in the same direction. During reassembling, lightly oil all the O-rings.

7. Spare Parts

You can find components subject to wear and the spare parts on the spare part list in the internet at www.renfert.com/p918.

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

8. Guarantee

Renfert extends a **3 year guarantee** on all parts of the high performance filter. Excluded from this guarantee are parts that are subject to normal wear and tear (filter inserts: pre-filter, fine filter and active coal filter). The guarantee is null and void if claimed damage or defect is the result of improper use, neglect of the operational, cleaning, maintenance and connection regulations, of self-performed repairs or of repairs performed by unauthorized personnel, of use of spare parts from other suppliers or of un-usual or improper influences as determined by the regulations for use. Repairs effected under guarantee do not extend the guarantee period.

9. Technical Specifications

Connection pressure: max. 10 bar, 145 psi

Operational temperature: - 20 °C to + 65 °C

Pre-filter element: 40 µm

Fine filter element: 0.01 µm

Residue oil content: max. 0.003 ppm at 21 °C

Air quality, ISO 8573-1: class 1.7.1

Throughput: 420 l/min (at 6.3 bar)

Dimensions: 210 x 280 mm

Weight: 2.6 kg

10. Parts supplied

- 1 High performance filter unit
- 1 Mounted distribution block
- 1 Connection hose with nipple connection
- 1 Wall mounting unit
- 1 User's manual

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Renfert GmbH • Industriegebiet • 78247 Hilzingen/Germany