
making work easy

- DE
- EN
- FR
- IT
- ES
- PT
- TR
- RU
- PL
- ZH
- JA
- KO
- AR
- UK

Sicherheitsbrenner eco

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Made in Germany

21-6715 06112019

Sicherheitsbrenner eco

934 0100 / 935 0100 / 935 0200

DEUTSCH

1. Einleitung

1.1 Beschreibung

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1. Introduction

1.1 Description

This safety burner eco, is a stable and reliable laboratory burner.

As soon as the flame goes out, within 90 seconds a safety valve activates to create an automatic gas lock. This stops further gas from flowing. The gas consumption of this burner is low because you are able to control the gas flow with a gas valve according to your requirements.

1.2 Symbols used

In these instructions for use and on the burner itself, you will find these symbols with the following meanings:

Danger

Immediate danger. Observe the accompanying documentation!

Warning

By not following the instructions the burner could become damaged.

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► **List, particular attention should be paid:**

- List
- List according to importance

⇒ Handling instructions / required action / input / operation sequence:

You will be asked to carry out the given task in the specified order.

◆ Result of the task / reaction of the unit / reaction of the programme:

The unit or the programme does not react to your action, or due to another incident.

Other symbols are described during use.

2. Safety

2.1 Intended use

Only for use with:

Liquid gas:	Burner No. 934-0100
Natural gas E:	Burner No. 935-0100
Natural gas E:	Burner No. 935-0200

Note:
Previously: Natural gas N (H)
Today: Natural gas E

The burner has been certified and released in conformance with the DVGW:
 Product identification number: DG-2411CO 0243

2.2 Hazard Statements and Warnings

If the burner is not used according to the supplied instructions for use, then the intended safety feature can no longer be guaranteed.

- ▶ The burner must be used under constant supervision.
- ▶ Even when the gas valve (3) on the burner is closed, the gas flow is not securely shut off.
- ▶ After use, the gas supply should be turned off immediately via the gas shut-off valve at the laboratory workbench (gas isolation valve) or by separating the gas connection safety valve.
- ▶ In classrooms, laboratory burners may only be used after turning on all the lined up shut-off systems (central room shut-off, group shut-off system, shut-off valve).
 Close attention must be paid to the DVGW worksheet G 621 (German legislation) when installing and using gas apparatus in laboratories and science classrooms.
- ▶ Should a strong smell of gas be detected, close the gas shut-off valve at the laboratory workbench immediately and initiate the corresponding safety measures (see DVGW law and valid national regulations).
- ▶ Always pay attention to the accident prevention regulations stipulated by the respective trade organisations!
- ▶ Further information concerning connection and operational regulations can be obtained from your local gas supply company and from your local gas installation specialist.
- ▶ Never clean the burner using boiling water or with compressed air!
- ▶ The burner flame can cause burns or surroundings to catch fire, keep your distance!
- ▶ The burner / mixing tube is hot!
- ▶ Information on REACH and SVHC is available on our website www.renfert.com, in the Support area.

2.3 Authorized Persons

Operation and maintenance of the unit may only be carried out by trained personnel.

2.4 Disclaimer

Renfert GmbH disclaims any responsibility for compensation or warranty if:

- ▶ the product is used for any purpose other than the intended use stipulated in the instructions.
- ▶ the product is altered in any way – other than the changes described in the instructions for use.
- ▶ the product is not repaired by a specialist or if third-party spare parts are installed, not original Renfert parts.
- ▶ the product is used further despite obvious safety deficiencies or damage.
- ▶ the product has been damaged by mechanical impact or if the product has been dropped.

3. Product Description

3.1 Component groups and functioning elements

- 1 Mixing tube
- 2 Sensor
- 3 Gas valve
- 4 Knurled nut
- 5 Hose spout
- 6 Actuation knob for ignition safeguard
- 7 Collection pan (only 934 0100 and 935 0100)
- 8 Sensor holder
- 9 Clamping screw
- 10 Allen key


Fig. 1

3.2 Supplied in Delivery

- 1 Safety Burner eco
- 1 Instructions for Use
- 1 Collection pan (only 934 0100 and 935 0100)
- 1 Allen key

4. Initial Operation

The initial operation of the burner may only be carried out by qualified personnel according to DVGW stipulations or valid national regulations.

Please observe the following:

- ⇒ Check if the type of gas and the connection pressure of the burner correspond to the gas supply.
- ⇒ Check the fit of the connection pipe and that it is not damaged.
- ⇒ Close the gas shut-off valve at the laboratory workbench (gas isolation valve).
- ⇒ Slide the certified gas hose (according to DIN 30664 or DVGW G 5501-P) up to the last groove of the hose spout (5) on the burner. If necessary, secure the hose in position with a clamp or ring clasp (see DVGW worksheet G 621 or relevant national specifications).

5. Operation

5.1 Igniting the flame

Caution / risk of burning!

When igniting the burner, turn the mixing tube away from the user and other objects.

Keep your distance!

The flame could be big, depending on the position of the gas valve!

⇒ Open the gas shut-off valve at the laboratory workbench (gas isolation valve).

⇒ Open the gas valve (3) to approx. 1 turn.

⇒ Press the actuation knob (6) in, so that you can hear the gas flowing.

⇒ Ignite the flame immediately!

⇒ After igniting the flame, keep the actuation knob pressed for another 10 seconds until the flame is burning independently. If the flame goes out, repeat the procedure.

The burner is now ready for use.

5.2 Adjusting the flame

Where? knurled nut (4) on the mixing tube.

How? Soft flame: turn the knurled nut (4) in an upward direction (air supply reduced).

Strong flame: the knurled nut (4) in a downward direction (air supply increased).

The size of the flame can be slightly adjusted by opening / closing the gas valve (3).

6. Cleaning / Maintenance

6.1 Cleaning

Caution:

Never clean the burner by boiling or steam cleaning!

Never use acetone when cleaning!

Do not blow compressed air into the hose connection or the mixing tube on the burner! This could damage the safety valve.

Only the following parts may be cleaned with boiling water or steam cleaner once they have been dismantled:

- ▶ **Mixing tube (1, Fig. 1)**
- ▶ **Sensor holder (8, Fig. 1)**
- ▶ **Collection pan (7, Fig. 1)**

When cleaning the external parts, please observe the following points:

- Only clean the burner once it is completely cool.
- Before cleaning, close the gas shut-off valve at the laboratory workbench (gas isolation valve).
- Do not use chemicals for cleaning.
- Never drill or unscrew the burner nozzle.
- The maximum cleaning temperature is 50 °C (122 °F).

6.2 Collection pan

(only 934 0100 and 935 0100)

The collection pan (7, Fig. 1) can be disassembled for cleaning.

6.3 Cleaning the mixing tube / sensor holder

Disassembly:

- ⇒ Using an allen key (10, Fig. 1), loosen the clamping screw (9) and sensor holder (8).
- ⇒ Pull the sensor holder upwards.
- ⇒ Unscrew the mixing tube (1).

Cleaning:

- ⇒ Clean the mixing tube and sensor holder with boiling water or the steam cleaner.

Assembly:

- ⇒ Screw the mixing tube back on.
- ⇒ Push the sensor holder over the mixing tube and sensor so that they sit approximately in the middle of the upper part of the mixing tube (Fig. 2).
- ⇒ Position the sensor so that the upper edge of the sensor's brass sleeve and the upper edge of the mixing tube are in line (Fig. 2).
- ⇒ Tighten the clamping screw in the sensor holder.
- ⇒ Carry out a function test.


Fig. 2

The allen key (10) can be stored inside the housing in the burner foot.


Fig. 3

6.4 Cleaning the burner nozzle

If the burner nozzle is blocked, it can be dismantled from the mixing tube and carefully cleaned with a soft cloth.

If cleaning in this manner is not possible, then the burner nozzle can only be cleaned and changed by an authorized specialist.

6.5 Spare parts

You can find components subject to wear and the spare parts on the spare part list in the internet at www.renfert.com/p918

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

Serial number and date of manufacturing are shown on the type plate of the unit.

6.6 Maintenance

The burner is maintenance-free.

7. Troubleshooting Guide

Fault	Cause	Solution
Burner will not ignite despite open shut-off valve at the laboratory workbench (gas isolation valve), no audible gas flow.	<ul style="list-style-type: none"> The main gas supply has not been switched on / gas bottle is empty. The actuation knob for automatic ignition has not been pressed. Nozzle is blocked. 	<ul style="list-style-type: none"> Switch the main gas supply on / exchange gas bottle. Keep the actuation knob pressed in. Unscrew the mixing tube; clean the nozzle carefully with a soft cloth. Do not unscrew the nozzle!
Burner will not ignite despite clearly audible gas flow.	<ul style="list-style-type: none"> Incorrect connection values. Caution: faulty pressure reducer. Defect in the sensor. 	<ul style="list-style-type: none"> Check the connection values and type of gas. Exchange pressure reducer. Send the burner in for repairs.
Smell of gas when burner is ignited.	<ul style="list-style-type: none"> Leaky connections. Leaky burner. 	<ul style="list-style-type: none"> Check connections. If the burner is leaky, do not repair yourself! Send the burner to be examined.

8. Technical Data

Type	934 0100	935 0100 / 935 0200
Gas type:	Liquid gas G 30 / 50 mbar	Natural gas E G 20 / 20 mbar
Dimensions (W x H x D):	141 x 100 x 121 mm [5.6 x 3.9 x 4.8 inch]	
Weight, approx.:	460 g [16.2 oz]	
Connection hose:	according to DIN 30664 or DVGW G5501-P	
Connection pressure:	50 mbar [0.725 psi]	20 mbar [0.290 psi]
Consumption:	50 g/h [1.76 oz/h]	48 l/h [1.7 ft ³ /h]
Watt:	650 W	450 W
Nozzle Ø:	0,34 mm [13.4 mil]	0,55 mm [21.7 mil]
Operating mode:	Continuous operation	
max. ambient temperature:	40 °C [104 °F]	

9. To change the gas type

The burner corresponds to the CE 30665 part 1 requirements and can be converted to use a different type of gas by changing the nozzle.

The type of gas and performance can be found in the technical data.

- The liquid gas version has “**Liquid gas 50 mbar**” written on the identification plate.
- The natural gas version has “**Natural gas E 20 mbar**” written on the identification plate.

When the nozzle is changed, the identification plate must also be changed.

The conversion sets can be found in the spare part list.

The nozzle may only be changed by qualified personnel according to the DVGW or valid national regulations.

10. Guarantee

All burners are checked before distribution. If you discover any functioning fault, please send the burner in for repair with a description of the fault. Do not carry out any repairs yourself. The burner nozzle and mixing tube have been carefully aligned and are compatible with one another; they must not be altered as they would then no longer correspond with the regulations.

When the equipment is used as intended and stipulated, Renfert warrants all parts of the unit with a **guarantee of 3 years**.

An original sales invoice from a specialist dealer is a precondition for a warranty claim.

Parts which are susceptible to regular wear and tear, such as consumables, are excluded from this warranty.

The guarantee is invalid after improper use, incompliance with the instructions for use regarding operation, cleaning, maintenance and connection, after repairs which were not carried out by specialists, if spare parts from other manufacturers have been installed and unpredictable or unpermitted external influences.

Guarantee claim do not lengthen the original guarantee duration.

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Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany