# Waxlectric I + II

No. 2139 / 2140

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ENGLISH

## 1. Introduction

We are pleased that you have decided to buy the *Waxlectric I/II*. This unit sets new standards with regard to function, capability and safety. To guarantee that you enjoy using this product for many years please follow these instructions carefully.

## 2. Description / Applications

*Waxlectric* is an electronically controlled waxing-up instrument. It is available in a 1-channel or 2-channel version both of which are suitable for almost all modelling and wax-up tasks.

- **An external power adaptor** supplies the waxing-up instrument with a reduced voltage of 12 V. This increases the user's safety during everyday use of the *Waxlectric*.
- The individual functions can be actuated via the wear-and-tear free, sealed **touch control panel**.
- The specific benefit of the *Waxlectric* is the **heating element**. This is contained within the modelling tips. The tips therefore respond more quickly to any change in temperature. The strong heat-conducting alloy also increases the thermal dynamics.
- The electronic waxing-up instrument is extremely **easy to maintain**. The modelling tips do not have an external finish. They can be ground according to your requirements. Should any defect with the heating element occur you need only replace the tip.
- The modelling tips can reach a **maximum temperature of 200°C [392°F]**. The working temperature can be controlled linearly from 50°C to 200°C [122°F to 392°F]. Be careful not to set the temperature too high, as some wax is very sensitive to overheating.

## 2.1 Ambient conditions

The device may only be operated:

- in interior areas,
- up to an altitude of 2,000 m above sea level,
- at an ambient temperature range between 5 - 40°C (41 - 104°F),
- at a maximum relative humidity of 80% at 31°C (87.8°F), dropping in a linear manner to 50% relative humidity at 40°C (104°F),
- at mains power fluctuations where the voltage fluctuations do not exceed 10% of the nominal value,
- over voltage category II,
- at pollution degree 2.

## 3. Installation / Operation

→ see illustrations 1-4


## 4. Handling

### 4.1 ON / OFF function:


### 4.2 Setting the temperature:

Setting the unit:



Switch between °C ↔ °F:  
Only possible when switching on !



Please observe that all 230 V units are pre-set to centigrade (°C) and all 120 V units are pre-set to Fahrenheit (°F).



**Setting the temperature:**


**Speed setting to maximum temperature (Waxlectric II version only possible at the right channel):**

**Speed setting to minimum temperature (Waxlectric II version only possible at the right channel):**

**4.3 Changing the tips:**

Make sure the O-ring is present (ill. 6) - it is the seal against wax penetration !

→ see illustrations 5, 6, 7, 1

**4.4 Changing the cork:**

→ see illustrations 5 + 8

## 5. Areas of Caution

- Before switching on check that the unit's voltage matches the main power supply.
- Units with faulty leads or other defects should not be used.
- Do not cover the power units when in use.
- Unroll power supply cable (overheating may occur in event of short circuit).
- The tips of the waxing-up instrument reach a working temperature of 200°C [392°F].  
> *Improper use may lead to a danger of burning!*
- The red marked end of the handle of the instrument gets very hot while working.  
> *Improper use may lead to a danger of burning!*

- Always insert the handpiece with the tip in the predetermined holder on the unit in order to avoid a potential fire.
  - The waxing-up instrument is solely for processing modelling waxes.
- > **Processing other materials can endanger the user!**

#### **5.1 Exclusion of liability**

Renfert GmbH rejects any claims for compensation as well as warranty claims in case of:

- the product is used for any other purpose than those mentioned in the instruction manual
- the product is changed or amended in any way beyond those mentioned in the instruction manual
- the product is repaired by a non-authorized party or the product is not fitted with original Renfert replacement parts
- the product is used continuously despite visible safety faults
- the product is subjected to mechanical jarring or is dropped.

### **6. Cleaning / Maintenance**

**Clean the unit from time to time with a soft cloth.** Never wash the unit! Rub the tips clean after use. This will prevent the wax burning into the tips and the tips losing their optimum heating capacity.

### **7. Spare Parts**

You will find the order numbers for wear-and-tear and exchange parts in the enclosed spare parts list.

### **8. Guarantee**

When used properly Renfert guarantees all parts of the *Waxlectric* – with the exception of wear-and-tear parts – for 3 years. Parts given to natural wear are excluded from this guarantee (e.g., all tips, cork fittings, etc.). The guarantee is void for damage caused by misuse, disregard of instructions for use, cleaning, maintenance and connection, as well as for self-repaired units or repair by untrained personnel, for use of spare parts other than from the manufacturer and influences outside the normal practice or unauthorised in the instruction manual. Services under guarantee do not extend the guarantee.

### **9. Technical Specifications**

Max. power requirement:  
5 W (*Waxlectric I*), 10 W (*Waxlectric II*)  
Min. working temperature: 50°C [122°F]  
Max. working temperature: 200°C [392°C]  
Power supply:  
12 V DC unregulated (check polarity!)  
Power supply - transformer:  
230 V AC, 50 Hz, 12 V DC  
120 V AC, 60 Hz, 12 V DC  
240 V AC, 50 Hz, 12 V DC

### **10. Parts supplied**

- 1 control unit with integrated holder for heating handle and tips
- 1 handpiece (or 2 handpieces with *Waxlectric II*)
- 1 modelling tip no. 2141-0102 (*Waxlectric I* version) resp. tips no. 2141-0101 + 2141-0103 (*Waxlectric II* version)
- 1 power supply
- 1 instruction manual with spare parts list


### **11. Units available**

No. 2139-0000 *Waxlectric I*, 230 V  
No. 2139-1000 *Waxlectric I*, 120 V  
No. 2139-2000 *Waxlectric I*, 240 V

No. 2140-0000 *Waxlectric II*, 230 V  
No. 2140-1000 *Waxlectric II*, 120 V  
No. 2140-2000 *Waxlectric II*, 240 V

## 12. Accessories

No. 2141-0101	1 probe, small	No. 2141-0107	1 Hollenback, horizontal
No. 2141-0102	1 probe, medium	No. 2141-0108	1 Hollenback, vertical
No. 2141-0103	1 probe, large	No. 2141-0109	1 mini probe – angled
No. 2141-0104	1 blade, narrow	No. 2141-0110	1 lance shaped blade
No. 2141-0105	1 blade, wide	No. 2141-0111	1 beaver tail blade
No. 2141-0106	1 blade, chamfer		

## 13. Troubleshooting Guide

Fault	Cause	Remedy
<b>No display, no function.</b>	<ul style="list-style-type: none"> <li>• Power supply interrupted.</li> <li>• Defect control unit.</li> <li>• Defect power adaptor.</li> <li>• Power supply too low.</li> <li>• Short circuit in the tip.</li> </ul>	<ul style="list-style-type: none"> <li>• Check power supply and reconnect.</li> <li>• Exchange via service.</li> <li>• Exchange with replacement part (see spare parts list).</li> <li>• Connect the right power adaptor (see spare parts list).</li> <li>• Clean the tips and handle (blow out).</li> <li>• Exchange the tip.</li> <li>• Heating handle defect --&gt; exchange with spare part (see spare parts list).</li> </ul>
<b>Unit functions but no heat in the tip.</b>	<ul style="list-style-type: none"> <li>• Tip is not screwed on correctly.</li> <li>• The contact between tip and handle is dirty.</li> <li>• The tip is defect.</li> <li>• The handpiece cable is not connected correctly to the unit.</li> </ul>	<ul style="list-style-type: none"> <li>• Screw in the tip hand-tight.</li> <li>• Clean the contact (blow out).</li> <li>• Replace the O-ring (see spare parts list).</li> <li>• Replace with the new tip (see point 12 / Accessories).</li> <li>• Insert the handpiece cable into the unit until it clicks.</li> </ul>
<b>Temperature memory not functioning.</b>	<ul style="list-style-type: none"> <li>• The memory chip is defect.</li> </ul>	<ul style="list-style-type: none"> <li>• Send for repair.</li> </ul>

***Tips & Tricks:***

- When "dipping" a thin probe with a rounded back into cold wax the wax will melt more quickly. A larger amount of wax can be taken up and also reduces the mechanical load of the relatively sensitive probe.
- Speedier and more effective modelling is achieved (30% saving on time!) when you pre-warm the wax in the *Vario E*. The wax can be used more sparingly and cools down quicker. This reduces shrinking to a minimum.
- As the tip inserts are not coated they can – starting from their basic shape – be ground, bent and polished as you require.
- The cable of the handpiece can be guided (e.g. over a specific fixed point (hook etc.) so that the handpiece hangs down from above and thus is at the ideal holding level. Using this cable/handpiece hanger a speedier grasp or change of handpiece is achieved.
- If you set the temperature to 50°C [122°F] you can use a slightly warmed *Waxlectric* blade to rub the wax crown border onto the die.

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