
making work easy

- DE
- EN
- FR
- IT
- ES
- PT
- TR
- RU
- PL
- CS
- AR

# Sicherheitsbrenner

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Made in Germany

21-9021 14012020

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# 1 Description / Applications

By purchasing a Renfert *safety burner* you now own a sturdy laboratory burner. A safety valve automatically turns the gas off within 90 sec. of the flame being extinguished. Thus no more gas can flow out. In addition, a tipping lever enables you to run the burner on main or pilot flame thus guaranteeing low gas consumption.

► **The burner is supplied in the following versions:**

No. 933-0100 > burner for liquid gas

No. 932-0100 > burner for natural gas E

► **Please note:**

Old

New

natural gas N (H)

natural gas E

The burner is type-tested and approved by the DVGW (German Association of Gas and Water Installers).  
Registr. number: NG-2211AO0738

## 2 Installation / Operation

The burner may only be installed by qualified staff in accordance with the DVGW regulations.

- Insert and fix the gas lever supplied with the burner into the screw thread hole on the burner's operating element.
- Check whether the type of gas and operating pressure required for the burner are compatible with your gas supply.
- Inspect the connection pipe to ensure it is correctly fitted and in perfect condition
- Shut off the gas tap on the laboratory bench.
- Push the gas hose (DIN 30664) onto the inlet nozzle of the *safety burner* right up to the last groove.
- Angle the burner's lever, preferably to position A, (pilot flame - see diagram) and open the gas cock for the lab bench.
- Press the gas lever in the middle of the burner so that the gas flow can be heard. Ignite the flame immediately!
- After the flame has been lit, hold the gas lever down for approx. 10 seconds until the flame burns independently. If the flame extinguishes repeat the process.

**Caution!**

**Always ignite the gas burner from the side!**

The burner is ready for use.

## 3 Areas of Caution

- This laboratory burner does not have an „off“ position.
- The laboratory burner must only be used under constant supervision. After use the gas supply should be turned off immediately via the gas shut-off valve or by separating the gas connection safety valve. In classrooms, laboratory burners may only be used after turning on all the lined up shut-off systems (central room shut-off, group shut-off system, shut-off valve). Close attention must be paid to the DVGW Work-sheet G621 (German legislation) when installing and using gas apparatus in laboratories and science classrooms.
- Should a strong smell of gas be detected, close the gas shut-off valve immediately and initiate the corresponding safety measures (see DVGW regulations).
- Always pay attention to the accident prevention regulations laid down by the respective trade organisations!
- Further information concerning connection and operational regulations can be obtained from your local gas supply company as well as from your local gas installation specialist.
- Never use water or compressed air to rinse or blow out burners! This could flush out / blow out the sealant grease causing the burner to leak.


### 3.1 Exclusion of liability

Renfert GmbH rejects any claims for compensation as well as warranty claims in case of:

- the product is used for any other purpose than those mentioned in the instruction manual
- the product is changed or amended in any way beyond those mentioned in the instruction manual
- the product is repaired by a non-authorized party or the product is not fitted with original Renfert replacement parts
- the product is used continuously despite visible safety faults

## 4 Handling

You can choose the size of the *safety burner's* flame using the tipping lever.

- Position A: Ignition or pilot flame burns  
 Position B: Main flame and ignition flame burn  
 Position C: Main flame burns

When choosing either position A or C always **turn the tipping lever to its limit**. By doing this you can identify quickly and unmistakably which position the burner was in should the flame ever go out. This will also prevent a flash-back in the burner's nozzle.

### Regulating the main flame:

- Where? The knurled nut on the large burner pipe.
- How? **Gentle flame:** turn the knurled nut upwards (reducing air supply).  
**Harsh flame:** turn the knurled nut downwards (increasing air supply).


### Regulating the pilot flame:

- Where? The sleeve over the air supply holes on the small burner pipe.
- How? **Gentle flame:** slide the sleeve over the air supply holes.  
**Harsh flame:** slide the sleeve away from the air supply holes (unstable flame).

The ignition flame should be as stable, and therefore as smooth, as possible in order to ensure a reliable burner operation. Adjust the flame so that it cannot be extinguished unintentionally.

### Adjusting the burner position:

The burner can be positioned at three different angles:

To adjust the burner to the angle which suits you the best, simply lift it out of its base and reinsert it in the required position. The burner is thus equally well suited for left-handed and right-handed users.

The burner should not be used in the 90° position when working with wax, since dripping wax may block the jets and cause the burner to malfunction.

## 5 Troubleshooting Guide

Fault	Cause	Remedy
<b>Burner will not ignite despite open shut-off valve, no audible flow.</b>	<ul style="list-style-type: none"> <li>• The main gas supply has not been switched on / gas bottle is empty.</li> <li>• Gas lever not pressed.</li> <li>• Nozzle is blocked.</li> </ul>	<ul style="list-style-type: none"> <li>• Operational error.</li> <li>• Operational error.</li> <li>• Screw off burner tubes and carefully clean the jets with a soft cloth.  <b>Do not screw out the jets!</b></li> </ul>
<b>Burner will not ignite despite clearly audible gasflow.</b>	<ul style="list-style-type: none"> <li>• Incompatible connecting valves.</li> <li>• <b>Caution:</b> faulty pressure reducer.</li> <li>• Faulty thermoelement.</li> </ul>	<ul style="list-style-type: none"> <li>• Check connection pressure and gas type.</li> <li>• Replace pressure reducer.</li> <li>• Send the burner in for repair.</li> </ul>
<b>Smell of gas when burner is ignited.</b>	<ul style="list-style-type: none"> <li>• Leaky connections.</li> <li>• Leaky burner.</li> </ul>	<ul style="list-style-type: none"> <li>• Examine connections.</li> <li>• Do not repair yourself. Send the burner to be examined.</li> </ul>

## 6 Cleaning / Maintenance

The *safety burner* is maintenance-free.

To remove external dirt please note the following points:

- Close the shut-off valve before cleaning.

### Attention:

**Never use water or compressed air to rinse or blow out burners! This could flush out / blow out the sealant grease causing the burner to leak. Always comply with the operating instructions!**

- Do not use chemical cleaning fluids.
- Use only a damp leather cloth to clean the frame.
- Never bore or unscrew the burner nozzles.

### Attention:

**do not blow compressed air down the burner's connector pipe!**

**The grease sealant will be removed or a sealing ring may be dislodged and the burner will leak.**

- The maximum temperature for cleaning is 50 °C (122 °F).

### Dismantling for cleaning:

#### Attention:

**The burner pipes and contacts are connected by a metal joint. This must not be bent!**

- Unscrew the large and small burner pipes.
- Pull the contacts out of the joint from below.
- The burner pipes can now be boiled out separately.

### Reassembly after cleaning:

- Slide the contacts into the joint from below.
- Screw on the burner pipes.
- **Attention:**  
reposition the contacts in their original position. Do not bend!
- Control the functioning.

- ▶ **If the burner nozzle becomes plugged, disassemble the tubes than hold the burner with the nozzle facing down. Use a hair dryer or other device generating hot air, melt the wax plug and allow the molten wax to completely run out.**
- ▶ **If this method fails to clean the nozzle, it must be either cleaned or replaced by an authorized technician.**

## 7 Spare Parts

You can find components subject to wear and the spare parts on the spare part list in the internet at [www.renfert.com/p918](http://www.renfert.com/p918).

Enter the following item number: 9330100 or 9320100

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

Serial number and date of manufacturing are shown on the type plate of the unit.

## 8 Guarantee

All *safety burners* are examined carefully by us before dispatch. Should any malfunctions occur please send in your burner for repair with a description of the problem. Do not try and repair the unit yourself!

The burner's nozzles and pipes are harmonised and should not be changed as they will then no longer comply with regulations.

When used properly Renfert **guarantees** all parts of the *safety burner* for **3 years**.

The original sales invoice of the specialist depot is required for a claim under guarantee. Parts given to natural wear and tear are excluded from this guarantee.

The guarantee is void for damage caused by misuse, disregard of the instructions for use, cleaning, maintenance and connection, for self-repaired units or repair by untrained personnel, for use of spare parts from other manufacturers and on influences outside the normal practice or unauthorised in the instruction manual.

Services under guarantee do not extend the guarantee.

## 9 Switching to a different type of gas

The *safety burner* conforms with the regulations set out in DIN 30665 Part 1 and enables a switch to another type of gas by changing the set of burner nozzles.

**Type of gas and power:**

Type of gas	Rated power					
	Small nozzle Ø	Large nozzle Ø	Load	Load large tube	Load small tube	Product ref. mark
Liquid gas G 30 / 50 mbar	0,18 mm [7.1 mil]	0,34 mm [13.4 mil]	0,75 kW	0,60 kW	0,15 kW	933-0100
Natural gas E G 20 / 20 mbar	0,32 mm [12.6 mil]	0,55 mm [21.7 mil]	0,65 kW	0,50 kW	0,15 kW	932-0100

- The natural gas E version is distinguished with the sticker „Erdgas E 20 mbar“.
- The liquid gas version is distinguished with the sticker „Flüssiggas 50 mbar“.

The predetermined lab burner can be switched over to natural gas E and liquid gas by changing the set of nozzles and its respective sticker. The relevant conversion sets are given in the table above.

Each set consists of a small nozzle, a large nozzle and a sticker naming the respective type of gas. The conversion sets are featured in the spare parts list.

**The changeover of the nozzles must only be carried out by qualified specialists in accordance with DVGW regulations (German legislation - check local regulations).**

## 10 Technical data

	Type	
	933-0100 Liquid gas	932-0100 Natural gas E
Height (depending on the angle):	113 / 122 mm [4.45 / 4.80 inch]	
Weight (with burner base):	640 g [22.4 oz]	
Connection hose:	According to DIN 30664	
Supply pressure:	50 mbar	20 mbar
Consumption:	48 g/h [1.69 oz/h]	63 l/h [2.2 ft <sup>3</sup> /h]
Output:	750 W	650 W
Mode of operation:	Continuous operation	
Ambient temperature, max.:	40 °C [104 °F]	
Dimensions (W x H x D):	160 x 122 x 60 mm [6.30 x 4.80 x 2.36 inch]	

## 11 Parts supplied

- 1 Burner
- 1 Gas lever set
- 1 Burner base
- 1 Instruction manual

## 12 REACH

Information on REACH and SVHC is available on our website [www.renfert.com](http://www.renfert.com), in the Support area.

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Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany