making work easy

- de
- en
- fr
- it
- es
- pt
- tr
- ru
- pl
- sv
- cs
- zh
- ko
- ar
- uk

# SILENT TS2

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Made in Germany

21-6670 30072024

# **SILENT TS2**

## **Nr. 2930 0050 / 2930 1050**

DEUTSCH

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## **Inhalt**

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ENGLISH

## Content

1. Introduction.....	2
1.1 Symbols .....	2
2. Safety .....	3
2.1 Intended Use .....	3
2.2 Improper Use.....	3
2.3 Ambient Conditions for Safe Operation.....	3
2.4 Ambient Conditions for Storage and Transport.....	3
2.5 Hazard and Warning Information.....	4
2.5.1 General Information .....	4
2.5.2 Specific Information .....	4
2.6 Authorized Persons.....	5
2.7 Disclaimer.....	5
3. Product Description.....	5
3.1 General Description .....	5
3.1.1 Areas of Application .....	5
3.2 Components and Functional Elements .....	6
3.3 Scope of Delivery .....	7
3.4 Accessories.....	7
4. Setting Up .....	8
4.1 Unpacking.....	8
4.2 Setup .....	8
4.2.1 External Exhaust Air Route .....	8
4.3 Electrical Connection.....	8
4.4 Compressed Air Connection.....	8
4.5 Installing the Pinch Valves .....	9
4.5.1 Suction Pipe Alteration .....	9
4.6 Connection to the Extraction Point.....	10
4.7 Electrical Connection.....	10
5. Operation .....	10
5.1 Switching the Unit On .....	10
5.1.1 Stand-By .....	11
5.2 Operation.....	11
5.2.1 Extraction .....	11
5.2.2 Extraction force .....	11
5.3 Programming .....	12
5.3.1 Self Diagnosis .....	12
5.3.2 Audible Signal .....	13
5.3.3 Dust Bag Full.....	13
5.3.4 Start-Up Threshold .....	14
5.3.5 Shut-Off Delay (3 - 30 sec.) .....	14
5.4 Operation with the SILENT flow sensor .....	15

6. Cleaning / Maintenance.....	16
6.1 Replacing the Dust Bag .....	16
6.2 Filters.....	17
6.2.1 Replacing the Fine Particle Filter .....	17
6.2.2 Replacing the Exhaust Filter .....	17
6.2.3 Replacing the Electronics Filter .....	17
6.2.4 Change Compressed filter .....	18
6.3 Self-Diagnosis.....	18
6.4 Safety Mechanism .....	18
6.5 Factory Set Parameter.....	18
6.6 Spare Parts .....	19
7. Troubleshooting.....	19
7.1 Error codes .....	20
8. Technical Data .....	21
9. Warranty.....	21
10. Disposal Information .....	22
10.1 Disposing of Consumables .....	22
10.2 Disposing of the Unit .....	22
10.2.1 Disposal Instructions for countries in the EU.....	22

en

# 1. Introduction

## 1.1 Symbols

In the instructions for use and on the unit itself you will find these symbols with the following meanings:

**Danger**  
This indicates a direct risk of injury. Consult accompanying documents!

**Electrical current**  
This indicates a risk of hazard due to an electrical current.

**Attention**  
Disregarding this warning may result in damage to equipment.

**Note**  
This provides the operator with useful information to improve and ease use.

The device complies with the requirements of the applicable EU directives.

This product complies with the relevant UK legislation.  
See UKCA Declaration of Conformity in the Internet under [www.renfert.com](http://www.renfert.com).

The device is subject to the EU directive 2002/96/EG (WEEE directive).

► **List, particular attention should be paid**

- List

⇒ Instructions / appropriate action / input / operational sequence:

You will be asked to carry out the action in a specified order.

- ◆ Result of an action / reaction of the device / reaction of the program:

The unit or program reacts as a result of your actions or when a specific incident occurs.

Other symbols are explained as they occur.

## 2. Safety

### 2.1 Intended Use

This device is designed to extract dry, non-explosive dust.

This device is intended solely for commercial use in a dental laboratory.

The intended use also includes compliance with the instructions specified by the manufacturer concerning operation, servicing and maintenance.

The extraction unit can be used in combination with a SILENT flow sensor within the scope of measures to protect against exposure to dust \*) as defined by the German ordinance on hazardous substances ("GefStoffV") / respective employers' liability insurance association.

For this purpose, the information on low-dust disposal of the extracted material (chapter "Emptying the dust drawer") and on low-dust replacement of the fine filter (chapter "Replacing the fine filter") in particular must be observed.

**In doing so, both the instructions for use of the extraction unit as well as those of the SILENT flow sensor must be observed.**

\*) Note for customers in Germany: A validated and recognized dust protection measure in accordance with testing principle GS-IFA-M20 (issued by the German Institute for Work and Health IFA) shall only be deemed to have been provided when all components, i.e. the flow rate monitoring unit such as the SILENT flow sensor, the extraction unit, and the detection unit, are validated and recognized.

In the case of multi-workbench extraction units and the use of Y-junctions or extractor switches, every suction channel / every suction point must be fitted with a flow rate monitoring unit.

### 2.2 Improper Use

Fire promoting, easily flammable or explosive materials must not be extracted with the SILENT TS2.

The extraction of liquids, smoldering or burning materials is prohibited.

This device is not intended for private, household use.

Any use other than specified in these instructions is deemed improper and constitutes a misuse of the device. The manufacturer shall not be liable for damages caused by improper use.

Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the device, increase the risk of serious injury and lead to damage to the environment or the device itself.

### 2.3 Ambient Conditions for Safe Operation

The device may only be operated:

- Indoors
- Up to an altitude of 2,000 m above sea level,
- At an ambient temperature of between 5 – 40 °C [41 – 104 °F] \*),
- At a maximum relative humidity of 80 % at 31 °C [87.8 °F], dropping to a linear of up to 50 % relative humidity at 40 °C [104 °F] \*),
- With mains power where the voltage fluctuations do not exceed 10 % of the nominal value,
- Under contamination level 2 conditions,
- Under over-voltage category II conditions.

\*) Between 5 - 30 °C [41 - 86 °F] the device can be operated at a relative humidity of up to 80 %. At temperatures between 31 - 40 °C [87.8 - 104 °F] the humidity must decrease proportionally in order to ensure operational readiness (e.g. at 35 °C [95 °F] = 65 % humidity, at 40 °C [104 °F] = 50 % humidity). The device may not be operated at temperatures above 40 °C [104 °F].

### 2.4 Ambient Conditions for Storage and Transport

For storage and transport the following specifications to ambient conditions apply:

- Ambient temperature - 20 – + 60 °C [- 4 – + 140 °F].
- Maximum relative humidity 80 %.

## 2.5 Hazard and Warning Information

### 2.5.1 General Information

- ▶ If the device is not used in compliance with the supplied instructions, the safety of the device can no longer be guaranteed.
- ▶ The device may only be operated using a mains cable with the country-specific plug system. Any necessary alterations must be carried out by a qualified electrician.
- ▶ The device may only be operated if the information on the identification plate conforms to the specifications of your local mains power supply.
- ▶ The device may only be plugged into outlets which are connected to the protective conductor system.
- ▶ The mains plug must be easily accessible.
- ▶ Disconnect the device from the mains before carrying out work on the electrical parts.
- ▶ Check connection cables (such as power supply cords), tubes and housing (i.e. the key-pad) regularly for damage (i.e. kinks, cracks and porosity) or signs of ageing. Devices with damaged connection cables, tubes or housing parts or other defects must not be operated!
- ▶ Defective devices must be put out of service immediately. Remove the mains plug and ensure the device is not used. Send the device for repair!
- ▶ Only operate the device under supervision.
- ▶ Please observe the national accident prevention regulations!
- ▶ It is the responsibility of the operator that national regulations during operation and regarding a repeated safety inspection of electrical equipment are complied with. For Germany these are the regulation 3 by DGUV (German Statutory Accident Insurance) in relation with VDE 0701-0702 (Association for Electrical & Electronic Technology) .
- ▶ Information on REACH and SVHC is available on our website [www.renfert.com](http://www.renfert.com), in the Support area.

### 2.5.2 Specific Information

- ▶ The German ordinance on hazardous substances (“GefStoffV”) or equivalent national regulations must be observed during operation, cleaning, and maintenance.
- ▶ The mains socket on the device is only designed for purposes specified in the operating instructions. Connecting other devices may cause material damage.
- ▶ Before connecting another electrical appliance to the mains socket, switch the unit off.
- ▶ Read the operating instructions of the other appliance and comply with the safety instructions contained in the document.
- ▶ Please observe the national regulations and permitted exposure to dust in a working environment. Please ask the “National Institute for Occupational Safety and Health” or other responsible authority.
- ▶ Please observe the Safety Data Sheets of the materials to be extracted.
- ▶ Always wear protective gear, when extracting hazardous materials.
- ▶ It is necessary to wear suitable personal protective equipment when emptying the dust drawer or cleaning, depending on the type of extracted material.
- ▶ When disposing of the extracted material or used filter, please observe the local specifications and accident prevention regulations!
- ▶ Make sure the dust drawer is fully closed during operation.
- ▶ Do not operate without a suction hose.
- ▶ Do not extract flammable or explosive gasses, fumes or dust.
- ▶ The following applications involve increased risks and are therefore not permitted:  
When cleaning milling chambers and SLM printing systems by suction, a sufficiently high concentration and purity (i.e. unmixed with other dental dusts such as gypsum, resins) can result in spontaneous ignition or the milling dust exploding or printing powder because of a potential exothermic reaction (e.g. due to oxidation).  
An alternative cleaning method must be selected for the following milling or printing media in particular (such as manual cleaning):
  - Wood
  - Titanium / titanium-aluminum
  - Light metals and light metal alloys (e.g. aluminum, magnesium)
  - Cobalt chrome powder (e.g. for use in SLM systems)If large quantities of light metals such as titanium alloy are processed (e.g. with sandpaper), producing very fine grinding dust, a sufficiently high concentration and purity may result in spontaneous combustion because of a potential exothermic reaction.

- ▶ **Do not extract hot materials.**
- ▶ **Do not extract liquids.**
- ▶ **If the dust extractor is employed to extract hazardous materials, appropriate personal protective gear must be worn and steps must be taken to ensure that the exhaust air is properly ventilated. Please refer to the associated safety data sheets for specific requirements.**
- ▶ **Dispose of extracted material according to local statutory regulations.**

## 2.6 Authorized Persons

Operation and maintenance of the device may only be performed by qualified personnel. Minors and pregnant women may only operate and service the device if they are wearing appropriate protective gear, in particular if the device is being used to extract hazardous materials. Any repairs not specifically described in these operating instructions may only be carried out by a qualified electrician.

## 2.7 Disclaimer

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- ▶ **The product is employed for any purposes other than those specified in the operating instructions.**
- ▶ **The product is altered in any way other than those alterations described in the operating instructions.**
- ▶ **The product is not repaired by an authorized facility or if non-original Renfert parts are implemented.**
- ▶ **The product continues to be used despite obvious safety faults or damage.**
- ▶ **The product is subjected to mechanical impacts or is dropped.**

# 3. Product Description

## 3.1 General Description

This device is a workbench extraction unit for the extraction of dusts generated in a dental laboratory. The extraction unit compiles two suction hoses to which one extraction point each may be connected. The extraction unit can be operated both manually and automatically, depending on the connected, electrical, dust-generating equipment.

### 3.1.1 Areas of Application

In the dental laboratory different areas of application are possible. Depending on the amount of dust to be extracted and the extraction performance available, both suction hoses can either operate together, or if preferred, each tube alternately.

Area of application	
Standard Workbench (WB)	Grinding, separating, polishing, etc. of different dental materials such as alloys, acrylics, porcelains, plasters, etc. in normal amounts using the technician's hand piece.
Work Preparation (WP)	Mainly grinding plasters and acrylics in large amounts using the technician's hand piece
Device (D)	Dust intensive devices such as the fine / recyclable sandblaster, etc.

#### Recommended use of the suction hoses:

- WB + WB      Suction hoses A and B simultaneously
- WB + WP      Suction hoses A and B simultaneously
- WP + WP      Suction hose A or B alternatively
- D + D          Suction hose A or B alternatively
- WB + D        Suction hose A or B alternatively
- WP + D        Suction hose A or B alternatively

## 3.2 Components and Functional Elements

- |    |                              |    |  |
|----|------------------------------|----|--|
| 1  | SILENT TS2                   | 12 | Compressed air connection with compressed air filter (Compressed air tube is sealed) |
| 2  | Key-pad                      | 13 | Motor compartment cover  |
| 3  | On / off switch              | 14 | Exhaust air filter   |
| 4  | Front panel/ Fine filter     | 15 | Pinch valve with extractor connection fitting (rotatable)                            |
| 5  | Front panel/ Dust drawer     | 16 | Suction pipe   |
| 6  | Fine filter                  | 17 | Connection set   |
| 7  | Dust drawer                  | 18 | Suction hose   |
| 8  | Mains cable                  |    |  |
| 9  | Device coupler socket (A, B) |    |  |
| 10 | Device protection switch     |    |  |
| 11 | Electronic filter            |    |  |


Fig. 1

- 3 On / off switch
- 30 Menu key
- 31 Activate self-diagnosis
- 32 Audible signal on / off
- 33 Programming dust drawer full
- 34 Programming start-up threshold
- 35 Programming shut off delay
- 36 7 segment display
- 37 Suction level display
- 38 Suction hose display
- 39 Suction hose keys
- 40 Enter button, save input / display enter button
- 41 ▢ Button / display function button
- 42 Display factory setting
- 43 ⊕ Button / display function button


Fig. 2

### 3.3 Scope of Delivery

- 1 SILENT TS2
- 1 Operating instructions
- 1 Quick reference card
- 1 Dust bag (ready for use in the dust drawer)
- 1 Compressed air tube (sealed)
- 2 Pinch valves with suction port
- 1 Connection set
- 2 Shockproof plug (only for 2930-0050)
- 1 4 m Suction Hose
- 4 PVC end bushings

### 3.4 Accessories

- 2921 0002 Dust bags (5 pieces)
- 2921 0001 Set of rollers
- 2925 0000 Extractor clamp (without glass screen)
- 2925 1000 Glass screen with holder (for the extractor clamp)
- 2921 0003 Set of mufflers, 2 pieces
- 2921 0004 External exhaust air control
- 2934 0007 90° angled connector Silent
- 90003 4305 Tube support adapter
- 90003 4240 Suction hose (3 m)
- 90003 4826 Suction hose, antistatic, 3 m, incl. 2 end mufflers
- 90115 0823 Suction hose (6 m)
- 90215 0823 Suction hose (9 m)
- 90003 4410 H+Hepa filter
- 90003 4430 Suction hose adapter universal

## 4. Setting Up

### 4.1 Unpacking

- ⇒ Remove the device and all the accessories from the delivery package.
- ⇒ Check the delivery for completeness (refer to the “Scope of Delivery” section).

### 4.2 Setup

The extraction unit is a free-standing appliance, intended to be set-up on the floor (e.g. under the work-bench).

It is possible to connect up to two suction points.

Position the extraction device so that:

- The exhaust vent (14, Fig. 1) is not blocked.
- The distance to the suction points is approximately even.
- The front of the device is easily accessible for removal of dust drawer.

If the device is to be installed in a cabinet, an opening for the exhaust air with the following dimensions must be provided, in addition to the openings for the suction hoses:

- Circular opening: Min. 120 mm diameters
- Rectangular opening: Min. 170 x 65 mm

#### 4.2.1 External Exhaust Air Route

An external exhaust air route (see accessories) allows the extracted air to leave the laboratory.

The installation details are supplied with the external exhaust air route.

**When the extraction unit is used in conjunction with an external ventilation system, a significant quantity of air is extracted from the room per hour.**

**This can create negative pressure within the room which, when using an air dependent naked flame fed by gas, liquid or solid fuel, can cause poisonous gasses (e.g. carbon monoxide) to be drawn into the working area.**

**It is therefore essential to ensure that the fresh air supply is sufficient and that the environmental air pressure is maintained, this should then be monitored by an authorized specialist (e.g. a certified Gas Service Engineer).**

### 4.3 Electrical Connection

**Before connecting the device, ensure that the voltage information on the identification plate corresponds with your local power supply.**

**Arrange the conducting parts (plug sockets, plugs and couplings) and install the extension cord so that the protection class is retained.**

⇒ Switch the device off at the on / off switch (3, Fig. 1).

⇒ Insert the mains cable into the plug socket.

### 4.4 Compressed Air Connection

The suction hoses are opened and closed via the pneumatic pinch valves (15, Fig. 3).

For this, compressed air is necessary. On the device, the compressed air tube is sealed.

Connection to the compressed air supply:

⇒ Select the appropriate tube connection piece from the connection set (17, Fig. 1) and attach to the end of the tube.

⇒ Connect to the compressed air.

**Observe the minimal / maximum connection pressure, see technical data!**

**Compressed air should be clean and dry (without condensed water). Wet compressed air can damage the unit.**

## 4.5 Installing the Pinch Valves

In order to prevent damage during transport, the following parts

- pinch valves (15)
- suction pipe (16)
- holder (20)

are not mounted on the device, but will be delivered in a pre-assembled state.

- ⇒ These must be installed onto the device before use:
- ⇒ Unscrew the knurled head screw (19, Fig. 4a).
- ⇒ Push the suction pipe (16) onto the suction support on the extractor and hang the holder (20) onto the holder screws (26, Fig. 4b).
- ⇒ Fix the holder in place with the knurled head screw (19, Fig. 4a).
- ⇒ Connect the compressed air tubes to the pinch valves:
  - Ensure the compressed air tubes (A, B) are connected to the correlating suction pipes (A, B).
  - Push the compressed air tube onto the pipe connection in the pinch valve as far as it will go.
 There will be a slight resistance to overcome when attaching.


Fig. 3

### 4.5.1 Suction Pipe Alteration

**!** During the alteration, make sure that the compressed air tube is not unintentionally ripped away from the pinch valves. If possible, loosen before the alteration (press the ring at the compressed air connection and pull the tube towards you) and tighten again after the work has finished.

- ⇒ Pull the pinch valves (15, Fig. 3) out of the suction pipe (16, Fig. 3).
  - ⇒ Unscrew the knurled head screw (19).
  - ⇒ Loosen the suction pipe and the holder (20) from the holding screws (26) by turning and then at the back, remove from the extractor connection fitting (A).
  - ⇒ Lever the suction pipe out of the holder (B), turn (C), and replace back into a new position in the holder (D).
  - ⇒ Insert the suction pipe into the extractor connection fitting on the extraction device, hang (E) the holder into the holding screws (26).
  - ⇒ Fix the holder with the knurled head screw (19).
  - ⇒ Push the pinch valve into the suction pipe as far as it will go.
- This process is made easier by turning slightly when pushing into place.


Fig. 4a


Fig. 4b

## 4.6 Connection to the Extraction Point

The extraction point is at the extractor connection fitting (15).

**When you turn the pinch valve, please ensure that there are no kinks in the compressed air tubes and that they are not unintentionally ripped off!**

The compressed air tubes are marked with A and B so that they can be easily allocated to the relevant suction hoses.

⇒ Shorten the suction pipes for both extraction points to the required length with a knife.

⇒ Screw the PVC end bushings onto the suction pipes (left handed thread).


Fig. 5

**Caution: Risk of injury!**

**When shortening the suction pipe please ensure to cut at a right angle.**

⇒ Insert the suction hose (18) into the extractor connection fitting (15). Please use a hose connection adapter (see accessories) if the diameter sizes are different, in order to minimize ambient noise and loss of extraction force.

⇒ Attach the suction hose to the extraction point (e.g. Dustex master plus, extractor clamp, etc.).

**Long suction hoses, tight bends and kinks will considerably reduce the extraction force at the extraction point.**

**If different suction hose lengths are used, the extraction force will be lower at the extraction point with the longer hose.**

**Make sure that the suction hoses are the same length and diameter. Excess hose can be cut to size or laid in a large arch.**

**Avoid steep pitches or hanging points along the hose path.**

## 4.7 Electrical Connection

⇒ Electrical appliances can be connected to the coupler sockets (9, Fig. 1) at the back of the unit.

Check the connections are correct - coupler socket (A, B) - suction hose / pinch valve (A, B).

**When connecting electrical appliances to the extraction unit, please ensure that the power consumption of all devices connected to the unit power outlet does not exceed the total maximum power for connected equipment (see chapter 8, "Tech. Data").**

**A two pole plug with protective contact is available for use as an adapter. This adapter may only be put together by an electrician!**

**The protective conductor system must not be interrupted by the adapter!**

## 5. Operation

**Before beginning to work, check that a dust bag has been inserted.**

**If operated without a dust bag, the intended protection against exposure to dust can no longer be guaranteed.**

The extractor unit is operated via the buttons on the key-pad (2, Fig. 1).

### 5.1 Switching the Unit On

The extractor is switched ON and OFF at the On / Off switch (3).

Once switched on (\*):

- it is possible to use the suction from one or both suction tubes (see chapter 5.2).
- Various settings can be programmed (see chapter 5.3).

(\* ) If the self-diagnosis function is activated, this will take place after the unit is switched on, before it is possible to work with the suction (see chapter 5.3.1).


Fig. 6

### 5.1.1 Stand-By

If the unit is not used for a long period of time, the extraction unit switches to stand-by mode. The display darkens.

**Change to the stand-by mode:**

- If no button has been pressed and no suction hose is open for a period of 2 min.

**Leave stand-by mode:**

- Press any button.
  - Switch on any electrical appliance connected to the extractor.
- Once switched on, the suction hose will open immediately and the extractor turbine will activate.

## 5.2 Operation

### 5.2.1 Extraction

When the extractor is switched on it goes to the stand-by mode.

When a suction hose is opened, the suction turbine switches to the last selected extraction force and the correlating pinch valve opens. When the second suction hose is then opened, the second pinch valve also opens.


Fig. 7

Opening / closing of the suction hose is as follows:

- Automatically:
  - The electrical appliance connected to the coupler socket (9, Fig. 1) is activated / deactivated.
- Manually:
  - By pressing a suction hose button (39).

The letters of the opened suction hose will be shown in the display (38).

If the electrical appliance is active, the relevant suction hose CANNOT be closed manually with the suction hose button, instead only by switching the electronic equipment off.

### 5.2.2 Extraction force

The extraction force can be adjusted in three levels; the current extraction force is shown in the display (37).

Set the extraction force:

- ⊕ Button (43), increases extraction force
- ⊖ Button (41), reduces extraction force

**Readjustment of the extraction force:**

The set extraction force applies for both suction hoses. When another suction hose is opened or closed, the power in the suction turbine is automatically readjusted so that the selected extraction force may be achieved in each opened suction hose.


Fig. 8

**To deactivate the extraction force readjustment:**

In the stand-by mode:

⇒ Press ⊕ / ⊖ buttons (43/41) together for 3 seconds.

- ◆ The indicator above the button ⊖ (41) shows “OFF”.

The deactivated extraction force readjustment is shown by the “OFF” indicator above the button ⊖ (41) permanently.

**To activate the extraction force readjustment:**

In the stand-by mode:

⇒ Press ⊕ / ⊖ buttons (43/41) together for 3 seconds.

- ◆ The indicator above the button ⊖ (41) goes “OFF”.


**i** When the extraction force readjustment is deactivated, the extraction power on one suction hose will reduce considerably if a second suction hose is then also activated.

**i** On extraction force level "3" with both hoses open, only a reduced extraction power is available as the suction turbine cannot be increased any further.

**i** The fuller the dust bag, the lower the extraction performance.

## 5.3 Programming

The menu function enables various extractor settings to be preset and saved.

In order to set different parameter, press the menu key (30) repeatedly.

To confirm and save the adjusted settings press the enter button (40).


An audible signal confirms that the settings have been successfully saved.

If you do not wish to save the adjustments, press the menu key again (30) and the programming will be aborted.

If the word "default" is shown in the display (42) press the buttons ⊕ and ⊖ (41, 43) together to revert to the standard factory settings (see table in chapter 6.5).

### 5.3.1 Self Diagnosis

The diagnosis symbol (31) indicates that self-diagnosis is performed:

⇒ Press menu key (30) 1 x .

◆ The diagnosis symbol (31) blinks.

⇒ Press the enter button (40).

◆ The self-diagnosis is performed 1 x.

For devices with a software version up to 1.00.15 (\*), it is possible to set the self-diagnosis automatically after each switch-on.

To enable / disable automatic self-diagnosis:

⇒ Press menu key (30) 1 x .

◆ The diagnosis symbol (31) blinks.

⇒ Switch the self diagnosis ON or OFF using the ⊕ / ⊖ buttons (41, 43).

⇒ Press the enter button (40).

(\*) The software version is shown briefly on the display after switching on.

### 5.3.2 Audible Signal

Various settings and warnings are confirmed or alerted via an audible signal. This signal can be switched on or off. The audible signal symbol (32) shows whether the signal is on or off.


Fig. 11

- ⇒ Press menu key (30) 2 x.
  - ◆ The audible signal symbol (32) blinks.
- ⇒ Switch the self diagnosis ON or OFF using the ⊕ / ⊖ buttons (41, 43).
- ⇒ Press the enter button (40).

### 5.3.3 Dust Bag Full

A full dust bag is indicated by the dust drawer symbol (33).

*The fill level of the dust bag depends on the type of dust accumulated.*


Fig. 12

The sensitivity of the “dust bag full” indication can be changed. This will also change the fill level within the dust bag, or duration until the signal “dust bag full” is shown.

- ⇒ Press menu key (30) 3 x.
  - ◆ The dust drawer symbol (33) blinks.
  - ◆ The set sensitivity level is shown.
- ⇒ Change the level using the / buttons (41, 43):
  - ◆ 1: full
  - ◆ 9: low fill level.
- ⇒ Press enter button (40).

### 5.3.4 Start-Up Threshold


Fig. 13

The start-up threshold, the point at which the connected electrical appliance opens the suction hose, can be individually adjusted for each hose.

- ⇒ Press menu key (30) 4 x.
  - ◆ The start-up threshold symbol (34) blinks.
  - ◆ The letters A and B blink.
- ⇒ Press the suction hose button (39) of the suction which is to be adjusted
  - ◆ The letter of the selected hose is indicated.
  - ◆ The display shows „OFF“.
- ⇒ Switch the electrical appliance off. Units equipped with a stand-by mode should be switched to stand-by (e.g. in the case of a hand piece, only switch the controller on without activating the hand piece).
- ⇒ Press the enter button (40).
  - ◆ The display shows “ON”.


Fig. 14

- ⇒ Switch the electrical appliance on, e.g. activate the hand piece at the speed at which you require the suction hose to open.
- ⇒ Press the enter button (40).
- ⇒ Select the next suction hose which you would like to adjust, or quit this menu by pressing the menu key (30).

### 5.3.5 Shut-Off Delay (3 - 30 sec.)

Shut-off delay = the amount of time between switching the connected electrical appliance off and the closure of the respective suction hose (= pinch valve).


Two types of shut-off delay times can be distinguished, each of which can be set independently of one another:

- Only one of the two suction hoses, A OR B, is open:
  - Set the time, from when the pinch valve in the suction hose closes and the suction turbine shuts off (e.g. for targeted suction of residue dust).
- Both suction hoses, A AND B, are open:
  - Set the duration of time, until the pinch valve in the suction hose closes and the suction power for the other suction hose has readjusted.

**Select a longer amount of time in order to prevent short interruptions (e.g. when working with a hand piece motor) from switching the extractor on and off too frequently, therefore causing the suction power to be continuously be readjusted.**

- ⇒ Press the menu key (30) 5 x.
- ◆ The shut-off delay symbol (35) blinks.

**Set shut-off delay A**

- ◆ The letter A is indicated.
  - ◆ The display shows the current shut-off delay setting for suction hose A.
- ⇒ Set the required shut-off delay with the ⊕ / ⊖ buttons (41, 43).  
 ⇒ Press enter button (40).

**Set shut-off delay B**

- ◆ The letter B is indicated.
  - ◆ The display shows the current shut-off delay setting for suction hose B.
- ⇒ Set the required shut-off delay with the ⊕ / ⊖ buttons (41, 43).  
 ⇒ Press enter button (40).

**Set shut off delay A and B**

- ◆ The letters A and B are indicated.
  - ◆ The display shows the current shut-off delay setting.
- ⇒ Set the required shut-off delay with the ⊕ / ⊖ buttons (41, 43).  
 ⇒ Press enter button (40).

The shut-off delay adjustment can be quit at any time before or after the value is saved, by pressing the menu key (30).

**i** *Whilst setting the shut-off delay, press the ⊕ / ⊖ buttons together to return to the factory set shut-off delay parameter.*

**5.4 Operation with the SILENT flow sensor**

If the extraction unit is operated in combination with a flow rate monitoring unit such as the SILENT flow sensor, the following must be observed:

- ▶ **If the flow rate is less than the defined threshold value for a certain period, this is detected by the SILENT flow sensor, the corresponding indicator begins to flash, and a warning tone is output every 30 seconds.**

**If the flow rate is too low at the suction point being monitored, safe operation is no longer guaranteed and the level of exposure to dust can be harmful to health!**

In this case, the following measures must be taken:

- ⇒ Stop working at all suction points of the affected extraction unit. This applies to the individual suction channels of any multi-workbench extraction unit as well as to all suction points operated via a Y-junction or via extractor switches.
- ⇒ Determine why the flow rate is too low and resolve the issue.

Cause	Measure
Suction level is too low	• Choose a higher suction level
Filter blockage	• Clean the filter.
Cleaning the filter does not result in the expected improvement.	• Replace the fine filter. • For low-dust replacement of the fine filter, observe Chapter 6.3.
Blockage in the suction point (suction mouth)	• Disconnect the suction hose at the suction point and check whether there is a blockage in the suction point; remove where appropriate.
Blockage in the suction line	• Check the suction line at all disconnection points / service openings section by section for blockages.
Leakage in the section of the suction line between the SILENT flow sensor and the extraction unit.	• Check the installation of the SILENT flow sensor in accordance with its instructions.

Determining the cause and measures can be carried out in the order shown.

After every measure, you should switch on the extraction unit again and check whether the minimum flow rate is again exceeded so that correct operation is possible again.

⇒ Work should only resume after the fault has been rectified and the minimum flow rate is exceeded again.

## 6. Cleaning / Maintenance

Inside the appliance there are no parts which require maintenance.

Opening the device, other than for the processes described below, is not permissible!

### 6.1 Replacing the Dust Bag

The set dust bag fill level has been reached:

- ◆ The dust drawer symbol (33, Fig. 2) blinks.
- ◆ An audible signal will be emitted 3 x.
- ◆ The display blinks "Cod01".

Emptying the dust drawer:

- ⇒ Switch the extractor off.
- ⇒ Pull the front panel forward and off (5, Fig. 1).
- ⇒ Pull the dust drawer (7) out towards the front.
- ⇒ Remove the protective film from the sticky tab and seal the dust bag opening adhesively to prevent dust escaping!
- ⇒ Bring the dust bag in the dust drawer to the disposal location.
- ⇒ Dispose of the dust bag properly.
- ⇒ Insert a new dust bag into the dust drawer. Ensure that the new dust bag is properly aligned in the guide grooves (21) and the glue flap faces upwards.


Fig. 16


Fig. 17

**!** Use only original Renfert dust bags (see spare parts list).

- ⇒ Inspect the dust drawer seal for damage and replace it, as required (24).
- ⇒ Insert the dust drawer, making sure it is straight, and push it in as far as it will go. Make sure that the dust drawer is seated correctly on the guide (23, Fig. 16).
- ⇒ Install the front panel (5, Fig. 1) on the bottom hooks then push it into the upper lock bolts (22, Fig. 16), pressing firmly until it locks in place.


Fig. 18

## 6.2 Filters

The extractor contains a 2 level filter system comprising:

- dust bag
- fine particle filter

**NEVER operate the dust extractor without the complete filtration system.**

**Use only the original Renfert filter (see spare parts list).**

### 6.2.1 Replacing the Fine Particle Filter

Examine and replace the fine particle filter (6) on a regular basis, at least once a year. It must always be replaced if the dust drawer symbol (33, Fig. 2) blinks, despite a new dust bag having just been replaced. There is a risk of increased exposure to dust when replacing the fine filter.

That's why the following measures must be taken:

- ⇒ Keep the plastic bag ready that contains the new filter.
- ⇒ Keep an additional dust extraction device or suction tube from another workbench extraction unit ready.
- ⇒ Wear at least an FFP2 mask for respiratory protection.
- ⇒ Pull the upper front panel (4, Fig. 1) forward and off.
- ⇒ Turn the fine particle filter in an anti-clockwise direction to release it.
- ⇒ Pull the fine particle filter straight forwards.
- ⇒ Place the fine filters immediately in the plastic bag and seal the bag tightly to prevent dust escaping.
- ⇒ Ensure proper disposal.
- ⇒ Install a new fine filter and fully insert it over the holder bar (25).
- ⇒ Turn the fine particle filter in a clockwise direction and hand-tighten it in place.
- ⇒ Inspect the seal on the front panel for damage and replace it, as required.
- ⇒ Snap the front panel back in place.
- ⇒ Use the additional dust extraction device to suction up dust that has been released.


Fig. 19

### 6.2.2 Replacing the Exhaust Filter

The exhaust filter (26) primarily traps particles rubbed off the suction turbine graphite brushes. The filter should be replaced annually:

- ⇒ Remove the suction pipe (see chapter 4.5.1).
- ⇒ Press the lock tabs (4 pcs.) on the filter cassette (14) and pull the filter cassette off.
- ⇒ Dispose of the filter and filter cassette properly.
- ⇒ Insert a new filter mat in the filter cassette so that the smooth, compressed side of the filter mat faces down or outwards when the filter is installed.
- ⇒ Replace the new filter cassette with the new filter.
- ⇒ Make sure the filter cassette is properly seated and locked in place on all lock tabs.


Fig. 20

### 6.2.3 Replacing the Electronics Filter

Cooling air for the unit's electronics passes through the electronics filter (11). The type of material trapped in the filter depends on the ambient laboratory conditions.

- ⇒ The filter should be examined annually and replaced as necessary:
- ⇒ Press on the lock tabs (2 pcs.) on the filter cassette (11) and pull the cassette off.
- ⇒ Dispose of the filter and filter cassette properly.
- ⇒ Insert a new filter in the filter cassette.
- ⇒ Re-install the new filter cassette with the new filter.
- ⇒ Make sure the filter cassette is properly seated and locked in place on both lock tabs.


Fig. 21

## 6.2.4 Change Compressed filter

- ⇒ Disconnect the appliance from the compressed air supply.
- ⇒ Roll away the O-ring seal.
- ⇒ Unscrew filter housing (hold up at the filter housing with the screwdriver attached).
- ⇒ Using an Allan key (SW 5 mm) unscrew the sleeve (26).
- ⇒ Install a new piece of non-woven material (28) and sieve (27), making sure to adhere to the correct.
- ⇒ Screw the sleeve back in and hand-tighten.
- ⇒ Screw filter housing together, replace the O-ring seal.


Fig. 22

## 6.3 Self-Diagnosis

During self-diagnosis, various device functions are checked. To carry out a self-diagnosis, see chapter 5.3.1. This process takes approximately 30 seconds.

**A blocked or defect suction hose can only be detected by the self-diagnosis and not when the device is in action.**

If an error is detected during the self-diagnosis, a fault code (see chapter 7.1) will be displayed. Please continue as instructed.

## 6.4 Safety Mechanism

The protection of the suction device is achieved by two device protection switches (10).

If an overload switch is released, it can be reset by pressing the button back in.

**Repeated triggering of a device protection switch is caused by a defect in the suction device. Send the device in for repair!**


Fig. 23

## 6.5 Factory Set Parameter

When the function “reset to factory set parameter” is activated, all individually set information is lost and overwritten by factory set values.

- ⇒ Extraction off
- ⇒ ⊕ and ⊖ Press button (41, 43) and keep pressed.
- ⇒ Extraction on
- ⇒ ⊕ and ⊖ Press button (41, 43) and keep pressed
  - ◆ An audible signal will confirm the factory set parameter is saved and the display will indicate „rES“.


Fig. 24

**Factory set parameter:**

Function / Feature	Adjustment Range	Factory Parameter
Extraction force, one suction hose is open	1 - 3	2
Extraction force, both suction hoses are open	1 - 3	2
Sensitivity of dust bag fill level Recognition	1 - 9	5
Audible Signal	On / Off	On
Automatic start-up threshold	1 - 100 W	8 W
Shut-off delay A or B, one suction hose is open	3 - 30 Sec.	3 Sec.
Shut-off delay A / B, both Suction hoses are open	3 - 30 Sec.	15 Sec.

Only for devices with a software version up to 1.00.15 (see chapter. 5.3.1):

Diagnosis - mode (after switching on)	On / Off	Off
---------------------------------------	----------	-----

## 6.6 Spare Parts

You can find components subject to wear and the spare parts on the spare part list in the internet at [www.renfert.com/p918](http://www.renfert.com/p918).

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

Serial number and date of manufacturing are shown on the type plate of the unit.

## 7. Troubleshooting

Error	Possible Cause	Corrective Action
<b>Extractor stops suddenly for no apparent reason and a error code is indicated in the display.</b>	<ul style="list-style-type: none"> <li>An error has been detected.</li> </ul>	<ul style="list-style-type: none"> <li>See table "Error Codes"</li> <li>If the instructions are of no help to you, or the error code is not listed:                             <ul style="list-style-type: none"> <li>- Make a note of the error code.</li> <li>- Send the device in for repair.</li> <li>- Inform the repair service which error code was noted.</li> </ul> </li> </ul>
<b>Suction power too low.</b>	<ul style="list-style-type: none"> <li>Dust bag is full.</li> <li>Changing the dust bag does not result in the expected improvement.</li> </ul>	<ul style="list-style-type: none"> <li>Change the dust bag.</li> <li>Replace the fine filter.</li> <li>For low-dust replacement of the fine filter, observe Chapter 6.2.</li> </ul>
<b>Dust bag over full.</b>	<ul style="list-style-type: none"> <li>Sensitivity of "dust bag full" signal is too large.</li> </ul>	<ul style="list-style-type: none"> <li>Set a smaller value (See chapter 5.3.3).</li> </ul>
<b>"Dust bag full" indicator goes on, even though the bag is not full.</b>	<ul style="list-style-type: none"> <li>Sensitivity of "dust bag full" signal is too low.</li> </ul>	<ul style="list-style-type: none"> <li>Set a larger value (See chapter 5.3.3).</li> </ul>
<b>Suction hose cannot be closed by pressing the suction hose button.</b>	<ul style="list-style-type: none"> <li>An electrical appliance which is connected to the suction hose is still in operation.</li> </ul>	<ul style="list-style-type: none"> <li>Switch the connected electrical appliance off.</li> </ul>
<b>"Dust bag full" indicator appears shortly after having replaced the dust bag.</b>	<ul style="list-style-type: none"> <li>Fine particle filter is blocked.</li> <li>Suction hose blocked.</li> <li>Cross-section of the dust generating device's hose fitting is too small.</li> </ul>	<ul style="list-style-type: none"> <li>Replace the fine particle filter (Chapter 6.2.1).</li> <li>Remove the obstruction in the suction hose.</li> <li>Reset the sensitivity of the "dust bag full" indicator (Chapter 5.3.3).</li> </ul>
<b>Extractor fails to start when a dust generating device is being used.</b>	<ul style="list-style-type: none"> <li>Extractor is not switched on.</li> <li>Device fuse has blown.</li> <li>Dust generating device is not connected to the extractor socket.</li> <li>Start-up threshold of the suction hose is too high.</li> </ul>	<ul style="list-style-type: none"> <li>Switch the extractor on (Chapter 5.1).</li> <li>Check and reinstall the fuse if necessary (Chapter 6.4).</li> <li>Connect the dust generating device to the coupler socket (Chapter 4.7).</li> <li>Adjust start-up threshold (Chapter 5.3.4).</li> </ul>
<b>Extractor starts, although a connected device is not being used.</b>	<ul style="list-style-type: none"> <li>Start-up threshold of one suction hoses is too low.</li> </ul>	<ul style="list-style-type: none"> <li>Check, which suction hose is open, the relevant letter will be indicated in the display.</li> <li>Adjust the start-up threshold of the suction hose (Chapter 5.3.4).</li> </ul>
<b>Extractor fails to stop when a connected device is switched off.</b>	<ul style="list-style-type: none"> <li>The dust generating device is not plugged in to the coupler socket on the extractor and the suction hose button is switched on.</li> <li>Start-up threshold is too low.</li> <li>Shut-off delay is too long.</li> </ul>	<ul style="list-style-type: none"> <li>Switch the suction hose button on the extractor off, and connect the dust generating device to the coupler socket on the extractor (Chapter 4.7).</li> <li>Adjust start-up threshold (Chapter 5.3.4).</li> <li>Set shorter shut-off delay (Chapter 5.3.5).</li> </ul>
<b>No satisfactory setting for the automatic start-up feature can be adjusted.</b>	<ul style="list-style-type: none"> <li>Some older hand pieces do not generate a sufficiently strong signal for the automatic start-up feature.</li> </ul>	<ul style="list-style-type: none"> <li>Switch to higher revs per minute (e.g. 10,000 rpm.).</li> </ul>

Error	Possible Cause	Corrective Action
The suction switches off and on in an uncontrollable manner when the hand piece is activated.	• The hand piece is being used at a lower speed than the rpm to which the automatic switch-on is set.	• The hand piece must at least be used at the speed at which the automatic switch-on is set.
Panels covering the fine particle filter and the dust drawer are difficult to open.	• Lock bolts are dirty.	• Clean lock bolts and lubricate them slightly if necessary.
The extraction force is not readjusted when a second suction hose is switched on or off.	• The automatic extraction force readjustment is deactivated.	• Activate the automatic extraction force readjustment (See chapter 5.2.2).
The display indicates "OFF" above the button  (41).	• The automatic extraction force readjustment is deactivated.	• Activate the automatic extraction force readjustment if required (See chapter 5.2.2).
There is no audible signal.	• The audible signal is deactivated.	• Activate audible signal (Chapter 5.3.2).
Pinch valves close / open slowly.	• The filter in the compressed air connection is blocked.	• Change the filter in the compressed air connection (see chapter 6.2.4).

## 7.1 Error codes

If the following error code appears during operation or self-diagnosis, please continue as follows:

- ⇒ Switch device off.
- ⇒ Deal with the fault as described in the table.
- ⇒ Switch the device back on.

For all other error codes, switch the unit off and then back on again.

If the error code is repeatedly indicated:

- ⇒ Make a note of the error code.
- ⇒ Contact the repair service.

Error Code	Cause / Meaning	Corrective Action
<b>Cod01</b>	• The set dust bag fill level has been reached.	• Change dust bag (Chapter 6.1).
<b>Err01</b>	<ul style="list-style-type: none"> <li>• Suction flow is blocked.</li> <li>• Insufficient compressed air pressure.</li> </ul>	<ul style="list-style-type: none"> <li>• Check the suction hoses for blockages, examine the dust bag fill level, and change if necessary (See chapter 6.1), examine the fine particle filter for contamination, and change if necessary (see chapter 6.2.1).</li> <li>• If the device is supplied with compressed air from a main laboratory switch together with the air compressor, then at first the compressed air will be insufficient and the self-diagnosis will indicate a fault. Switch the device ON at the On / Off switch (3, Fig. 1).</li> </ul>
<b>Err05</b>	• Overheated electronics.	• Switch the extractor off and allow cooling. If the extractor is located in a closed cabinet, make sure the exhausted air is able to escape (See chapter 4.2).
<b>Err26</b>	<ul style="list-style-type: none"> <li>• Suction hose A or B permanently open.</li> <li>• The fine filter is too clouded.</li> <li>• Suction tube is blocked.</li> <li>• Installation with excessively narrow cross-section in the suction line.</li> </ul>	<ul style="list-style-type: none"> <li>• Make sure the pinch valve on the suction hose is not blocked that the pinch valve membrane is not defect.</li> <li>• Check to ensure that the pinch valve is assembled correctly.</li> <li>• Change fine filter (see chapter 6.2.1).</li> <li>• Check suction channel from suction point to pinch valve for blockage and remove if necessary.</li> <li>• Check the suction tube from the suction point to the pinch valve for cross-section restrictions and correct if possible. If necessary, use tube connector adapter (see accessories).</li> </ul>
<b>Err51</b>	• Error occurred whilst adjusting the start-up threshold.	• Cancel the adjustment of the start-up threshold using the menu key (30, Fig. 2) and enter the required setting afresh as described in Chapter 5.3.4.

## 8. Technical Data

Working Voltage	230 V	120 V
Permissible mains voltage:	220 - 240 V	120 V
Mains frequency:	50 / 60 Hz	60 Hz
Suction turbine power: *)	1400 W	1400 W
max. connecting power: *)	2000 W	360 W
Total connected power: *)	3400 W	1800 W
Mains input fuse:	2 x 15 A (T)	
max. connection pressure:	8 bar [116 psi]	
min. connection pressure:	4.5 bar [65 psi]	
LpA **) (at max. volume flow):	56,9 dB(A)	
Number of suction hoses:	2	
Flow rate, max. ***):	3300 l/min	
Negative pressure, max. ****):	200 hPa [2.9 psi]	
Fine particle filter:		
Filter surface area, approx.:	0.8 m <sup>2</sup> [1240 sq inch]	
Filter quality:	Class M according to EN60335-2-69	
Fill volume, dust bag, approx.:	7.5 l [2 US gal]	
Dimensions (Width x Height x Depth):	224 x 591 x 584 mm [8.8 x 23.3 x 23.0 inch]	
Weight (empty), approx.:	27 kg [59.5 lbs]	
Ø Suction fittings:		
internal:	45 mm [1.77 inch]	
external:	50 mm [1.97 inch]	

- \*) Working voltage power
- \*\*) Sound pressure levels according to EN ISO 11202
- \*\*\*) Free-blowing turbine at nominal turbine voltage
- \*\*\*\*) At nominal turbine voltage

## 9. Warranty

Provided the unit is properly used, **Renfert warrants all components for 3 years.**

**The suction motor has a guarantee of 3 years, with a maximum of duration of 1000 operating hours (motor running time).**

Warranty claims may only be made upon presentation of the original sales receipt from the authorized dealer.

Parts which are subject to natural wear and tear (wear parts) and consumables are excluded from the guarantee. These parts are marked in the spare part list.

The warranty is voided in the case of improper use; failure to observe the operating, cleaning, maintenance and connection instructions; in case of independent repairs or repairs by unauthorized personnel; if spare parts from other manufacturers are employed, or in case of unusual influences or influences not in compliance with the utilization instructions.

Warranty service shall not extend the original warranty.

# 10. Disposal Information

## 10.1 Disposing of Consumables

Full dust bags and filters must be disposed of under compliance with locally applicable regulations. Depending on the material trapped by the filters, protective gear may need to be worn during disposal.

## 10.2 Disposing of the Unit

The unit must be disposed of by an authorized recycling operation. The selected firm must be informed of all possible health hazardous residues in the unit.

### 10.2.1 Disposal Instructions for countries in the EU

To conserve and protect the environment, prevent environmental pollution and improve the recycling of raw materials, the European Commission adopted a directive that requires the manufacturer to accept the return of electrical and electronic units for proper disposal or recycling.

**Within the European Union, units with this symbol should not therefore be disposed of in unsorted domestic waste.**

Please contact your local authorities for more information on proper disposal.

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declare under our sole responsibility, that the product(s)

Name	Part No.
SILENT TS	2921 0050
SILENT TS2	2930 0050
SILENT TC	2935 0000 / 2935 0000 03
SILENT TC2	2936 0000 / 2936 0000 03
SILENT EC2	2937 0000 / 2937 0000 03
SILENT compact	2934 0000 / 2934 0000 03
SILENT compactCAM	2934 2000 / 2934 2000 03 2934 2500 / 2934 2500 03
SILENT powerCAM EC	2939 0000
SILENT powerCAM TC	2938 0000
iVAC silent	36 2934 2000

is (are) in conformity with the relevant regulatory requirements by compliance with the UK designated standards.

#### UK legislation

*Electromagnetic Compatibility Regulations 2016*

*Supply of Machinery (Safety) Regulations 2008*

*The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012*

#### Designated standards

EN 61326-1:2013

EN ISO12100:2010

EN 61010-1:2010/A1:2019/AC:2019-04

EN IEC 63000:2018

Tilo Burgbacher,  
Engineering Director

Hilzingen, 02.11.2022

The Technical Documentation has been retained by Renfert GmbH

**Hochaktuell und ausführlich auf ...  
Up to date and in detail at ...  
Actualisé et détaillé sous ...  
Aggiornato e dettagliato su ...  
La máxima actualidad y detalle en ...  
Актуально и подробно на ...**

**[www.renfert.com](http://www.renfert.com)**

Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany  
Tel.: +49 7731 82 08-0 • Fax: +49 7731 82 08-70  
[www.renfert.com](http://www.renfert.com) • [info@renfert.com](mailto:info@renfert.com)

Renfert USA • 3718 Illinois Avenue • St. Charles IL 60174/USA  
Tel.: +1 6307 62 18 03 • Fax: +1 6307 62 97 87  
[www.renfert.com](http://www.renfert.com) • [info@renfertusa.com](mailto:info@renfertusa.com)  
USA: Free call 800 336 7422

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**Renfert GmbH, Industriegebiet, 78247 Hilzingen / Germany**