making work easy

MT plus

FAQ - frequently asked questions

Made in Germany

FAQ - Deutsch

FAQ - English

FAQ - Français

FAQ - Italiano

FAQ - Español

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FAQ - English

- i** For information on operating elements, see [Fig. 1](#) on page 2 of this document.
- i** For chapter indications, refer to the complete instruction manual on the internet.

Fault	Possible cause	Corrective action
Device fails to start.	<ul style="list-style-type: none">Power supply not established.Mains fuse faulty.Front cover not properly closed.Motor overheated.	<ul style="list-style-type: none">Check the mains power supply.Check the mains fuse.Close the front cover (Fig. 11).Allow the motor to cool down, then restart at the main power switch.
Trimmer disc wobbles.	<ul style="list-style-type: none">Disc is not lying flush on the flange.Disc not tightened sufficiently.	<ul style="list-style-type: none">Reinstall the disc according to the directions. Pay attention that the disc is mounted exactly on the lock pin (Fig. 30).Remove all deposits from the back of the disc and from the stopping face of the flange.Hand-tighten the disc (Fig. 29).
Excessive noise while the trimmer disc is running.	<ul style="list-style-type: none">Disc not adequately tightened.	<ul style="list-style-type: none">Hand-tighten the disc (Fig. 29), or reinstall, as required (see operating instructions, chapter 9.1).
Trimmer table wobbling.	<ul style="list-style-type: none">Trimmer table incorrectly installed.	<ul style="list-style-type: none">Reinstall the table according to the directions (Fig. 40).

Wet Trimming:

Fault	Possible cause	Corrective action
Disc buildup.	<ul style="list-style-type: none"> • Water volume incorrectly adjusted at the device. • Water supply not fully open. • Spray tube plugged. • Inlet strainer plugged. • Inlet water pressure too low. 	<ul style="list-style-type: none"> • Increase the water volume at the trimmer (Fig. 9). • Open the water supply (Fig. 6). • Clean / decalcify the spray tube (see operating instructions, chapter 9.3). • Remove the strainer with pliers and clean / decalcify it. • The diameter of the inlet pipe is too small resp. too many users are connected to a small pipe >>> Check if water pressure is at least 1 bar (ask your plumber). If necessary the water pressure can artificially increased by a domestic water work (electric pump and pressure tank). • Ask your plumber.
Water sprays out of the grinding opening.	<ul style="list-style-type: none"> • Too much water. • Water drainage impeded. 	<ul style="list-style-type: none"> • Reduce the water volume at the trimmer (Fig. 9). • Check/properly lay out drain tube >>> slope (Fig. 3a + 3b).
No rinse water.	<ul style="list-style-type: none"> • Water connection not properly established. • Spray tube soiled / calcified. 	<ul style="list-style-type: none"> • Check the water connection (see operating instructions, chapter 8.2 steps 4 - 6). • Clean/decalcify the spray tube (see operating instructions, chapter 9.3).
Water is dripping out between the front door and the housing.	<ul style="list-style-type: none"> • Gypsum powder or plaster grains on the sealing surface. 	<ul style="list-style-type: none"> • Clean the sealing surface and the sealing ring.

Dry Trimming:

Fault	Possible cause	Corrective action
Disc buildup.	<ul style="list-style-type: none">Models are too wet.	<ul style="list-style-type: none">Clean the disc with a brush under running water or with plaster solvent (Renfert GO-2011 Art. No. 2011-0000).
Excessive dust being blown out of the grinding opening.	<ul style="list-style-type: none">Extraction not switched on.Insufficient suction.Extractor or dust bag full.Extractor hose blocked.	<ul style="list-style-type: none">Switch the extraction on.Increase the suction.Replace the dust bag.Check / clean the hose.

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