making work easy

de

en

fr

it

es

# SIMPLEX model designer

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Made in Germany

21-9053 31102024 / Vers. 1.10

making work easy

# SIMPLEX model designer

# Inhalt

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# 1 Introduction

## 1.1 Symbols used

In these Instructions For Use (IFU) or on the equipment, symbols with the following meanings are used:

### Caution

**Risk of damage to the equipment if the instruction is not followed.**

### Note

**Provides the operator with useful information that facilitates use.**

- ▶ List, special attention should be paid to

- List
- Subordinate list

⇒ Instructions / required action / input / operational sequence:

You are asked to carry out the specified action in the order indicated.

- ◆ Result of an action / equipment response / software response:

The equipment or software responds as a result of your actions or when a specific incident occurs.

Other symbols are explained as they occur.

## 2 Safety

### 2.1 *Intended use*

SIMPLEX model designer is a CAD software for creating digital working models and planning models from intraoral scan (IOS) data for orthodontic use.

Extending the scope of use beyond the above-mentioned indications is outside the intended use of this software and is therefore not permitted.

## 3 Installation

To download the software, you must log in to “myRenfert” and register your SIMPLEX model designer.

The link to a Quick Start Guide, how to log in to myRenfert and how to register the software can be found on the dongle packaging.

After downloading, you must install the software on your computer:

⇒ Open the download directory on your PC.

⇒ Unzip the ZIP file.

⇒ Click on “modeldesigner\_.....\_installer.exe” and follow the instructions.

- ◆ The SIMPLEX model designer software will be installed on your PC.

### 3.1 Unlocking the software

The software is protected with a dongle.

It can only be used when the dongle is plugged into the PC:

⇒ Plug the dongle into a free USB port of your PC.

## 4 User interface

When you start the software, you will see the following user interface

- 1 Menu bar:  
The menu bar contains all program functions.
- 2 Toolbars (can be enabled/disabled):  
For easier and faster access, various functions can also be selected using the icons in the toolbar. Some functions can also be selected using shortcuts.
- 3 Action menu:  
This area indicates additional actions that can be performed with the mouse and keyboard in the current work step.
- 4 Status bar (can be enabled/disabled)


Fig. 1

**Some software functions are not available!**

For the intended use of the software, i. e. creation of planning and working model from oral scan data, several software functions are not required and therefore cannot be called up.

Although they are listed here for the sake of completeness, they are disabled in these Instructions For Use and either grayed out or marked accordingly.

## 5 Software functions / operation

### 5.1 Assignment of mouse buttons

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| Select / highlight | Left mouse button                                   |
| Rotate             | Hold down the right mouse button and move the mouse |
| Move/drag          | Press mouse wheel and move mouse                    |
| Zoom               | Rotate the mouse wheel                              |

While editing data, the currently available mouse/keyboard interactions are displayed in the action menu (top right corner of the screen, 3, Fig. 1).

It is also possible to use the 3D mouse (3Dconnexion Spacemouse) for operating the software.

### 5.2 File menu

| Icon                                                                              | Shortcut | Menu item          | Description                                                                    |
|-----------------------------------------------------------------------------------|----------|--------------------|--------------------------------------------------------------------------------|
|  |          | New Order          | Creates a new order and starts the wizard.                                     |
|  |          | Open Order         | Opens the database with all created orders.                                    |
|                                                                                   |          | Import Geometry *) | This function permits to load implemented standard data such as fixation pins. |
|                                                                                   |          | Print Preview      | Displays the print preview.                                                    |
|                                                                                   | Ctrl+P   | Print              | Prints the current view.                                                       |
|                                                                                   |          | Print Setup        | Changes the printer and the print options.                                     |
|                                                                                   |          | Recent File *)     | Back to previous file                                                          |
|                                                                                   |          | Exit               | Exits the software                                                             |

\*) Not available in this software release.

### 5.3 Edit menu


| Icon                                                                                | Shortcut  | Menu item         | Description                           |
|-------------------------------------------------------------------------------------|-----------|-------------------|---------------------------------------|
|  | Ctrl+Z    | Undo              | Undoes the last step.                 |
|                                                                                     |           | Delete Objects    | Deletes one or more selected objects. |
|  | SHIFT+DEL | Clear All Objects | Deletes all loaded objects.           |

## 5.4 Object menu

| Icon                                                                                | Shortcut | Menu item                | Description                                                      |
|-------------------------------------------------------------------------------------|----------|--------------------------|------------------------------------------------------------------|
|                                                                                     | CTRL+I   | Object Info              | Displays information about the order and any existing objects.   |
|                                                                                     |          | Object Color *)          | Opens a color selection dialog to change the color of an object. |
|                                                                                     |          | Mosaic                   | Displays each object in a different color.                       |
|    |          | Points                   | Various display modes for wireframes.                            |
|    |          | Wireframe                |                                                                  |
|    |          | Flat-Shaded              |                                                                  |
|    |          | Gouraud-Shaded           |                                                                  |
|    |          | Transparent              |                                                                  |
|    |          | Flat + Wireframe         |                                                                  |
|                                                                                     |          | Show Objects *)          | Shows selected objects again                                     |
|                                                                                     |          | Hide Objects *)          | Hides selected objects                                           |
|                                                                                     | S        | Toggle Show/Hide *)      | Toggles between visible and hidden objects                       |
|                                                                                     | W        | Toggle Wax Layer         | Show or hide the selected item.                                  |
|                                                                                     | P        | Toggle Occlusal Plane    | Show or hide the selected item.                                  |
|                                                                                     | R        | Toggle Median Plane      | Show or hide the selected item.                                  |
|                                                                                     | B        | Toggle Tuber Plane       | Show or hide the selected item.                                  |
|                                                                                     | S        | Toggle Socket            | Show or hide the selected item.                                  |
|                                                                                     | G        | Toggle Gingiva Mask      | Show or hide the selected item.                                  |
|                                                                                     | Z        | Toggle Section Plane     | Show or hide the selected item.                                  |
|                                                                                     | Y        | Toggle Fixation          | Show or hide the selected item.                                  |
|                                                                                     | T        | Toggle Articulator *)    | Show or hide the selected item.                                  |
|                                                                                     | A        | Toggle Annotation        | Show or hide the selected item.                                  |
|  | U        | Toggle Upper Jaw         | Show or hide the selected item.                                  |
|  | L        | Toggle Lower Jaw         | Show or hide the selected item.                                  |
|  | O        | Toggle Pre-op Scan       | Show or hide the selected item.                                  |
|                                                                                     | X        | Toggle Upper / Lower Jaw | Changes the view.                                                |

\*) Not available in this software release.

## 5.5 Analysis menu

| Icon                                                                              | Shortcut | Menu item          | Description                                                                                                                                                                                                                                                                                  |
|-----------------------------------------------------------------------------------|----------|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  |          | Point Distance 3D  | Starts "Point Distance 3D", a function that allows you to select points on the user interface by left-clicking the mouse and to display their 3D distance to each other.                                                                                                                     |
|                                                                                   |          | Point Distance 2D  | Starts "Point Distance 2D", a function that allows you to select points on the user interface by left-clicking the mouse and to display their 2D distance to each other. As you move the mouse, the 2D distance between this point and the current mouse position is displayed in real time. |
|  | Ctrl+O   | Occlusion          | Show or hide the selected item.                                                                                                                                                                                                                                                              |
|  |          | Material Thickness | Another useful tool for analyzing the created model is to output the material thickness. The different material thicknesses are indicated in different colors. Move the mouse over the user interface to display the exact values at the respective points.                                  |

## 5.6 Display menu

| Icon                                                                                | Shortcut | Menu item                                       | Description                                                                        |
|-------------------------------------------------------------------------------------|----------|-------------------------------------------------|------------------------------------------------------------------------------------|
|                                                                                     |          | Rotation<br>X-Axis<br>Y- Axis<br>Z-Axis<br>Free | Sets the rotation parameters based on individual axes or without any restrictions. |
|                                                                                     |          | Predefined Views                                | Sets the view to one of seven predefined perspectives.                             |
|    | Ctrl+F2  | Top                                             |                                                                                    |
|    | Ctrl+F3  | Bottom                                          |                                                                                    |
|    | Ctrl+F4  | Left                                            |                                                                                    |
|    | Ctrl+F5  | Right                                           |                                                                                    |
|    | Ctrl+F6  | Front                                           |                                                                                    |
|    | Ctrl+F7  | Back                                            |                                                                                    |
|    | Ctrl+F8  | Isometric                                       |                                                                                    |
|   |          | Reset View                                      | Resets the 3D Viewer to default (centered from top).                               |
|  | F6       | Center View                                     | Display models in the center of the screen.                                        |
|                                                                                     |          | Ignore Colors                                   | Ignores the automatic detection of color markers when defining the borders.        |
|  | Ctrl+G   | Draw Grid                                       | Show or hide the selected item.                                                    |
|                                                                                     |          | Set Rotation Point                              | Sets the rotation point on the 3D model for future view rotations.                 |
|  |          | Rectangular Zoom                                | Set the zoom range by moving the mouse while holding down the left mouse button.   |
|                                                                                     |          | Store Viewport                                  | Saves the current view parameters for later use.                                   |
|                                                                                     |          | Restore Viewport                                | Sets the view according to the previously saved parameters.                        |
|                                                                                     | Ctrl+C   | Snap Screen                                     | Copies screenshot to clipboard.                                                    |

## 5.7 View menu

| Icon | Shortcut | Menu item           | Description                                                                                            |
|------|----------|---------------------|--------------------------------------------------------------------------------------------------------|
|      |          | Toggle "Toolbars"   | Show or hide the selected item.                                                                        |
|      |          | Toggle "Status bar" | Show or hide the selected item.                                                                        |
|      |          | Customizing         | Customer-specific adaptations can be loaded here, e. g. additional desired geometry shapes / adapters. |
|      |          | Model types *)      | Default values for designing a model type can be entered here.                                         |
|      |          | Options (Settings)  | Displays general settings / options.                                                                   |

\*) Not available in this software release.

## 5.7.1 Model Types

Create and define model types that can be selected when creating a new order.

The data of a model type stored here will be displayed in different steps as default settings.

## 5.7.2 Options (Settings)

- ▶ **Reendering:**  
Various settings for optimizing screenshots.
- ▶ **Database Folder:**  
Folder where orders are stored.  
If a new folder is selected, the previous jobs can be copied to the new folder by clicking “Copy database”.
- ▶ **Other folders:**  
Folders for importing scan data and exporting model data, as well as the folder where the SIMPLEX sliceware software is installed.
- ▶ **Language:**  
Setting of the selected interface language. The language will not be changed until the software is restarted.
- ▶ **Construction:**  
Various settings for model design and data import and export.
  - **Remember Dimensions:**  
The last manual settings will be used for the next design.
  - **Reset Dimensions:**  
Resets the settings to default values.
  - **Automatic Saving:**  
Data is saved automatically when you click “Next” in the wizard or at regular intervals.
  - **Basal Chamfer:**  
Basal chamfering of the dental arch.
  - **Lower Fixation:**  
For working models, the created articulator supports are automatically extended to the bottom of the model.
  - **Import Reduction:**  
Reduces the scan data based on the selected accuracy value.
  - **Export Reduction:**  
Reduces the model data based on the selected accuracy value.
  - **Minimum Height:**  
Extends the scan data to this value. If set to zero, the bottom of the model is set to the lowest scan point.
  - **Customize Fixation hole:**  
The hole diameter for the fixation pins can be increased for specific printers.
  - **Occlusal Plane:**  
Defines the jaw data set on which the occlusal plane is created.
- ▶ **Result Rotation:**  
Adjustable orientation of the design after clicking the “Export Model” button in the “Finalize” area for correct transmission to the 3D printer.


Fig. 2

## 6 Planning Model – Design sequence

A wizard guides the user through the process of creating a planning model in 5 steps.

It is only possible to move to the next step after completing the previous step with all required entries and actions.

Before that, the “Next” button is disabled.

A planning model is created in the following 5 steps:

1. Order
2. Import
3. Direction
4. Model
5. Finalize

The wizard is started via “File – New Order” or “File – Open Order”

⇒ File – New Order:

- ◆ A new order is created.

⇒ File – Open Order:

- ◆ This function permits to continue editing an existing order. A selection window will open showing the orders that have already been created. Select the desired order and load it with “OK”.

**A reloaded order will always be opened in the step in which it was closed in the previous session.**

**To make changes again, you can go back to the previous step from any step using the “Back” button.**

### 6.1 Planning model – Order

Enter the patient data and select whether to create a working or a planning model.

- 1 Patient: The patient's name (optional).
- 2 Type: Selection of previously created model types including their parameters, see also chapter 5.7.1.
- 3 Technician: Name of the technician (optional).
- 4 Date: Assigned by the system, cannot be changed.
- 5 Comment: Enter additional, individual details about the order (optional).
- 6 Task: Select whether to create a Planning Model or a Working Model.

⇒ Select “Planning Model”!

- 7 Finalize this step and open the next step.
- 8 Cache your entries.

Fig. 3

### 6.2 Planning Model – Import

Load the scan data. You can load the following file formats:

- STL
  - PLY with color information
- 1 Import Upper Jaw (maxilla).
  - 2 Import Lower Jaw (mandible).
  - 3 Import Pre-Op Scan (optional).
  - 4 Trim Scan (refer to 6.2.1).
  - 5 Show Penetration (refer to 6.2.2).
  - 6 Finalize this step and open the next step.
  - 7 Cache your entries.

Fig. 4

**“Show Penetration” can only be selected if maxilla and mandible are loaded.**

### 6.2.1 Trim Scan

Click the left mouse button to drag out a box (1) and then delete the contents.

The available mouse and keyboard actions are displayed in the action menu (2).


Fig. 5

### 6.2.2 Show Penetration

**i** This function is only available if both, maxilla and mandible, were loaded.

This view shows the premature contacts in static occlusion.

⇒ Press the Enter key.

- ◆ A penetration log is created as a PDF file.

⇒ Press the ESC key.

- ◆ The view is closed.


Fig. 6

### 6.3 Planning Model – Direction

- 1 Define Occlusal Plane (refer to 6.3.1).
- 2 Define Median Plane (refer to 6.3.2).
- 3 Define Tuber Plane (refer to 6.3.3).
- 4 Finalize this step and open the next step.
- 5 Cache your entries.

**i** The definition of the planes must then be confirmed by pressing the “Enter” key.


Fig. 7

#### 6.3.1 Define Occlusal Plane

The mouse cursor changes to a crosshair.

If maxilla and mandible are loaded, the jaw defined in the occlusal plane settings will be displayed.

The three highest points (usually in the middle of the incisors and on the large molars) must be selected by clicking on them with the left mouse button.

These three points are used to define the occlusal plane.


Fig. 8

The occlusal plane can be selected by moving the arrow in the centre (1), as well as each of the three points (2).

The action menu (3) indicates which actions are possible with the mouse and keyboard.

## 6.3.2 Define Median Plane

The mouse cursor changes to a crosshair.

Left-click on the first median point according to the model analysis (e. g., papilla incisiva).

Move the mouse to draw a straight line from this point. The median plane is determined by clicking on another reference point (e. g. palatal suture).

The median plane can be precisely adjusted by clicking and dragging the corner points.


Fig. 9

## 6.3.3 Define Tuber Plane

The tuber plane is automatically created perpendicular to the occlusal plane and the median plane.

The exact position can be changed by rotating the mouse wheel.

When the tuber plane is in the desired position, press Enter to confirm all planes.


Fig. 10

**If the occlusal plane is subsequently changed, the median and tuber planes are deleted and must be repositioned.**

**If the median plane is subsequently changed, the tuber plane is deleted and must be recreated.**

## 6.4 Planning Model – Model

- 1 Socket (Model base) / Form (Shape).
- 2 Finalize this step and open the next step.
- 3 Cache your entries.

Check that the selected model base size matches the model.

If the model base is too small or the if model is placed outside the base, an error message is displayed.

The selected model bases are displayed on the model.

The position of the model bases can be changed.

⇒ Move the mouse cursor to one of the model bases:

- ◆ Forward/backward arrows appear on screen.
- ◆ Use the mouse wheel to move both model bases forwards / backwards.

⇒ Move the mouse cursor to one of the red control points of the forward/back arrows:

- ◆ The arrow changes to an up/down arrow.
- ◆ This jaw can be moved up/down with the mouse wheel.


Fig. 11


Fig. 12

The "...variable" model base types can be scaled by pressing the mouse button and Ctrl key.

The yellow line on the jaw model defines the boundary line for the model base shape.

This line can be moved by clicking on individual points and dragging the mouse.

If a point is activated with the mouse and the mouse wheel is rotated, all points move up/down.

## 6.5 Planning Model – Finalize

The models are calculated with the entered parameters. Various adjustments can now be made.

**The model to which the function is to be applied can be determined by selecting “Upper Jaw / Lower Jaw” (maxilla/mandible).**

- 1 Engrave Model (refer to 6.5.1).
- 2 Finalize Model (refer to 6.5.2).
- 3 Export Model (refer to 6.5.3).
- 4 Open Sliceware (refer to 6.5.4).
- 5 Upper Jaw / Lower Jaw (maxilla/mandible): Select the model to which the following actions will be applied.
- 6 “Finish” closes the wizard.
- 7 Cache your entries.


Fig. 13

### 6.5.1 Engrave Model

A menu window opens in which the following entries can be made:

- Label (Text):  
Enter the desired engraving text; the patient's name is automatically suggested.
- Sunken (Embossed):  
Select whether the lettering should be embossed or raised.
- Depth:  
The depth of the embossment.

The mouse cursor changes to a crosshair.

The action menu (3, Fig. 1) indicates which actions can be performed with the mouse.

Close the menu window with “OK” to execute the action.

**If you select “Engrave Model” multiple times, multiple different engravings will be applied.**


Fig. 14

### 6.5.2 Finalize Model

A menu window opens with the following tabs:

- Freeforming
- Constructing

The mouse cursor changes to a selection tool.

This function permits to add or remove material to/from the model using several methods.

The action menu (3, Fig. 1) indicates which actions can be performed with the mouse.


Fig. 15

### 6.5.3 Export Model

Click this button to open the Windows Explorer.

If a path has been defined under “Settings – Export Folder”, it will be opened.

You can also navigate to the desired export folder and save the data records directly there as open STL files.

### 6.5.4 Open sliceware

Click this button to start the SIMPLEX sliceware software and to automatically transfer the model to the software.

# 7 Working Model – Design sequence

A working model is created in 6 steps:

1. Order
2. Import
3. Direction
4. Block out
5. Socket (Model base)
6. Finalize

The wizard is started via “File – New Order” or “File – Open Order”

⇒ File – New Order:

- ◆ A new order is created.

⇒ File – Open Order:

- ◆ This function permits to continue editing an existing order. A selection window will open showing the orders that have already been created. Select the desired order and load it with “OK”.

**A reloaded order will always be opened in the step in which it was closed in the previous session.**

## 7.1 Working Model – Order

When selecting whether a planning model or working model is to be created, select “Working model”!

Enter the patient data and select whether to create a working or a planning model.

- 1 Patient: The patient's name (optional).
- 2 Type: Selection of previously created model types including their parameters, see also chapter 5.7.1. \*).
- 3 Technician: Name of the technician (optional).
- 4 Date: Assigned by the system, cannot be changed.
- 5 Comment: Enter additional, individual details about the order (optional).
- 6 Task: Select whether to create a planning model or a working model.

⇒ Select “Working Model”!

7 Finalize this step and open the next step.

8 Cache your entries.

Fig. 16

## 7.2 Working model – Import

Load the scan data. You can load the following file formats:

- STL
- PLY with color information

- 1 Import Upper Jaw (maxilla).
- 2 Import Lower Jaw (mandible).
- 3 Import Pre-Op Scan (optional).
- 4 Trim Scan (refer to 7.2.1).
- 5 Show Penetration (refer to 7.2.2).
- 6 Finalize this step and open the next step.
- 7 Cache your entries.

Fig. 17

### 7.2.1 Trim Scan

Click the left mouse button to drag out a box (1) and then delete the contents.

The available mouse and keyboard actions are displayed in the action menu (2).


Fig. 18

### 7.2.2 Show Penetration

**i** This function is only available if both, maxilla and mandible, were loaded.

This view shows the premature contacts in static occlusion.

⇒ Press the Enter key.

- ◆ A penetration log is created as a PDF file.

⇒ Press the ESC key.

- ◆ The view is closed.


Fig. 19

### 7.3 Working Model – Direction

- 1 Define Occlusal Plane (refer to 7.3.1).
- 2 Define Median Plane (refer to 7.3.2).
- 3 Block Out (refer to 7.3.3).
- 4 Upper Jaw / Lower Jaw (maxilla/mandible): Select the model to which the following actions will be applied.
- 5 Finalize this step and open the next step.
- 6 Cache your entries.

**i** The definition of the planes must then be confirmed by pressing the Enter key.


Fig. 20

#### 7.3.1 Define Occlusal Plane

The mouse cursor changes to a crosshair.

If maxilla and mandible are loaded, the jaw defined in the occlusal plane settings will be displayed.

The three highest points (usually in the middle of the incisors and on the large molars) must be selected by clicking on them with the left mouse button.

These three points are used to define the occlusal plane.


Fig. 21

The occlusal plane can be selected by moving the arrow in the centre (1), as well as each of the three points (2).

The action menu (3) indicates which actions are possible with the mouse and keyboard.

### 7.3.2 Define Median Plane

The mouse cursor changes to a crosshair.

Left-click on the first median point according to the model analysis (e.g., papilla incisiva).

Move the mouse to draw a straight line from this point. The median plane is determined by clicking on another reference point (e.g. palatal suture).

The median plane can be precisely adjusted by clicking and dragging the corner points.


Fig. 22

### 7.3.3 Insertion Direction

Blocking out in the insertion direction is possible by clicking on the "Insertion Direction" button.

This divides the software window into three different areas. By rotating the model in the large left-hand window (hold down the mouse cursor and left mouse button), the undercuts change according to the viewing direction.

The changed undercut areas are visualized in the small windows on the right.

The views in the right windows can be freely defined by moving the mouse cursor and holding down the left mouse button without changing the selected insertion direction. If the mouse cursor is moved over undercut areas, a numerical value appears which describes the depth of the undercut at this point.


Fig. 23

⇒ Enter key:

- ◆ The newly set insertion direction is accepted.

⇒ ESC key:

- ◆ Cancels the selection.

⇒ DEL key:

- ◆ Deletes the insertion direction already selected.

“Angle” can be used to define an angle for the block-out wax in relation to the insertion direction.

## 7.4 Working Model – Block-out

If an insertion direction has been entered, all undercut areas are automatically blocked out.

- 1 Alter Wax Layer (refer to 7.4.1).
- 2 Upper Jaw / Lower Jaw (maxilla/mandible): Selection of the model.
- 3 Finalize this step and open the next step.
- 4 Cache your entries.

The model to which the function is to be applied can be determined by selecting “Upper Jaw / Lower Jaw” (maxilla/mandible).


Fig. 24

### 7.4.1 Alter Wax Layer

Additional application, removal and smoothing of block-out wax.  
The mouse cursor changes to a selection tool.


Fig. 25

The size can be adjusted by holding down the Shift key and rotating the mouse wheel.  
The displayed parameters are applied to the selection tool based on the selected action (apply material, remove material, smooth surface).

#### Intensity / Thickness

**Intensity:** Can be set as a percentage (%) and reacts to the duration of the mouse click.

**Thickness:** Can be adjusted during application and removal, allowing for uniform application/removal in mm.

### 7.5 Working Model – Socket (Model Base)

You can now define base plates and, if required, articulator supports.

- 1 Model base ("Socket"; refer to 7.5.1).  
- Fixation \*)
- 2 Upper Jaw / Lower Jaw (maxilla/mandible): Selection of the model.
- 3 Finalize this step and open the next step.
- 4 Save your entries.

**i** The model to which the function is to be applied can be determined by selecting "Upper Jaw / Lower Jaw" (maxilla/mandible).

\*) Not available in this software release.


Fig. 26

## 7.5.1 Socket (Model base)

### Form (Shape):

Various model base shapes are available. They can be selected via the drop-down menu.

### Thickness:

Determines the thickness of the base plate.

If the mouse cursor is moved over the base plate, the color changes and the object is enabled.

If you hold down the left mouse button, you can move the base plate in XY direction.

Rotate the mouse wheel to scale the base plate.

Move the mouse to the end of the base plate to enable a crosshair.

Move the mouse cursor to the corresponding axis and rotate the mouse wheel to scale the base plates separately in XY extension or to adjust the height.


Fig. 27

## 7.6 Working Model – Finalize

The models are calculated with the entered parameters. Various adjustments can now be made.

**i** The model to which the function is to be applied can be determined by selecting “Upper Jaw / Lower Jaw” (maxilla/mandible).

- 1 Engrave Model (refer to 6.5.1).
- 2 Finalize Model (refer to 6.5.2).
- 3 Export Models (refer to 6.5.3).
- 4 Open Sliceware (refer to 6.5.4).
- 5 Upper Jaw / Lower Jaw (maxilla/mandible): Selection of the model.
- 6 Save your entries.
- 7 “Finish” closes the wizard.


Fig. 28

### 7.6.1 Engrave Model

A menu window opens in which the following entries can be made:

- Label (Text):  
Enter the desired engraving text; the patient's name is automatically suggested.
- Sunken (Embossed):  
Select whether the lettering should be embossed or raised.
- Depth:  
The depth of the embossment.

The mouse cursor changes to a crosshair.

The action menu (3, Fig. 1) indicates which actions can be performed with the mouse.

Close the menu window with “OK” to execute the action.

**i** If you select “Engrave Model” multiple times, multiple different engravings will be applied.


Fig. 29

## 7.6.2 Finalize Model

A menu window opens with the following tabs:

- Freeforming
- Constructing

The mouse cursor changes to a selection tool.

This function permits to add or remove material to/from the model using several methods.

The action menu (3, Fig. 1) indicates which actions can be performed with the mouse.


Fig. 30

## 7.6.3 Export Model

Click this button to open the Windows Explorer.

If a path has been defined under “Settings – Export Folder”, it will be opened.

You can also navigate to the desired export folder and save the data records directly there as open STL files.

## 7.6.4 Open sliceware

Click this button to start the SIMPLEX sliceware software and to automatically transfer the model to the SIMPLEX sliceware software.

**Hochaktuell und ausführlich auf ...  
Up to date and in detail at ...  
Actualisé et détaillé sous ...  
Aggiornato e dettagliato su ...  
La máxima actualidad y detalle en ...  
Актуально и подробно на ...**

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Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany