making work easy







- DE
- EN
- FR
- IT
- ES
- PT
- TR
- RU
- PL
- CS
- ZH
- JA
- KO
- AR
- UK






Millo / Millo pro

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Made in Germany

21-6543 21112022

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- 1 Introduction..... 2
- 2 Symbols..... 2
- 3 Information for the purchaser..... 2
- 4 Application..... 2
 - 4.1 Correct use of the unit..... 2
 - 4.2 Unapproved use..... 2
 - 4.3 Ambient conditions (according to DIN EN 61010-1)..... 3
- 5 Hazards and warnings..... 3
- 6 Approved personnel..... 4
- 7 Disclaimer..... 4
- 8 Package contents..... 4
- 9 Models..... 4
- 10 Accessories..... 4
- 11 Prior to installation..... 4
 - 11.1 Connecting the extractor..... 5
- 12 Setting up and installation..... 5
 - 12.1 Setting up..... 5
 - 12.2 Electrical connection..... 5
 - 12.3 Absaugung anschließen..... 5
- 13 Operation..... 5
 - 13.1 Controls..... 5
 - 13.2 Switching on / off..... 6
 - 13.3 Operating instructions..... 6
 - 13.4 Adjusting the height of the cutter..... 6
 - 13.5 Changing the cutter..... 6
- 14 Cleaning / Servicing..... 7
 - 14.1 Cleaning the housing..... 7
 - 14.2 Cleaning the work table..... 7
 - 14.3 Changing the fuse..... 7
- 15 Replacement parts..... 8
- 16 Repair..... 8
- 17 Troubleshooting..... 8
- 18 Technical data..... 9
- 19 Guarantee..... 9
- 20 Disposal instructions..... 9
 - 20.1 Disposal instructions for countries in the EU..... 9

1 Introduction

Thank you for purchasing the Millo or Millo pro dental arch trimmer. This unit sets new standards in terms of functionality, efficiency and ergonomics.

Please read the following operating instructions carefully and adhere to the safety instructions to ensure long-term, problem-free operation of the unit.

2 Symbols

The following symbols contained in these instructions or on the unit indicate:

Hazard
Immediate risk of injury.

Dangerous voltage
Risk of electric shock.

Caution
Risk of damage to the unit if the instruction is not followed.

Information
Useful information to facilitate operation and handling.

The device complies with the requirements of the applicable EU directives.

This product complies with the relevant UK legislation.
See UKCA Declaration of Conformity in the Internet under www.renfert.com.

The device is subject to the EU directive 2002/96/EG (WEEE directive).

Caution!
Wear protective glasses.

Caution!
Wear a face mask.

Only to be used indoors.

Unplug the unit before opening it.

Other symbols are explained as required.

3 Information for the purchaser

Instruct operators on the application, possible hazards during use and operation of the Millo / Millo pro arch trimmer based on the user information.

Ensure that this user information is available to operators.

4 Application

4.1 Correct use of the unit

Only use the unit for trimming dental stone or resin dental arches

4.2 Unapproved use

Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the device, increase the risk of serious injury and lead to damage to the environment or the device itself.

4.3 Ambient conditions (according to DIN EN 61010-1)

The unit should only be operated:

- indoors,
- up to 2,000 metres above sea level,
- at an ambient temperature of 5 °C - 40 °C [41 °F - 104 °F] *),
- at a maximum relative humidity of 80 % at 31 °C [87.8 °F], linear reduction up to 50 % relative humidity at 40 °C [104 °F] *),
- with mains electricity supply provided that the voltage fluctuation is within 10 % of the rated value,
- with Pollution Degree 2,
- with Overvoltage Category II.

*) The unit can be used at a temperature of 5 °C - 30 °C [41 °F - 86 °F] and at a humidity of up to 80 %. At temperatures of 31 °C - 40 °C [87.8 °F - 104 °F] the humidity must reduce proportionately to ensure that the unit can be operated (e.g. at 35 °C [95 °F] = 65 % humidity, at 40 °C [104 °F] = 50 % humidity). The unit should not be operated at temperatures above 40 °C [104 °F].

5 Hazards and warnings

Granted protection is not given any longer in case the unit is not operated according to the instruction manual on hand.

Only to be used indoors. The unit is only intended for dry operation and should not be used or stored outdoors or in wet conditions.

The unit should only be operated using a mains cable with a country-specific plug system. Any modification required should only be carried out by an electrician.

The unit should only be operated if the data on the rating plate corresponds with the data of regional mains voltage.

The unit should only be plugged into earthed sockets.

The mains socket should be easily accessible.

Regularly check connecting cables and hoses (e.g. mains cable) for damage (e.g. kinks, tears, porosity) or deterioration. Do not operate units with damaged mains wires, hoses or other defects.

The unit should always be operated under supervision.

The plug-in socket is only suitable for connecting an extractor. The recommended connected load should not be exceeded.

Risk of injury!

Risk of injury if the components used are not approved. Always use original Renfert components.

Wear protective glasses!

Wear a face mask!

Depending on the work piece involved, the operator's work behavior (e.g., application pressure), the setup location and - in particular - the duration of operation, the use of this device can result in hearing damage.

The operator is obligated to have the workplace related assessment level checked in accordance with local requirements and to institute appropriate hearing protection measures.

Wear appropriate personal protective equipment when working with materials that are hazardous to health.

Hazardous dust or fumes may be created when trimming with the unit. An extractor that is appropriate for the materials released should be connected.

Switch off the unit when trimming is complete.

Unplug the unit before carrying out repairs or maintenance on the electrical components.

Switch off and unplug the unit before cleaning or servicing it.

It is the responsibility of the operator that national regulations during operation and regarding a repeated safety inspection of electrical equipment are complied with. For Germany these are the regulation 3 by DGUV (German Statutory Accident Insurance) in relation with VDE 0701-0702 (Association for Electrical & Electronic Technology).

Information on REACH and SVHC is available on our website www.renfert.com, in the Support area.

6 Approved personnel

Operation and servicing of the Millo / Millo pro dental arch trimmer should only be carried out by trained personnel.

7 Disclaimer

Renfert GmbH is not liable for claims for compensation or claims under guarantee if:

- the product is used for purposes other than those stated in the operating instructions.
- the product is modified in any way – apart from modifications described in the operating instructions.
- the product has not been repaired by a specialist firm or original Renfert replacement parts have not been used.
- there is continued use of the product despite obvious safety defects or damage.
- the product has been subjected to mechanical knocks or has been dropped.

8 Package contents

- 1 Millo / Millo pro dental arch trimmer
- 1 Cutter, tapered 6-blades (coarse), NOT suitable for plastic models.
- 1 Allen key
- 1 Operating instructions
- 1 Grounded wall plug (No. 1805-0000 only)

9 Models

1804-0000	Millo, 230 V, 50 / 60 Hz
1804-1000	Millo, 120 V, 60 Hz
1805-0000	Millo pro, 230 V, 50 / 60 Hz
1805-1000	Millo pro, 120 V, 60 Hz

10 Accessories

No. 1806-0001	Cutter, conical cross-cut (fine)
No. 1806-0002	Cutter, tapered 6-blades (coarse), NOT suitable for plastic models.
No. 2921-0050	Silent TS, 230 V, 50 / 60 Hz
No. 2921-1050	Silent TS, 100 - 120 V, 50 / 60 Hz
No. 90003-4305	Universal hose connector adapter
No. 2926-0000	Y-junction
No. 15-0823	PVC Flex super elastic suction hose (sold by the metre)

11 Prior to installation

Compare the data on the rating plate with the data of the regional mains voltage before installation.

The unit should only be plugged into earthed sockets.

11.1 Connecting the extractor

A corresponding connector (Figure 12) is included in the standard delivery of 230 V units with outlets in accordance with DIN 49441 (Figure 2, E). This connector can be used to make an adapter for your local plug configuration.

**This adapter may only be made by a qualified electrician!
The adapter must not interrupt the ground circuit system!**

12 Setting up and installation

12.1 Setting up

- Place the unit on a stable, flat workbench.
- Position the unit so that the mains plug is easily accessible.

12.2 Electrical connection

- Ensure that the voltage corresponds to the voltage on the rating plate of the unit.
- Plug in the unit.

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The Millo / Millo pro must be operated with an extractor.

The minimum air flow volume of the extractor must be 2040 l/h.

Attach the extractor hose to the connector on the unit (Fig. 6).

12.3.1 Option Millo pro

The Millo pro has a plug-in socket.

This plug-in socket has an electrical current when the Millo pro is switched on so that an extractor can be operated directly with it.

Do not exceed the connected load for the plug-in socket.

**The plug-in socket should only be used for connecting an extractor.
Connecting other electrical units to the plug-in socket could damage the Millo pro.**

- Plug the extractor into the plug-in socket on the back of the Millo pro (Fig. 7).

Section 11.1 at the end of these instructions offers information concerning electrical systems employing different power plug designs.

The Millo / Millo pro is now ready for operation.

13 Operation

13.1 Controls

Millo pro (Fig. 1, 2) / Millo (Fig. 3, 4, 5)

- A ON / OFF button
- B ON / OFF LED display
- C Height adjustment (only Millo pro)
- D Connector for extractor
- E Plug-in socket for connecting an extractor (only Millo pro)
- F Safety cut-out (Millo pro, Fig. 2)
- G Fuse links (Millo, Fig. 5)
- H Mains cable

13.2 Switching on / off

Switching on the unit:

- Press the ON / OFF button (A, Fig. 1 & 3):
 - LED display lights up;
 - cutter rotates.
 - Only Millo pro:
 - Electricity is supplied to the plug-in socket and the connected extractor, which is already switched on, starts.

Switching off the unit:

- Press the ON / OFF button (A, Fig. 1 & 3):
 - LED display goes out;
 - cutter comes to a stop.
 - Only Millo pro:
 - The electricity supply to the plug-in socket is switched off and the connected extractor comes to a stop.

13.3 Operating instructions

**Sections of the dental arch may break off during trimming.
Wear protective glasses!**

**The dust created during trimming can be hazardous to health.
Wear a face mask!**

**Risk of injury to hands!
The cutters of the trimmer are sharp.**

**Risk of injury!
When working with the tapered coarse cutter on plastic dental arches risk of injury by hooking into and flying around.
Use the tapered cross-cut cutter for plastic dental arches.**

**Depending on the pressure applied, sound pressure levels of over 85 db(A) can occur.
Wear ear protection!**

**Risk of injury by becoming entangled in the unit!
Do not operate with unprotected long hair or loose clothing.**

- Control the dental arch on the work table using both hands.
- Always keep hold of the dental arch when it is in contact with the cutter.
- Move the dental arch counter to the rotational direction of the cutter:
 - optimal removal
 - reduces the risk of the cutter catching.

13.4 Adjusting the height of the cutter

Only the Millo pro (Fig. 11)

Raising the cutter:

- Turn the adjustment wheel to the left.

Lowering the cutter:

- Turn the adjustment wheel to the right.

13.5 Changing the cutter

**Risk of injury!
The cutters are sharp. Take great care when changing the cutters.**

**Risk of injury!
The cutters heat up during trimming. Allow them to cool before changing them.**

Always use Renfert cutters.

Changing the cutters:

- Switch off the unit.
- Unplug the unit.
- Allow the cutter to cool.

- Remove the work table (Fig. 8).
- Remove the protective cap (Fig. 9).

i The protective cap is more easily removed by pressing it together at the flat surfaces of the clamping ring (Fig. 9).

- Loosen the clamping ring with the Allen key supplied (Fig. 10).

i The Allen key is kept in a holder under the work table.

- Remove the cutter.
- Insert the new cutter fully into the clamping ring.
- Tighten the clamping ring with the Allen key (Fig. 10).

Ensure that the clamping ring has been adequately tightened.

- Replace the protective cap.
- Place the Allen key in the holder.
- Replace the work table.

14 Cleaning / Servicing

Switch off and unplug the unit before cleaning or servicing it.

14.1 Cleaning the housing

Never clean the unit under running water or immerse it in water.

Do not steam clean the unit.

Risk of injury!
Cover sharp cutters if necessary.

Do NOT use cleaning agents that are caustic or abrasive or that contain solvent.

- Switch off the unit.
- Unplug the unit.
- Wipe the housing with a damp cloth.
- Remove the work table (Fig. 8).
- Vacuum out the dust receptacle and wipe with a damp cloth.
- Dry the unit with a soft, clean cloth.

14.2 Cleaning the work table

- Remove the work table (Fig. 8).
- Clean the work table, e.g. steam clean.
- Dry the work table and replace it on the unit.

14.3 Changing the fuse

14.3.1 Millo

Before replacing the fuse, disconnect the unit from the mains power supply by unplugging the power cord from the wall outlet.

Never use fuses with a higher rating.

- Unplug the unit.
- Place the unit on its side.
- Open the fuse holder (G, Fig. 5) and change the fuse.

14.3.2 Millo pro

- Unplug the unit.
- Press in the safety cut-out (F, Fig. 2) on the side of the unit.

15 Replacement parts

You can find components subject to wear and the spare parts on the spare part list in the internet at www.renfert.com/p918.

Enter the following item number: 18040000 / 18050000

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

Serial number and date of manufacturing are shown on the type plate of the unit.

16 Repair

Repairs should only be carried out by qualified personnel.

Repairs to the electrical equipment, which are not specified in these operating instructions, should only be carried out by an electrician.

17 Troubleshooting

Problem	Cause	Remedy
Unit does not operate after switching it on.	<ul style="list-style-type: none"> • Mains cable not plugged in. • No electricity supply to the mains socket. • Unit connected to an extractor that is switched off. 	<ul style="list-style-type: none"> • Plug in the mains cable. • Check the mains socket. • Check the extractor, refer to the extractor instructions if necessary.
Unit switches off during operation.	<ul style="list-style-type: none"> • Motor fuse triggered, as the motor temperature is too high. • Fuse triggered. 	<ul style="list-style-type: none"> • Unplug the unit and allow it to cool for approx. 1.5 hrs. • Use an extractor (see Section 12.3). • Switch on the safety cut-out / replace fuse (see Section 14.3). • If the fuse triggers again, return the unit for repair.
Heavy dust accumulation on the work table.	<ul style="list-style-type: none"> • Extractor not used. • Extraction power too low. 	<ul style="list-style-type: none"> • Use an extractor (see Section 12.3). • Check the filter of the extractor. • Check that the suction hose is not bent. • Check that the suction hose is fitted securely to the connector (D, Figure 2, 4). • Use a higher extraction setting.
Fuse triggers.	<ul style="list-style-type: none"> • Load of the connected extractor too high, adhere to the max. connected load (see technical data). 	<ul style="list-style-type: none"> • Connect the extractor separately to the electricity supply.
Reduced removal or application of greater pressure required.	<ul style="list-style-type: none"> • Cutter blunt 	<ul style="list-style-type: none"> • Change the cutter (see Section 13.5).
Cutter slips out during operation.	<ul style="list-style-type: none"> • Cutter not inserted far enough into the clamping ring. • Clamping ring not tightened enough. 	<ul style="list-style-type: none"> • Insert the cutter fully into the clamping ring. • Tighten the clamping ring further.
Cutter cannot be inserted.	<ul style="list-style-type: none"> • Original Renfert cutter not used. 	<ul style="list-style-type: none"> • Always use original Renfert cutters.
Cutter cannot be removed.	<ul style="list-style-type: none"> • Clamping screw not loosened enough. • Corrosion on the cutter shank. 	<ul style="list-style-type: none"> • Loosen the clamping screw further. • Return the unit for repair.
Clamping screw cannot be loosened.	<ul style="list-style-type: none"> • Screw overtightened. 	<ul style="list-style-type: none"> • Return the unit for repair.
Connected extractor does not operate.	<ul style="list-style-type: none"> • Extractor not switched on. • Extractor is in automatic operation mode. • No electricity supply to the plug-in socket. 	<ul style="list-style-type: none"> • Switch on the extractor. • Set the extractor to continuous operation. • Return the unit for repair.

18 Technical data

Rated voltage:	220 - 240 V, 50 / 60 Hz 120 V, 60 Hz
Power consumption:	120 VA (230 V) / 130 VA (120 V)
max. permitted connected load of the plug-in socket (only Millo pro):	2200 VA (230 V) / 1300 VA (120 V)
Total connected load (only Millo pro) max.:	2320 VA (230 V) / 1430 VA (120 V)
Unit fuse (only Millo):	T1AL, 250 VAC (230 V) / T2AL, 250 VAC (120 V)
Rated speed:	2770 / 3350 U/min (50 / 60 Hz)
LpA *) (idling):	< 70 db(A)
Dimensions, without cutter (W x H x D):	214 x 207 x 252 mm [8.4 x 8.2 x 9.9 inch]
Weight, approx.:	5,7 kg [12.6 lbs]

*) sound pressure level in accord. with EN ISO 11202

19 Guarantee

Renfert gives a **3-year guarantee** on all parts of the Millo / Millo pro provided that the product is used in accordance with the operating instructions. The original sales invoice of the specialist depot is required for a claim under guarantee.

Parts which are subject to natural wear and tear (wear parts) and consumables are excluded from the guarantee. These parts are marked in the spare part list.

The guarantee is rendered void in the case of incorrect use, non-adherence to the operating, cleaning, servicing or connection instructions, repairs carried out by the owner or repairs that are not carried out by a specialist firm, use of replacement parts from another manufacturer or unusual or unapproved uses not specified in the user guidelines.

Successful claims under guarantee do not extend the guarantee period.

20 Disposal instructions

The unit should be disposed of by a specialist firm. The specialist firm should be informed of any hazardous residue in the unit.

20.1 Disposal instructions for countries in the EU

To conserve and protect the environment, prevent environmental pollution and improve the recycling of raw materials, the European Commission adopted a directive that requires the manufacturer to accept the return of electrical and electronic units for proper disposal or recycling.

Within the European Union units with this symbol should not therefore be disposed of in unsorted domestic waste:

For more information regarding proper disposal please apply at your local authorities.

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declare under our sole responsibility, that the product(s)

Name	Part No.
Millo pro	1805 0000 / 1805 0000 03
Millo	1804 0000 / 1804 0000 03

is (are) in conformity with the relevant regulatory requirements by compliance with the UK designated standards.

UK legislation

Electromagnetic Compatibility Regulations 2016

Supply of Machinery (Safety) Regulations 2008

The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

Designated standards

EN 61326-1:2013

EN ISO 12100:2010

EN 61010-1:2010/A1:2019/AC:2019-04

EN IEC 63000:2018

Tilo Burgbacher,
Engineering Director

Hilzingen, 02.11.2022

The Technical Documentation has been retained by Renfert GmbH

**Hochaktuell und ausführlich auf ...
Up to date and in detail at ...
Actualisé et détaillé sous ...
Aggiornato e dettagliato su ...
La máxima actualidad y detalle en ...
Актуально и подробно на ...**

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Renfert GmbH, Industriegebiet, 78247 Hilzingen / Germany
Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany