

**Renfert**

making work easy



Genius evo

**NEW** Bionic Hair brush

# NEW Genius evo

## Bionic Hair brush

Genius evo expands the popular Genius brush range with a brush featuring Bionic Hair. Offering integrated protection in the handle, this premium quality brush impresses with its modern design. However, the Bionic Hair makes it truly special. The special molecular structure of the bristles allows moisture to be absorbed and is distinguished by outstanding longevity. The fineness, elasticity, and resilience of the brush tip allow targeted application of ceramic.

### Advantages

- Innovative brush tip protection thanks to a smooth rotating mechanism in the brush handle.
- Precision work and optimal moisture management due to the special molecular structure of the Bionic Hair fibers.
- Exceptionally long-lasting and cost-effective thanks to easily replaceable brush inserts.
- Consistent haptics for all handle positions.

### Details

- The ergonomic design and low weight ensure easy, relaxed handling during layering.
- The Color Guide ring not only highlights the brush size in color but also includes integrated rollaway protection.

### Tip

If the brush tip bends inadvertently (e.g., due to incorrect storage), it can be held briefly in hot water at approx. 90°C. This automatically causes the bristles of the brush to straighten so that the brush tip returns to its original shape.



An all time classic among the ceramic brushes with Bionic Hair

### Ordering information

Genius evo Set	with 4 brush sizes (2, 4, 6, 8)	No. 17242000
[1] Genius evo size 2	1 piece	No. 17242002
[2] Genius evo size 4	1 piece	No. 17242004
[3] Genius evo size 6	1 piece	No. 17242006
[4] Genius evo size 8	1 piece	No. 17242008
Genius evo refill, size 2	2 pieces	No. 17242112
Genius evo refill, size 4	2 pieces	No. 17242114
Genius evo refill, size 6	2 pieces	No. 17242116
Genius evo refill, size 8	2 pieces	No. 17242118



Integrated protective function thanks to rotating mechanism

