




Renfert Polishing Pastes

Overview of Renfert Polishing Pastes with recommended Polishing Tools

Item		Recommended Polishing Tool				
Polishing Paste	with Dia- mond Grit	Material/ Range of Application	Processing/ Preliminary Polishing	Polishing	High Lustre Polishing	
Renfert Polish hybrid materials No. 510 3000	✓	Hybrid materials				
Renfert Polish all-in-one No. 5100100	✓	Lithium disilicate, zirconium dioxide, veneering ceramics, high performance polymers				
Renfert Polish ZrO ₂ No. 510 1100	✓	Zirconium dioxide				
Renfert Polish LiSi ₂ No. 510 2100	✓	Lithium disilicate				
Kohinoor L No. 5160001	✓	Veneering ceramic				
Brinell L No. 5190001	✓	Precious PFM alloys, ceramic				
Opal L No. 5200001		Veneering composites				
Saphir No. 515 0000		Precious metal alloys				
Saphir No. 515 0000		Cobalt-chrome alloys				
Universal polishing paste No. 513 1000		PMMA, resins				

*Included in Brinell L No. 5190001 scope of delivery

Polishers for Handpiece



Polishers for Lathe

