

SILENT flow sensor

NEW Sensor module for workbench extraction units



Certified sensor module as an upgrade for workplace safety. A component to enhance safety for all the users not only those at NPM workplaces. The module prevents unnoticed loss of suction power and the resulting risk of health hazards.

Advantages

- Health protection ensured by monitoring the volume flow at the suction mouth.
- Acoustic and optical notification when the flow falls below the limit value, i.e. work safely at all times under optimum conditions!
- App-ready - updates and additional functions can be performed via the Renfert CONNECT app.

Details

- Easy to retrofit for certified Renfert extraction systems.
- Module for hazard analysis at workplaces (e.g. cobalt, non-precious metal).
- Easy cleaning option.



Compact mounting under the workbench

Technical data

Permissible mains voltage	100–240 V
Permissible mains frequency	50/60 Hz
Power consumption (<)	1 W // 1 VA
Weight	800 g // 28 oz
Dimensions (W x H x D)	107 x 200 x 95 mm // 4.2 x 7.9 x 3.7"
Ø suction fittings interior	45 mm // 2"
Ø suction fittings exterior	53 mm // 2"

Ordering information

SILENT flow sensor 100-240 V	No. 29140000
SILENT flow sensor 100-240 V Universal adapter	No. 29141000